

RE-800KD

Omron Compressor
Nebulizer for Kids



Target Product Specification

		RE-C800KD	Ref. NE-C801
Product Image			
Size	W x H x D	4.5" × 3.4" × 1.7" 115 × 85 × 43 [mm]	5.6" × 3.9" × 2.8" 142 × 98 × 72 [mm]
	Cubic Volume [W x H x D]	420 [cm³]	1002 [cm ³]
Weight		180 [g]	270 [g]
Sound Level		45 or less [dBA]	46 [dBA]
Worldwide Use		Yes (w/ AC Adapter)	Yes (w/ AC Adapter)
Nebulization Rate		0.3 [ml/min]	0.3 [ml/min]
Particle Size MMAD		Approx. 3 μm	Approx. 3 μm
Contents		Compressor × 1 Nebulizer Kit × 1 Air Tubing × 1 Mouthpiece × 1 Air Filter × 5 AC Adapter × 1 Storage Bag × 1 Child Mask × 1 Toy × 1	Compressor × 1 Nebulizer Kit × 1 Air Tubing × 1 Mouthpiece × 1 Air Filter × 5 AC Adapter × 1 Storage Bag × 1

Item	OHI Model #	RE-C800KD
Nebkit (Neb Cup+Tubing +Mouthpiece)	C803NEB	X
Mouthpiece	C910	X
Nosepiece	C911	X
Air Filter	C30FL	X
Adult Mask	9920	X
Child Mask	9921	X
AC Adapter	HEM-ADPTW5 (Same as BPM)	X

Toy and storage bag will be set up as a maintenance part

Pediatric Compressor Nebulizer (RE-C800KD)



Our smallest, lightest, most convenient tabletop nebulizer – ever.

Quieter than traditional compressor nebulizers

- Approximately 45dB \approx Sound level of your refrigerator

Less weight for easy transportation

- Approximately 6.3 oz (compressor only)
= 5X lighter than most other tabletop compressor models

Ideal particle size for efficient treatment

- MMAD is approximately 3 μ m (based on EN13544-1:2007)

3 years warranty

Reimbursement Pending

*** Available Mid-September 2020**

