

Pulmo-Aide[®] Compressor Nebulizer System

The Gold Standard for High Performance Medication Delivery



Pulmo-Aide Compressor
Nebulizer System
Model 5650D

The DeVilbiss Pulmo-Aide has been providing **effective aerosol treatments for more than thirty years**. Modern styling, durability and ease of use make this unit **perfect for in-home asthma treatments**. The Pulmo-Aide comes ready to use with all the attachments - filter, tubing, and nebulizer.

The **Pulmo-Aide is the only approved compressor for use with TOBI** (tobramycin solution for inhalation). The DeVilbiss Pulmo-Aide is unmatched in quality, value and performance.

Performance Features

- Maximum compressor pressure – 35 psig or greater
- Compressor free air flow – 9 lpm or greater
- 5-year limited warranty
- Longest lasting compressor on the market with an average product life of 10 years
- Unmatched quality and durability - Great for nursing homes and long-term users

Design Features

- Metal frame motor construction
- Unique 4-pole motor design extends life of the unit
- Large built-in compartment for storing nebulizers and tubing
- Built-in nebulizer holder is a convenient place to store your nebulizer if your treatment is interrupted.

Available Models

5650D	Pulmo-Aide Compressor Nebulizer System 115VAC, 2 prong plug
5650H	Pulmo-Aide Compressor Nebulizer System 230VAC, 3 prong plug
5650C	Pulmo-Aide Compressor Nebulizer System CSA approved

Parts and Accessories

DEVILBISS PART #	PRODUCT DESCRIPTION
800D	Reusable Nebulizer – blue (1 each)
800C	Reusable Nebulizer – purple/green (1 each)
4650D-621	Disposable Nebulizer (50 pack)
AAM12	Adult Aerosol Mask (12 pack)
AAM20	Adult Aerosol Mask (50 pack)
PAM10	Pediatric Aerosol Mask (50 pack)
DLI050	Pediatric Children's Character Mask (50 pack)
5650D-601	Compressor Filter (package of 25)
5650D-602	Compressor Filter (package of 5)
5650D-612	Filter Kits (6 each – caps and filters)
5650D-611	Gray Hose Tubing Connector (package of 6)

800 Reusable Jet Nebulizer

- MMAD 3.6 – 4 microns
- Approximately 75% of particles are below 5 microns in size
- Reduces treatment time nearly in half – 3cc treatment in 7 – 9 minutes
- 6cc capacity
- Dishwasher safe to promote frequent cleaning
- Ergonomic design fits comfortably in the hand – even with limited functionality
- 4 components for ease of assembly/disassembly
- Comfortable mouthpiece design
- Standard 22mm fitting, no adapters required
- Can be used inline with a ventilator
- 1 year warranty
- Tested for use with Albuterol, Pulmicort (budesonidel, Ipratropium bromide)
- Suitable for use with most common liquid medications

Technical Specifications

Built to meet the requirements of HCPCS code: E0570		Filter Replacement:	No tools required
Product Dimensions:	10.5H X 6.3W X 10.1D	Warranty:	5 Years
Product Weight:	7.1 lbs.	Safety Approvals:	UL1431
Maximum Air Pressure:	30+ psig	Motor/Pump Assembly:	Diaphragm, Primarily Metal Parts
Maximum Liter Flow:	9+ lpm	Neb Storage:	Yes
MMAD:	5 microns or less	Master Carton Quantity:	4
Sound:	51 dBA		



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com