

PulmoMate® Compressor Nebulizer System

Robust Performance in a Convenient Package



PulmoMate Compressor
Nebulizer System
Model 4650D

The DeVilbiss PulmoMate keeps pace with the world by providing an unbeatable combination of quality and value. The sleek case design fits easily into luggage, backpacks and diaper bags for easy transport and includes a convenient storage compartment for housing tubing and nebulizers

The DeVilbiss PulmoMate features an updated compressor motor for long-term durability and performance. This system delivers a consistent and optimal pressure flow to any nebulizer for the most effective aerosol treatment.

Performance Features

- Maximum compressor pressure – 35 psig or greater
- Compressor free air flow – 8 lpm
- Robust quality and durability
- Consistent performance over the life of the unit
- 5-year limited warranty

Design Features

- Metal frame motor construction
- Integrated handle for easy carrying
- Large built-in compartment for storing nebulizers and tubing
- Intuitive design is easy to use and maintain

Available Models

46502	PulmoMate Compressor Nebulizer System 115VAC, 2 prong plug
4650D	PulmoMate Compressor Nebulizer System 230VAC, 3 prong plug

Parts and Accessories

DEVILBISS PART #	PRODUCT DESCRIPTION
800D	Reusable Nebulizer – blue (1 each)
800C	Reusable Nebulizer – purple/green (1 each)
4650D-621	Disposable Nebulizer (50 pack)
AAM12	Adult Aerosol Mask (12 pack)
AAM20	Adult Aerosol Mask (50 pack)
PAM10	Pediatric Aerosol Mask (50 pack)
DLI050	Pediatric Children's Character Mask (50 pack)
5650D-601	Compressor Filter (package of 25)
5650D-602	Compressor Filter (package of 5)
5650D-612	Filter Kits (6 each – caps and filters)
5650D-611	Gray Hose Tubing Connector (package of 6)

800 Reusable Jet Nebulizer

- MMAD 3.6 – 4 microns
- Approximately 75% of particles are below 5 microns in size
- Reduces treatment time nearly in half – 3cc treatment in 7 – 9 minutes
- 6cc capacity
- Dishwasher safe to promote frequent cleaning
- Ergonomic design fits comfortably in the hand – even with limited functionality
- 4 components for ease of assembly/disassembly
- Comfortable mouthpiece design
- Standard 22mm fitting, no adapters required
- Can be used inline with a ventilator
- 1 year warranty
- Tested for use with Albuterol, Pulmicort (budesonidel, Ipratropium bromide)
- Suitable for use with most common liquid medications

Technical Specifications

Built to meet the requirements of HCPCS code: E0570		Filter Replacement:	No tools required
Product Dimensions:	3.8H X 7.0W X 13.0D	Warranty:	5 Years
Product Weight:	5.3 lbs.	Safety Approvals:	UL1431
Maximum Air Pressure:	35+ psig	Motor/Pump Assembly:	Primarily Metal Parts
Maximum Liter Flow:	8 lpm	Neb Storage:	Yes
MMAD:	5 microns or less	Master Carton Quantity:	6
Sound:	53 dBA		



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA
 800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com