

Plastic Cylinder Wrenches (E-Keys)


EK-1 25/pk
 Also known as "E Keys".
 Used on Post Valve Oxygen Cylinder to Open and Close.



E-Key

Oxygen Regulator Connector Swivel

XM-1 (XMAS Trees), 25/pk
 For Use on Concentrator or Regulator with DISS Fitting in Order to Attach Cannula.



XM-1

Cylinder Wrench

18501 For Use with Small Cylinders used in Most Portable Oxygen Set Ups, 10/pkg



18501

FEATURES:

- Has place to store Viton rubber/ brass sealing washer.
- Used to open post valves.
- Can be used to assist in install or removal of both large wing nut and small t-bar regulators.

Concentrator Water Trap

CON 900 25/bx



CON 900

FEATURES:

- Easy inline downstream installation with small bore tubing. Can be cleaned with detergent or in a vinegar bath.

Tubing Extension Connector

CON 400 Swivel, 50/bg
CON 500 Straight, 50/bg
CON 550 Ridged, 50/bg
CON 700 Y Connector, 10/bg



CON 400



CON 500



CON 550



CON 700

FEATURES:

- Rigid plastic tubing connector allows extension tubing to be connected to present cannula for greater length(s).

COMPARE TO:		
MEDX	Hudson®	Salter®
CON 400	N/A	1220
CON 500	N/A	2005
CON 550	1420	1215
XM-1	2555	1205



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.