NE-C801

OMRON CompAir[®] Compressor Tabletop Nebulizer System

The CompAir Compressor Nebulizer System (NE-C801) is a small, lightweight and convenient tabletop nebulizer that features OMRON advanced design and engineering.





Small, light tabletop nebulizer

weighing less than 9.5 oz. (3.5x lighter than most other compressor models)

Advanced engineering

allows it to administer medication as effectively as tabletop models three times its size

Fast, efficient operation

with a nebulization rate of 0.30mL/minute and particles as small as $3\mu m$ MMAD

Quiet motor

runs at less than 52 dB, similar to the sound of an electric toothbrush*

Affordable and efficient alternative

for anyone with asthma, COPD or other respiratory conditions

Medicine particles with an average size of 3 μ m,

which have a greater probability to penetrate into the deep $\ensuremath{\mathsf{lung}^{\dagger}}$

Product	OMRON Model Number	UPC	NDC	HCPCS Code
CompAir [®] Compressor Tabletop Nebulizer System	NE-C801	073796458010	73796-0458-01	E0570

* dBA comparisons. Center for Hearing and Communication. Retrieved March 13, 2023, from https://chchearing.org/noise/common-environmental-noise-levels/

† Sheth P, Stein S, Myrdal. Factors influencing aerodynamic particle size distribution of suspension pressurized metered dose inhalers. AAPS PharmSciTech, Vol. 16, No. 1, February 2015 (# 2014). DOI: 10.1208/s12249-014-0210-z.

OMRON