

RE-C800KD

# OMRON CompAir® Compressor Nebulizer with Kid's Accessory

The CompAir Compressor Nebulizer with Kid's Accessory is a small, lightweight and tabletop nebulizer designed specifically for kids. The OMRON advanced design features our lightest compressor nebulizer yet—only 6.3 oz—which makes this pediatric nebulizer easily portable.



## Friendly design

features an approachable, kid-friendly animal character

## Portable

our lightest compressor nebulizer yet, weighs in at only 6.3 oz

## Quiet operation

operates with a low noise level, a little louder than a quiet library (~45 dBA)\*

## Efficient delivery

of particles as small as 3µm MMAD, which helps medicine penetrate to the deep lung†

## Reliable

3-year warranty

Product	OMRON Model Number	UPC	NDC	HCPCS Code
CompAir® Compressor Nebulizer with Kid's Accessory	RE-C800KD	073796800147	73796-0800-14	E0570

\* dBA comparisons. Center for Hearing and Communication. Retrieved March 13, 2023, from <https://chchearing.org/noise/common-environmental-noise-levels/>

† Sheth P, Stein S, Myrdal. Factors influencing aerodynamic particle size distribution of suspension pressurized metered dose inhalers. AAPS PharmSciTech, Vol. 16, No. 1, February 2015 (# 2014). DOI: 10.1208/s12249-014-0210-z.