

HCPCS
A7000



610-12B

COMPARE TO

DeVilbiss® 7305D-632

CHAD
drive

NEW

1P
ACK

HCPCS

A7000

A7002

evo



22330

COMPARE TO

DeVilbiss® 73050-633

800cc Disposable Suction Canister

Item #	Product Description	UOM
610-12B	800cc Disposable Suction Canister	12/cs
610-48BP	800cc Disposable Suction Canister	48/cs

- 800cc disposable suction canister
- Raised etched cc increments
- Lid comes with auto shut-off float valve
- Universal ports and size for use with most suction machines.

800cc Disposable Suction Canister Kit

Item #	Product Description	UOM
22330	800cc Disposable Suction Canister Kit	1/ea
22330-12	800cc Disposable Suction Canister Kit	12/cs

- 800cc disposable suction canister kit.
- 10" tubing, elbow connector, 72" suction tubing, hydrophobic filter.

accessories		UOM
610-12B	800cc Disposable Suction Canister	12/cs
SUCP TUBING 10D	Soft, Flexible, Medical Grade Silicon Suction Tubing, 10"	1/ea
SUCP TUBING 72	Suction Tubing with Blue Tip, 72"	1/ea
18600-ELBOW	Universal Flexible Elbow for Suction Machines	20/bx
18600-FILTER	Hydrophobic Filter for Suction Machines	3/bx