

Pediatric Compressor Nebulizer

Item#

NEB300BEAR

Comes with Disposable Neb Kit, Pediatric Aerosol Mask, & Filters
HCPCS# E0570 + A7003 + A7015 + A7013

NEB300BEAR-R

Comes with Disposable and Reusable Neb Kit, Pediatric Aerosol Mask, & Filters
HCPCS# E0570 + A7003 + A7005 + A7015 + A7013

Neb kits include nebulizer bottle, t-piece, mouthpiece, 6in corrugated reservoir tube, and 7ft supply tube.



Great value

**FAST
TREATMENT
TIMES**

Specifications

Weight	2.75 lbs
Input	120V/60Hz
Current	1.5A
Medication Capacity	6ml
Particle Size	Respirable Fraction 0.5 to 5um 0.5um 80%
MMAD	2.0um
Average Nebulization Rate	0.3-0.4ml/min
Noise Level	Below 60dBA
Compressor Pressure Range	30-45psi
Operation Pressure Range	12-14.5psi (0.85-1.1 bar)
Liter Flow Range	8-10LPM
Operation Temperature Range	50°F-104°F (10°C-40°C)
Operation Humidity Range	10 to 95% RH
Storage Temperature Range	14°F-149°F (-10°C-65°C)
Storage Humidity Range	10 to 95% RH

Comes in an attractive retail box

