



CONSISTENT

MEDICATION DELIVERY

VORTEX[®]
Non Electrostatic Holding Chamber



The next generation of PARI VORTEX[®], now even more user friendly, helps minimize the technique and coordination needed for effective inhaler treatments. Patients of all ages are assured a more efficient and consistent medication dose from their MDI for optimal treatments.

www.pari.com/vortex

NON-ELECTROSTATIC HOLDING CHAMBER



PARI VORTEX® Pediatric



Felix Frog Mask
(age 3+ years)



Chloe Ladybug Mask
(age 1-2 years)

Configurations


	Part #	HCPCS #
VORTEX Non-Electrostatic Valved Holding Chamber	051F7000	A4627
VORTEX with Toddler Ladybug Mask	051F7100	A4627
VORTEX With Child Frog Mask	051F7200	A4627

Configurations


	Part #	HCPCS #
VORTEX Hospital Package <i>Case Qty. 50</i>	051H7400P50	A4627
Chloe Ladybug Mask	044F5110	A7015
Felix Frog Mask	044F5210	A7015
Adult Mask	044F7247	A7015

THE DIFFERENCE IS IN THE DELIVERY¹

Vortex significantly reduces oropharyngeal deposition, which may cause unwanted side effects, vs. MDI alone.



Vortex delivers a higher respirable drug dose than MDI alone and OptiChamber Diamond.



VORTEX Grows With Your Patients

From infant to adult, Vortex offers a series of facemasks for various age groups, eliminating the need to purchase a new holding chamber when changes in mask size are needed.



¹Measurements according to Ph. Eur. 2.9.18 and USP 601: Next Generation Impactor; const. flow 30l/min; 50 ± 5 % relative humidity and 23 ± 2 °C; 3 test devices each measured in triplicates (n=3x3=9); pMDI Sultano[®] (100µg Salbutamol, GlaxoSmithKline); *All registered marks and trademarks belong to their respective companies
²US Patent Number: 7,562,656



PARI Respiratory Equipment, Inc.
Tel: 1-800-FAST-NEB (327-8632)
product.info@pari.com
www.pari.com | Midlothian, VA 23112