



Carex Health Brands



**Mobility & Comfort for
your active lifestyle**

Crosstour™ Rolling Walker

Provides mobility assistance for
increased independence



1 - Part Check



Frame
(backrest and wheels attached by brake cable)



Front Wheels x2



Storage Bag

2 - Attach Wheels Using Pop-pins



3 - Unfold Walker Frame



4 - Attach Backrest



DO NOT CROSS BRAKE CABLES WHEN ATTACHING WHEELS AND BACKREST



5 - Adjust Handle Bars Using Pop-pins



6 - Attach the Storage Bag

