



**SLIP RESISTANT
PADS TO REDUCE
MOVEMENT DURING
TRANSFER**

**ASSISTS WITH LATERAL
SIDWAYS MOVEMENT
OF A USER IN SITTING
POSITION FROM ONE PIECE
OF FURNITURE TO ANOTHER**

FEATURES

- Transfer board assists with lateral sideways movement of a user in sitting position from one piece of furniture to another.
- Waterproof transfer board can be used in a wet area
- Two slip resistant pads to reduce movement during transfer
- The arc shape enables a variety of positioning options
- Can be used for transfers from bed to chair, wheelchair, shower commode, etc.
- 60° angle lead in top edge on transfer board

SPECIFICATIONS

Overall Dimensions:

L: 710mm x W: 220mm x H: 10mm

Product weight: **1.5kg**

Max. Load: **130kg** - when gap is 200 - 300mm

180kg - when gap is less than 200mm

Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180

E: orders@for-de.com.au

1 Business way Malaga Western Australia 6090

Part of the **For-de Group Pty Ltd**

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.