



Eagle  
Dental


Feel the Difference



Carbide Cutters

# TC CARBIDE CUTTERS

## Cut 10 E Standard Cross Cut




Cat. No.	0110.023HP	0210.023HP	0410.023HP	0610.023HP	0710.023HP	0910.023HP	1010.023HP	1610.014HP	1810.023HP
Size $\varnothing$ 1/16mm	023	023	023	023	023	023	023	014	023
Thickness mm	14.0	16.0	16.0	17.0	8.0	5.5	7.0	3.1	4.0

HP


**Cut 10 E:** For coarse abrasion on chrome-cobalt alloys, precious metal alloys, chrome-nickel alloys, model cast alloys, and acrylics.

Pack of 1



5110.060HP	5210.060HP	5310.060HP	5410.060HP	5510.060HP	5610.045HP	5710.040HP	5810.060HP	5910.060HP	6110.040HP	6210.040HP
060	060	060	060	060	045	040	060	060	040	040
12.2	12.2	12.2	13.9	15.0	12.9	14.1	12.5	13.9	8.9	7.8

Pack of 1




6310.060HP	6410.031HP	7010.060HP	7110.023HP	7110.027HP	7110.040HP	7110.050HP	7210.040HP	7210.060HP	7710.023HP	7710.029HP
060	031	060	023	027	040	050	040	060	023	029
9.5	11.5	11.0	2.0	2.4	3.5	4.5	11.5	14.2	5.0	5.0

Pack of 1

# TC CARBIDE CUTTERS

## Cut 20 EF Fine Cross Cut




Cat. No.	0120.023HP	0220.023HP	0420.023HP	1220.016HP	1620.014HP	1820.023HP	5220.060HP	5320.060HP
Size $\varnothing$ 1/16mm	023	023	023	016	014	023	060	060
Thickness mm	14.0	16.0	16.0	8.0	3.1	4.0	12.0	12.2

HP

**Cut 20 EF:** Suitable for all dental materials. Enables accurate operation on any surface.

Pack of 1




5420.060HP	5620.045HP	5720.040HP	5820.060HP	6120.040HP	6320.060HP	6420.031HP	7020.060HP	7120.023HP	7220.060HP	7720.023HP	7720.029HP
060	045	040	060	040	060	031	060	023	060	023	029
14.0	13.0	14.0	12.5	9.0	9.5	11.5	11.0	2.0	14.0	3.0	5.0

Pack of 1

## TC CARBIDE CUTTERS


### Cut 24 SPF Fine Spiral Cross Cut



Cat. No.	5424.060HP	6924.040HP	7224.040HP	7724.023HP	1324.008HP
Size $\varnothing$ $\frac{1}{16}$ mm	060	040	040	023	008
Thickness mm	13.9	13.7	11.5	5.0	4.2

Pack of 1


### Cut 25 FST Fine Straight Cross Cut



Cat. No.	0825.023HP	6925.040HP	7225.060HP
Size $\varnothing$ $\frac{1}{16}$ mm	023	040	060
Thickness mm	8.0	13.7	14.2

Pack of 1

### Cut 30 SF Superfine Cross Cut



Cat. No.	0130.023HP	0830.023HP	0930.023HP	1030.023HP	1130.023HP	1230.016HP	1830.023HP	5730.040HP	7730.023HP	7730.029HP
Size $\varnothing$ $\frac{1}{16}$ mm	023	023	023	023	023	016	023	040	023	029
Thickness mm	14.0	8.0	5.5	7.0	5.5	8.0	4.0	14.0	5.0	5.0

Pack of 1

HP

**Cut 24 SPF:** An exceptionally fine cutting surface on all alloys - especially on titanium since the special cut pattern avoids clogging.

HP

**Cut 25 FST:** Excellent cutting performance and facilitated chip removal on tough materials like titanium and NP alloys, P metals, plaster of Paris, model castings, veneer, and prosthetic resins.


**Cut 26 FST:** Suitable for titanium and titanium alloys.

HP

**Cut 30 SF:** For fine cutting of gold, alloys, and all versions of resins and composites. Due to the special cut, also ideally suited for ceramics to avoid occurrence of crackles.

## TC CARBIDE CUTTERS

### Cut 50 SG Coarse Cross Cut



Cat. No.	5250.060HP	5450.060HP	5550.060HP	5650.045HP	5850.060HP	7050.060HP
Size $\varnothing$ $\frac{1}{16}$ mm	060	060	060	045	060	060
Thickness mm	12.2	13.9	15.0	12.9	12.5	11.0

Pack of 1

HP

**Cut 50 SG:** For efficient reduction of resins including tray materials as well as plaster.


### Cut 51 GE Medium Coarse Cross Cut



Cat. No.	5251.060HP	5451.060HP
Size $\varnothing$ $\frac{1}{16}$ mm	060	060
Thickness mm	12.2	13.9

Pack of 1

### Cut 55 XT Super Coarse Cross Cut



Cat. No.	5455.060HP	6055.070HP
Size $\varnothing$ $\frac{1}{16}$ mm	060	070
Thickness mm	13.9	14.2


Pack of 1

HP

**Cut 51 GE:** For bulk reduction on dry plaster.

**Cut 55 XT:** For gross reduction of resins and plaster as well as resin-based tray materials.

### Cut 75 GSQ Coarse Straight Blade Cross Cut



Cat. No.	0175.023HP	7275.060HP	7775.040HP
Size $\varnothing$ $\frac{1}{16}$ mm	023	060	040
Thickness mm	14.0	14.2	9.0

Pack of 1

HP

**Cut 75 GSQ:** Plain tothing with transverse section. For trimming of soft acrylics and soft relinings.



Eagle  
Dental

Feel the Difference