





Eagle  
Dental

Feel The Difference



## DEGRANULATION KIT

Designed by Dr. Efraim Kfir



### Kit Features & Benefits:


- Designed morphologically to tear out all the remnants of soft tissue which are well attached into the bone and grab them out without harming the bone.
- The debridement process is done quickly and effectively.
- The drills are more precisely and do not grab the bone compare to tungsten burs.
- The drills are connected to physiodispenser working in sterile environment.
- The speed can be controlled comparing to high speed depends the kind of the procedure needed and preventing emphysema.
- Multi use unique suction which is designed exclusively for excellent direct visualization of the infected area (is always the best technique) to ensure optimal debridement.

	1 mm	2.5 mm	3 mm	3.5 mm
Code	Strauss G 1 mm	Strauss G 2.5 mm	Strauss G 3.0 mm	Strauss G 3.5 mm
Head size	1 mm	2.5 mm	3.0 mm	3.5 mm
Total length	31 mm	31 mm	31 mm	31 mm




### Advantages:

- Optimal debridement of the inflammatory tissue is achieved which cannot be done manually by surgical spoon curette.
- Multi task drills can be used in implantology, periodontal surgery, endodontic surgery, and root removal.
- Maximum safety when working near anatomical structure
- Multiuse drills (more than 100 procedures)




# SAFE SINUS LIFT KIT

Designed by Dr. Efraim Kfir



## Kit Features & Benefits:

- High durability burs – can be used for more than 100 cases.
- Safe sinus uplift procedure.
- Great visibility of the sinus membrane.
- No damage to the membrane.
- Fast and easy process.
- Accurate and repeatable results.



## Main characteristic:

- Can be connected to a physio-dispenser for maximum sterility.
- Can be used in implantation, periodontic surgery and endo removal.
- Special a-traumatic diamond burs with calibrated stopper.
- The special burs allows for a fast and accurate penetration of the hard tissue without harming the membrane.
- Unique suction which clears all fluids but doesn't remove bone substitutes.
- Custom made membrane separator allowing an easy detachments of the membrane.

## The Kit Contents:

- The Kit Contains 6 Specialized Burs: All burs are 3 mm in diameter, equipped with built-in calibrated stoppers allowing penetrating depths ranging from 5-10 mm (with 1 mm increments). The unique design enables the practitioner to drill safely through the antral floor, in order to expose the sinus (schneiderian) membrane.
- Dedicated Suction Syringe.
- Sinus Membrane separator: This unique curette is designed for circumferential antral membrane detachment and elevation all around the osteotomy site.



## Zekrya

A unique bur that has a cutting tip and 6 helical flutes making them ideal for sectioning an impacted wisdom tooth prior to its extraction, separating roots and removing.



Cat. No.	Zekrya23	Zekrya28
Head Size	016	016
Head Length	11	11
Total length	23	28
FG	✳	✳



**100% Control**

## Lindemann

The Lindemann bur is designed for aggressive cutting of bone structure. The bur is a one-piece Tungsten Carbide. The special geometry allows an efficient cutting of hard tissue or bone tissue while removing the debris.



Cat. No.	Lindemann161	Lindemann162
Head Size	016	016
Head Length	10	10
Total length	25	25
FG	✳	✳
RA	✳	✳
Blade Geometry	Plain Cut	Cross Cut


Available in HP upon request

## Trephine Bur

Surgical burs for bone grafting

Cat. No.	Size	Ø Outer	Ø Inner	Drilling Depth
TB3X4	040	4,0 mm	3,0 mm	max. 16,0 mm
TB4X5	050	5,0 mm	4,0 mm	max. 16,0 mm
TB5X6	060	6,0 mm	5,0 mm	max. 16,0 mm


- The burs are made from high quality surgical grade stainless steel
- Laser depth marking for easy identification
- Different shapes and sizes are available upon request





## Tissue Trimming

- Designed for soft-tissue surgery, the ceramic soft tissue trimmer bur is a rotary instrument that can be used as an alternative to using scalpel, electro-surgery and surgical blades.
- Using the tissue trimmer provides minimum bleeding.
- The biocompatible hard oxide points are made to cut tissue without burning.



Cat. No.	FG-CER	FG-CEROS
Head Size	015	015
Head Length	5	8
Total length	22.5	25
FG	✳	✳

**Tissue Trimmer**

## 6 Blades Operative & Surgical

Cat. No.	1/4	1/2	1	2	3	4	5	6	7	8
Head Size	005	007	008	010	012	014	016	018	021	023
Head Length	0.4	0.6	0.6	0.8	1	1.1	1.2	1.5	1.9	1.9
FGOS	✳	✳	✳	✳	✳	✳		✳		✳
HP	✳	✳	✳	✳	✳	✳	✳	✳	✳	✳

**Round**

Cat. No.	1702
Head Size	016
Head Length	4.1
FGOS	✳
HP	✳

**Round End Cross Cut Tapered Fissure**

Cat. No.	556	557	558	559	560
Head Size	009	010	012	014	016
Head Length	4	4.5	4.5	4.5	4.5
FGOS	✳	✳	✳	✳	
HP	✳	✳	✳	✳	✳

**Cross Cut Fissure**

Cat. No.	701	702	703
Head Size	012	016	021
Head Length	4.5	4.5	5
FGOS	✳	✳	✳
HP		✳	✳

**Cross Cut Tapered Fissure**

Cat. No.	CC001-040M	CC001-050M	CC001-060C
Shape	001	001	001
Size	040	050	060
Roughness	M	M	C

**Carbide Cutters**



## Magic Touch



Perfect Burs for fast & smooth grinding of Zirconia, E-Max® and porcelain  
The Safest, the Fastest, the Smoothest Crown & Bridge Adjustment



See demo video

Cat. No.	A3PF	A4PF	B4PF	K2PF	K2RPF	M3PF	C31PF	PR15PF	E12PF	PR2PF	E5PF	T2PF	E7PF
Code	801-018PF	801-022PF	807-018PF	368-023PF	379-023PF	856L-018PF	856-025PF	889-012PF	858-012PF	390-018PF	862-014PF	878K-022PF	885-012PF
ISO No.	001-018	001-022	225-018	257-023	277-023	198-018	197-025	170-012	172-012	254-018	249-014	298-022	130-012
Head Size	018	022	018	023	023	018	025	012	012	018	014	022	012
Head Length			5	5.3	4.4	9	8	4	8	3	8	8	8
Total Length	19	19	19	20	19	21	21.5	20	22	19	21	20.5	21
<b>Magic</b>													

Available in HP upon request



CEREK KIT 1  
For full description see page 8



## FG Diamond Burs



Cat. No.	A1L	A2L	A3L	A4L	A5L
Code	801L-012	801L-014	801L-018	801L-022	801L-032
ISO No.	001-012	001-014	001-018	001-022	001-032
Head Size	012	014	018	022	032
Total Length	24	24	24	24	24
Fine					▼
Medium	▽	▽	▽	▽	▽
Coarse	▽	▽	▽	▽	▽
S. Coarse		▼	▼	▼	▼

Surgical Length Round



Cat. No.	A2XL	A3XL	A4XL
Code	801XL-014	801XL-018	801XL-022
ISO No.	001-014	001-018	001-022
Head Size	014	018	022
Total Length	28	28	28
Fine			
Medium	▽	▽	▽
Coarse			
S. Coarse			

Extra Long Surgical

## Periodontic Burs

Cat. No.	MR10L2	MR10L3	MR015L	MR010L	C2411S	C244
Code	Special	Special	839L-015	839L-010	849-014	849-016
ISO No.	Special	Special	150L-015	150L-010	196-014	196-016
Head Size	012	012	015	010	014	016
Head Length	2	2			2	4
Total Length	24	24	24	24	23.5	23.5
M	▽	▽	▽	▽	▽	
XXF					▼	▼

Periodontic

special custom made burs with laser depth markings





# Eagle Dental



## Crown Lengthening Kit

Profile	Cat. No.	Application
	A4LC	round diamond bur
	A3LC	round diamond bur
	MR10L2M	side cutting diamond bur (narrow)
	MR10L3M	Side cutting diamond bur (narrow)
	FG957	End cutting Carbide bur
	MR015LM	End cutting diamond bur
	MR010LM	End cutting diamond bur
	C2411S/M	Tapered, side cutting diamond bur
	C2411S/XXF	Tapered, side cutting diamond bur
	C244/XXF	Tapered, final polishing diamond bur



CAT #: PERIOKIT1

The new Crown Lengthening kit is designed to be used in an orderly sequence and guide the clinician to achieve the desired outcome precisely, predictably and safely. It includes special custom made burs with laser depth markings. Designed by Dr. Cobi Landsberg and Prof. Nitzan Bichacho.

## Magic Touch Cerec Kit

A complete set for Cerec restoration E-MAX® and Zirconia adjustments.

Profile	Cat. No.	Application
	A3PF	Magic Touch - adjustments / endo access
	E5PF	Magic Touch - adjustments / crown cutting
	K2RPF	Magic Touch - adjustments
	M3PF	Magic Touch - adjustments / crown cutting
	RDRAZ14M	Zirconia smoothing polisher
	RDRAZ14F	Zirconia polishing polisher
	RDRAE14M	Lithium Disilicate smoothing polisher
	RDRAE14M	Lithium Disilicate polishing polisher



CAT#: CERECKIT1

Designed by Dr. Eran Dolev, Department of Oral Rehabilitation, School of Dental Medicine, Tel – Aviv University. ISCD certified CEREC trainer.