

# MAURO

## Ergo-line Relax Chair



# Mauro Ergo-line Relax Chair

Care Quip introduces the new Mauro Ergo-line Relax reclining chair. Featuring a comfortable and supportive recline mechanism that tilts into a zero gravity position, this chair has been specifically designed for the care sector.

## 01. Contoured Backrest

The seat shell has been designed following the shape of the human spine (S-shape). It offers an exceptional support at the level of the lower back and the shoulders. An additional neck support cushion is possible, but certainly not required.

## 02. Relax Mechanism

Less pressure points on the body and a good blood circulation thanks to the Ergo-line relax mechanism that inclines the backrest, seat and legrest simultaneously into a "Zero Gravity" position (\*).

## 03. Central Brake

The chair is extremely manoeuvrable and mobile, thanks to the large castors. The 2 castors at the back legs have central brakes and a directional lock.

## 04. Adjustable Armrests

The lowerable armrests can be adjusted in 4 different positions. In their lowest position, the armrests are on the same level as the seat which facilitates patient transfers.

## 05. Modularity

The relax chair can be composed the way you want it. Various accessories make the chair more complete, such as a removable headrest and a dinner tray.

## 06. Colours

ED7800 Graphite Grey  
ED7810 Brick Red  
ED7820 Midnight Blue



(\* Zero gravity position: this seat/lying position is based on NASA data. A neutral body position converts all muscle power into a biomedical balance. (Source: "From outer space to you", Omni, March 1994)

## Specifications

Length x Depth x Height:	720 x 780-1550 x 1250mm
Seat Height From Footrest:	460mm
Seat Depth:	510mm
Seat Width:	535mm
Height Armrests:	735-530mm
Safe Working Load:	150kg

