



Etac Fix and Cut food preparation

Award-winning ergonomically designed aids for daily living

 **etac**[®]




Creating Possibilities

Fix kitchen grip

This kitchen product keeps the provisions in place, making it possible to work with one hand only or with a double grip to do simple food preparation.

Fix simplifies household chores

Fix preparation board stands firmly on its four friction feet or the enclosed suction feet. The jaws hold a bowl, grater or piece of food in place with ease. There is also a detachable spiked plate for smaller pieces of food. This can be removed and placed with the spikes pointing downwards to obtain a flat surface.



Description	Order Code	Material	Weight
Fix Preparation Board	80501004	Chopping board: Polystyrene Metal parts: Stainless steel Jaws: PC/ABS	1.3 kg

Cut cutting aid

This brace keeps the bread in place, making it possible to work with one hand only or with a double grip to do simple food preparation.

Cut stands firmly and simplifies cutting

Cut is a functional cutting aid which simplifies kitchen work by standing firmly and holding in place whatever is to be cut. By turning the transverse lath, the thickness of the slices is easily changed. The knife fits in the guides, making it possible to cut safe and straight with impaired sight.



Description	Order Code	Material	Weight	Care
Cut Cutting Aid	80501005	Polystyrene	415 g	Dishwasher-safe



 **etac**[®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”