

“The SD2 Wheelchair is an upgrade in both design and features and replaces the old Drive Super Deluxe Wheelchair.”



## SD2 ALUMINIUM WHEELCHAIR

### Product features

- Elegant metallic black aluminium frame
- Double cross brace provides added support to the wheelchair
- Attendant or user accessible parking brakes on rear wheels ensure safety along with attendant cable brakes on transit models
- Improved footrests swing in and out to enable transfers on and off chair, whilst flip-up and removable armrests facilitate side transfers
- Push button adjustment facilitates height and length adjustable arm pads
- Incorporated as standard are deluxe seat and back cushions and a padded tension adjustable backrest for added comfort
- Rear stepper tubes with integrated height adjustable anti-tip wheels facilitate safe curb-climbing for both attendant and independent user
- Half folding back mechanism allows easy storage and transportation
- Quick release rear wheels on Self Propel models
- Black pneumatic rear tyres and puncture proof front tyres are durable and low maintenance
- Chair is fitted with reflectors to improve outdoor visibility at dusk or when dark
- One touch lap belt for added safety
- Features handy rear pocket for storing items
- Comes with calf strap as standard for additional leg support



Code	Width (Open)	Width (Closed)	Length	Height	Seat Width	Seat Length	Seat Height	Weight	Weight Capacity
SD2SP16BLKAU	61cm (24")	28cm (11")	105cm (41")	100cm (39")	40cm (16")	40cm (16")	52cm (20")	18kg (39lb)	135kg (21st)
SD2SP18BLKAU	66cm (26")	28cm (11")	105cm (41")	100cm (39")	45cm (18")	40cm (16")	52cm (20")	18kg (39lb)	135kg (21st)
SD2SP20BLKAU	71cm (28")	28cm (11")	105cm (41")	100cm (39")	52cm (20")	40cm (16")	52cm (20")	19kg (41lb)	135kg (21st)
SD2SP22BLKAU	76cm (30")	28cm (11")	105cm (41")	100cm (39")	56cm (22")	40cm (16")	52cm (20")	19kg (41lb)	135kg (21st)
SD2TS16BLKAU	57cm (22")	28cm (11")	102cm (40")	100cm (39")	40cm (16")	40cm (16")	52cm (20")	17kg (37lb)	135kg (21st)
SD2TS18BLKAU	62cm (25")	28cm (11")	102cm (40")	100cm (39")	45cm (18")	40cm (16")	52cm (20")	17kg (37lb)	135kg (21st)
SD2TS20BLKAU	67cm (26")	28cm (11")	102cm (40")	100cm (39")	52cm (20")	40cm (16")	52cm (20")	17kg (37lb)	135kg (21st)
SD2TS22BLKAU	72cm (28")	28cm (11")	102cm (40")	100cm (39")	56cm (22")	40cm (16")	52cm (20")	18kg (39lb)	135kg (21st)

“Our new and improved XS2 Wheelchair uses a similar frame to the XS Wheelchair but adds many new features”



## XS2 ALUMINIUM WHEELCHAIR

### Product features

- Sturdy, silver aluminium frame is lightweight for transporting, durable and aesthetically pleasing
- Simple folding mechanism, half-folding backrest, removable arms as well as quick-release footrests and self-propelling rear wheels allow space-saving storage and easy transportation
- Padded desk-style armrests provide additional comfort and allow unhindered access to work surfaces
- Footrests swing in and out to enable transfers on and off chair, whilst flip-up armrests facilitate side transfers
- A one-touch lap belt is also supplied on both self-propel and transit models, providing additional stability for user
- Incorporated as standard are a seat cushion, tension adjustable backrest and push-button length adjustment on armrests, allowing chair to be easily adapted to suit comfort preferences
- Long-reach brakes on rear wheels give easy access to attendant and user for increased safety, with transit model also featuring 'immediate-access' cable brakes to reduce speed or safely park chair
- Rear stepper tubes with integrated anti-tip wheels facilitate safe curb-climbing for both attendant and independent user
- Chair is fitted with reflectors to improve outdoor visibility at dusk or when dark
- All-round solid Polyurethane (PU) tyres offer durability and minimum maintenance
- Available in three different seat widths to suit various user sizes



Code	Width (Open)	Width (Closed)	Length	Height	Seat Width	Seat Length	Seat Height	Weight	Weight Capacity
XS2SP16SILAU	59cm (23")	26cm (10")	103cm (40.5")	93cm (37")	40cm (16")	43cm (17")	52cm (21")	17kg (37lb)	135kg (21st)
XS2SP18SILAU	63cm (25")	26cm (10")	103cm (40.5")	93cm (37")	45cm (18")	43cm (17")	52cm (21")	17kg (37lb)	135kg (21st)
XS2SP20SILAU	68cm (27")	26cm (10")	103cm (40.5")	93cm (37")	50cm (20")	43cm (17")	52cm (21")	17kg (37lb)	135kg (21st)
XS2TS16SILAU	59cm (23")	26cm (10")	92cm (36")	93cm (37")	40cm (16")	43cm (17")	52cm (21")	15kg (33lb)	135kg (21st)
XS2TS18SILAU	63cm (25")	26cm (10")	92cm (36")	93cm (37")	45cm (18")	43cm (17")	52cm (21")	15kg (33lb)	135kg (21st)
XS2TS20SIL	68cm (27")	26cm (10")	92cm (36")	93cm (37")	50cm (20")	43cm (17")	52cm (21")	15kg (33lb)	135kg (21st)

Leaflet code: DDVH\_031\_XS2\_&\_SD2\_wheelchair\_Leaflet