

Compact Fold for Portability



Folding Backrest



Swing-away Detachable Footrest



Slowing and Parking Brake



days[®] Escape Wheelchair

With a foldable aluminium frame, the Escape has been specifically designed for ease of transportation and storage. Suitable for both indoors and out, it gives you flexibility wherever you go.

In addition to being portable, these attendant or self-propelled chairs offer both comfort and practicality. With a full-length armrest, padded upholstery and detachable swing-away footrests, they provide all the basic comfort necessities. Peace of mind and low maintenance are provided by the puncture-proof tyres.

Key features

- Folding frame and backrest for compact storage
- Generously padded flame retardant upholstery
- Puncture-proof tyres
- Full length armrest
- Attendant-operated cable brakes
- Swing-away and detachable footrests
- Seat belt included
- Available in 2 colours: Silver Blue, Ruby Red



Specifications


Specifications	Attendant	Self Propelled
Seat Width (on the outside of the tubes)	450mm	450mm
Seat Width (on the inside of the tubes)	400mm	400mm
Seat Depth	400mm	400mm
Seat Height - Front (from the ground)	510mm	510mm
Seat Height - Rear (from the ground)	490mm	490mm
Width Between Arms	400mm	400mm
Backrest Height Above Seat	400mm	400mm
Footrest Length (from seat to the surface of the footrest)	350-450mm	350-450mm
Armrest Height Above Seat	230mm	230mm
External Width (Open)	530mm	640mm
External Width (Folded)	270mm	320mm
External Depth (w/o Footrests)	750mm	800mm
External Depth (w Footrests)	995mm	1056mm
Chair Weight	12kg	15kg
Max User Weight	100kg	115kg

Product Codes

Silver Blue
Ruby Red

PAT-091171727
PAT-091556893

PAT-091566249
PAT-091566256



Silver
Blue

Ruby
Red

1 Year
Product
Warranty

2 Years
Main Frame
Warranty

ATTENDANT
100
kg
Maximum
User Weight

SELF-PROPELLED
115
kg
Maximum
User Weight