

## EZ ADJUST BEDRAIL

BECAUSE IT'S OUR MOST VERSATILE BED RAIL
















## **FEATURES & BENEFITS**

- Bed Rail/Hand Rail Adjustable length allows for use as both a bed rail to prevent falling out of bed or as a hand rail to assist in getting out of bed.
- Adjustable Bed rail easily adjusts from 26" to 34" to 42" after it has been placed on the bed.
- Collapsible Rail folds down to side of bed to allow user to get out of bed or provide space for making bed.
- Reversible Can be placed on either side of the bed
- Dual Safety Strap Secures to bed frame with included safety strap.

## **SPECIFICATIONS**

Width of Rail: Locks in place at 26", 34" & 42" after installed on bed

Height of Rail: 23" high from the base of mattress

Gap Between Rails: 3"

Bed Type: Fits 12"- 18" mattress

Weight of Product: 14 lbs

Weight Capacity: 300 lbs

Assembly: Installs in minutes with 4 bolts and allen

wrench (included)

Package Dimensions: 28" × 21" × 1.5" Case of 4 Dimensions: 28" × 21" × 6"; 60 lbs

## WHAT IS IN THE BOX









WEIGHT CAPACITY



**EXTENDS**