



IPR Handsaws System



- Starter **Yellow** 0.06mm, Single-Sided | Extra-Fine Diamond
Start IPR safely and open the interproximal contact easily



- Opener **Red** 0.12mm, Double-Sided | Medium Diamond
Reduce enamel efficiently on two adjacent teeth



- Widener **Blue** 0.15mm, Double-Sided | Coarse Diamond
Widen the contact and reduce enamel effectively



- Extra-Widener **Green** 0.20mm, Double-Sided Coarse Diamond
Further widen the interproximal contact

IPR Handsaws Individual 8 Pack

8x Yellow 0.06mm

8x Red 0.12mm

8x Blue 0.15mm

8x Green 0.20mm



IPR Handsaws Assorted 8 Pack

1x Yellow 0.06mm

1x Red 0.12mm

1x Blue 0.15mm

1x Green 0.20mm

1x Clear 0.10mm

1x Cyan 0.12mm

1x Purple 0.25mm

1x Brown 0.30mm





IPR Handsaws System



- Single-Sided Opener **Clear** 0.10mm, Single-Sided | Medium-Fine Diamond
Serrated, further opening strip which affords great accuracy



- Single-Sided Widener **Cyan** 0.12mm, Single-Sided | Medium-Coarse Diamond
Initial widener specifically for one surface reduction



- Super-Widener **Purple** 0.25mm, Double-Sided | Coarse Diamond
IPR workhorse as most cases average 0.3mm reduction



- Mega-Widener **Brown** 0.30mm, Double-Sided | Coarse Diamond
More aggressive but with high degree of accuracy and control


IPR Handsaws Individual 8 Pack

8x Clear 0.10mm

8x Cyan 0.12mm

8x Purple 0.25mm

8x Brown 0.30mm



IPR Gauges

Keyring of 6

0.10mm

0.20mm

0.25mm

0.30mm

0.40mm

0.50mm

