

Disposable Dental

Product Guide

one 





The importance of providing and verifying an ultra-safe group of instruments and products for yourself and your staff when a potentially high risk patient needs treatment cannot be overstated.

From initial examinations and restorative procedures to complex crown and bridge procedures, One Time has you protected.

Similarly, if your sterilisation protocols are disrupted due to autoclave breakdown, simply turn to One Time products as your ever reliable backup.

A must have for every practice.



Disposable Handpiece

Quality super-sterilised disposable hand-piece for single use with high risk patients.

Low cost, high-speed handpiece as a backup to autoclave failure.

High torque, water and air cooled, high speed up to 370,000rpm and 65 db noise.

Simple coupling interchange into standard dental units, Midwest 4 hole.

Sold individually or in Packs of 3.

Pre-oil before use. For dental and single use only, do not repeat use and store at ambient temperature.


Visit us @ VIMEO

1122 | EOCA

for footage of our handpiece in use

Disposable Crown and Bridge Preparation Kits

Everything you need to prepare for a crown in the anterior or posterior region in one kit



Anterior Kit

Anterior C&B Preparation Bur Set
10 x Anterior C&B preparation burs.

Disposable Handpiece

Quality super-sterilised disposable handpiece for single use with high risk patients.

Examination Kit

1x Probe, 1x Mirror, 2x Masks,
1x Tweezers, 2x Cotton rolls,
1x Patient bib (with ties).

Posterior Kit

Posterior C&B Preparation Bur Set
10 x Posterior C&B preparation burs.

Disposable Handpiece

Quality super-sterilised disposable handpiece for single use with high risk patients.

Examination Kit

1x Probe, 1x Mirror, 2x Masks,
1x Tweezers, 2x Cotton rolls,
1x Patient bib (with ties).

Restorative Kit

Restorative C&B Preparation Bur Set
10 x Posterior C&B preparation burs.

Disposable Handpiece

Quality super-sterilised disposable handpiece for single use with high risk patients.

Examination Kit




1x Probe, 1x Mirror, 2x Masks,
1x Tweezers, 2x Cotton rolls,
1x Patient bib (with ties).

Restoration Instrument Kit

1x Ball Burnisher, 1x Amalgam
Plugger, 1x Flat Plastic.

Disposable Restoration Kits

Everything you need for every variety of restorative procedure



Disposable Bur Sets

For every case, crown and bridge or restorative. Simply use & discard. Maximum efficacy, minimal risk and cost efficient.



Anterior

Crown and Bridge Preparation Bur Set

10 x Anterior C&B preparation burs.

One Time Use



Posterior

Crown and Bridge Preparation Bur Set

10 x Posterior C&B preparation burs.

One Time Use



Restorative

Crown and Bridge Preparation Bur Set

8 x Restorative C&B preparation burs.

One Time Use



Disposable Restoration Instruments

Small or large, for all conventional restorative procedures: use and discard



Select from;

1. Ball Burnisher
2. Amalgam Plugger
3. Flat Plastic Metal Handle



4. Restoration Instrument Kit which includes all 3 instruments.

Disposable Examination Kit

A complete pack for every exam: completely disposable

- 1x Mirror
- 1x Probe
- 2x Masks
- 1x Tweezers
- 2x Cotton Rolls
- 1x Patient bib with ties

**Single use only, do not repeat use.
Do not use if package is damaged.**



To order or for more information
1122Corp.com

