



1



Insert smartphone  
Tighten thumbscrew

2



Center over camera  
Tighten thumbscrew

# dermLite Universal Smartphone Adapter

## 3Gen



 DermLite and 3Gen are registered trademarks of 3Gen, Inc. | 31521 Rancho Viejo Road | Suite 104 San Juan Capistrano, CA 92675 | U.S.A.  
Tel. +1-949-481-6384 | Fax +1-949-240-7492  
E-Mail: [info@3geninc.com](mailto:info@3geninc.com) | [www.dermLite.com](http://www.dermLite.com)  
©2017 by 3Gen, Inc. | All rights reserved.  
REF A DLUA-1701  
Patents/Trademarks: [www.dermLite.com/patents](http://www.dermLite.com/patents)

EC	REP
----	-----

M. Devices Group, Healthcare Education Centre,  
The Church, Portland Street, Southport, PR8 1HU, UK.  
Tel: +44 1704 544 944 | Fax: +44 1704 544 050



55 mm min.

98 mm max.