



Delivering medical supplies with the right quality, right quantity, at the right time, and at a fair price.

BV MEDICAL®



Customer Service

First-class customer service where people can reach people.



Small Business

We are a small business with a big mission.



Our Vision

To enhance customers' health through a committed working relationship.

About Us

Vision

To enhance customers' health through a committed working relationship.

Mission

We are a small business, with a big mission. We forge a relationship of trust and reliability with our customers through first-class customer service where people can reach people, not a faceless corporation. We rapidly deliver diagnostic, respiratory, and medical instruments with the right quality, right quantity, at the right time, at a fair price. We back that up with solid, generous warranties.

Values

Customer Service Excellence Teamwork **Ethical Focus**

Authorized Reseller

BV Medical® is an Authorized Reseller of Omron®, KaWe®, Kinsa®, PediaPals, Accuro, A&D Medical®, and TecniMed products. We are also a Reseller of Littmann® products.

















Please visit our website at www.bvmedical.com or call our Customer Service Team at (888) 822-8293

Business Hours:

Monday - Friday 8:00 am to 4:00 pm Central



Optimum Series

Optimum aneroid sphygmomanometer with chrome plated manometer and luminescent dial for easy reading. Not made with natural rubber latex.

Available in the following sizes:

Child 7.25" -10.5"

20-108-013-L

Adult 10" - 16"

20-108-014-L

L Adult 13.5" - 20"

20-108-015-L

Thigh 16" - 26"

20-108-016-L



Basic Series

Basic aneroid sphygmomanometer with nylon cuff markings and standard air release valve. Not made with natural rubber latex.

Available in the following sizes:

Child 7.25" -10.5"

20-100-013-L

Adult 10" - 16"

20-100-014-L

L Adult 13.5" - 20"

20-100-015-L

Thigh 16" - 26"

20-100-016-L



Classic Series

Classic aneroid sphygmomanometer with nylon cuff markings and carrying case.

Not made with natural rubber latex.

Available in the following sizes:

Child 7.25" -10.5"

20-208-013-L

Adult 10" - 16"

20-208-014-L

L Adult 13.5" - 20"

20-208-015-L

Thigh 16" - 26"

20-208-016-L

OMRON®





OMRON is a leader in home blood pressure monitors. Clinically proven for accuracy.

BP7900

Uses Bluetooth® technology to let you measure, store, and review 90 readings.

Lightweight, slim, and portable design.

UltraSilent Inflation, so measurements can be taken quietly and discreetly.

Available Upper Arm Units:

OMRON Evolv	B P 7 0 0 0
OMRON 3 Series	B P 7 1 0 0
OMRON 5 Series	B P 7 2 0 0
OMRON 5 Series	B P 7 2 5 0
OMRON 7 Series	B P 7 3 5 0
OMRON 10 Series	B P 7 4 5 0
OMRON BP + EKG	B P 7 9 0 0



OMRON 7 Series Wireless

OMRON is a leader in home blood pressure monitors. Clinically proven for accuracy.

BP6350

Uses Bluetooth® technology to let you measure, store, and review 90 readings.

Lightweight, slim, and portable design.

UltraSilent Inflation, so measurements can be taken quietly and discreetly.

Available Wrist Units:

OMRON 3 Series

BP6100

OMRON 7 Series

BP6350

A&D®



A&D Premium

For over 40 years, A&D has been clinically validating their units to ensure they are accurate and easy-to-use.

UA-789AC

Ideal for people with larger upper arms. The extra large cuff measures 16.5" to 23.6".

Includes an irregular heartbeat indicator.

Simple one-button operation and a large display for easy use.

Available Upper Arm Units:

A&D Premium	U A - 1030T
A&D Basic	U A - 611
A&D Essential	U A - 6 5 1
A&D Deluxe	U A - 6 5 1 B L E
A&D Advanced	U A - 7 0 5
A&D Premium	U A - 7 6 7 F
A&D Pro	U A - 767 P V S
A&D Premium	U A - 787EJ
A&D Premium	U A - 789 A C



A&D Premium

For over 40 years, A&D has been clinically validating their units to ensure they are accurate and easy-to-use.

UB-543

Multi-user devices that provides trusted accuracy and reliability.

Compact for storage and travel.

Includes an irregular heartbeat indicator and pressure rating indicator.

Available Wrist Units:



Classic Stainless Steel

Precisely machined from surgical stainless steel for optimum performance and durability. Dual-head chestpiece. Not made with natural rubber latex.

Available in the following colors:

Black 38-700-014

Burgundy 38-700-054

Hunter Green 38-700-124

Magenta 38-700-164

Navy Blue 38-700-184

Purple 38-700-234

Royal Blue 38-700-254



Dual Sprague Rappaport

Multi-functional dual-head that can be changed out for a pediatric diaphragm, adult diaphragm, pediatric bell, child bell, and adult bell.

Available in the following colors:

Black 38-416-014

Burgundy 38-416-054

Gray 38-416-104

Navy Blue 38-416-184



Single Sprague Rappaport

Multi-functional dual-head that can be changed out for a pediatric diaphragm, adult diaphragm, pediatric bell, child bell, and adult bell.

Available in the following colors:

Black 38-415-014

Burgundy 38-415-054

Gray 38-415-104

Navy Blue 38-415-184



Professional Series

The Professional Series dual-head model consists of an adjustable aluminum binaural. A lightweight and comfortable design.

Available in the following colors:

Black 38-412-014

Burgundy 38-412-054

Gray 38-412-104

Navy Blue 38-412-184



3M™ Littmann® Classic II™ Pediatric

The 3M™ Littmann® Classic II Pediatric Stethoscope is specially designed for use on smaller patients.

Available in the following colors:

Black

3 M 2113

Raspberry

3 M 2122



3M™ Littmann® Lightweight II S.E.

The 3M™ Littmann® Lightweight II S.E. Stethoscope is an ideal entry-level instrument with a teardrop shape for easy auscultation.

Available in the following colors:

Black

3 M 2 4 5 0

Burgundy

3 M 2 4 5 1

Caribbean Blue

3 M 2 4 5 2



3M™ Littmann® Cardiology IV™ Diagnostic

Sturdy and refined, the Littmann Cardiology IV helps isolate subtle sounds so you can hear in the most critical and challenging environments.

Available in the following colors:

Black

3 M 6 1 5 2

Rose Pink

3 M 6 1 5 9

3M™ Littmann® Classic III™ Stethoscope

Tunable diaphragm technology and resilient next-generation tubing. Includes snap-tight, soft-sealing ear tips. Made in the USA.

Available in the following colors:

Black

3 M 5 6 2 0

Navy Blue

3 M 5 6 2 2

Burgundy

3 M 5 6 2 7

Available with a mirror finish:

Ceil Blue

3 M 5 9 5 9

Plum

3 M 5 9 6 0

Pearl Pink

3 M 5 9 6 2

Navy Blue

3 M 5 8 6 3



Visiofocus

The Visiofocus line is the only noncontact thermometer in the world to project the temperature reading. The Visiofocus is our mid-range model, Visiofocus Smart is our value model, and the Visiofocus Pro is our professional model. Each model has unique features for intended use.

Available in the following models:

Visiofocus

40-640-000-L

Visiofocus Smart

40-600-000-L

Visiofocus Pro

40-840-000-L







BV Medical® AnimaTemps Dog Digital Thermoemter

A kid-friendly digital thermometer which can be used for oral, rectal, or underarm use. This thermometer features a kid-friendly dog design.

40-208-000

Replacement Parts:

Batteries

41-100-010

Probe Covers

41-009-000



BV Medical® Instant Ear Thermometer

Our BV Medical Instant Ear Thermometer provides temperature readings instantly. Includes a stand, probe covers, and more.

40-500-000

Replacement Parts:

Probe Covers

41-500-100



BV Medical® Digital Thermometer

A 10-second rigid tip digital thermometer which can be used for oral, rectal, or underarm use. Durable thermometer with various features.

40-101-000-L

Replacement Parts:

Batteries

41-100-010

Probe Covers

41-009-000



BV Medical® Basic Pulse Oximeter

The BV Medical® Basic Pulse oximeter measures saturation of peripheral oxygen and pulse rate. This model is lightweight, compact, and simple to use.

50-100-004-L



BV Medical® Pediatric Pulse Oximeter

The BV Medical® Pediatric Pulse Oximeter accurately measures oxygen saturation level (SpO2) and pulse rate for children. For 2 years and older.

50-102-003-L



BV Medical® Classic Pulse Oximeter

The BV Medical® Classic Pulse Oximeter has a bright organic LED digital display making it easy to read. Displays saturation of peripheral oxygen, pulse rate, and perfusion index.

50-104-004-L



BV Medical® Optimum Pulse Oximeter

The BV Medical® Optimum Pulse Oximeter accurately measures oxygen saturation level (SpO2), pulse rate, and perfusion index (PI). This adult model features a large green display to make measurements easy-to-read.

50-200-010-L

OMRON®



OMRON CompAir®

The CompAir Compressor Nebulizer System is a small, lightweight and convenient tabletop nebulizer that features OMRON advanced design and engineering.

Prescription Required for End User

NE-C801



OMRON CompAir® for Kids

The COMP A-I-R mini compressor nebulizer for kids with VVT technology, small design, near silent operation, and kids accessory. A kid-friendly nebulizer.

Prescription Required for End User

RE-C800KD

OMRON Portable MicroAir®

It's time you break free from your respiratory restrictions. It's time to enjoy a long day at the beach, explore a new country or take a ski trip up north. That's because with the portable MicroAir® nebulizer, any place is the right place. It's a handheld nebulizer like no other.

Prescription Required for End User

NE-U100













BV Medical® Otoscope Penlight

The BV Medical Otoscope Penlight is valuable for ear, nose, and throat use by medical professionals. The penlight can be detached from the otoscope piece and be used separately. The pre-focused high intensity bulb allows for bright illumination. Easily attach to your pocket with the metal clip.

60-400-011

BV Medical® Mini Otoscope

The BV Medical Conventional Mini Otoscope is easy to use and carry. The 2.5 Volt vacuum bulb provides standard illumination. The ear funnel connection is via a screw thread for secure use. The BV Medical Mini Otoscope is available in different illumination options and colors for clinic or personal use.

Available in the colors Black, Royal Blue, Pink, and White

Conventional 2.5V

60-840-250-L

Also available in the following models:

F.O. Xenon 2.5V

60-830-250

F.O. Xenon/LED 2.5V

60-831-250-L





KAWE®



KaWe® EUROLIGHT

The high-quality metal housing of the KaWe EUROLIGHT product line makes the compact instruments sturdy, so you can have both excellent quality and high durability.

Available in the following model:

C10

01.11110.001



KaWe® COMBILIGHT

The KaWe COMBILIGHT otoscope optimally illuminates with an excellent cost to performance ratio. They are robust and operate reliably.

Available in the following model:

C10

01.12110.001



KaWe® PICCOLIGHT

The KaWe PICCOLIGHT F.O. LED Standard otoscope features a fiber optic (F.O.) LED illumination with a connector for a pneumatic test. This product is shockproof and energy-saving. The KaWe PICCOLIGHT F.O. LED Standard otoscope is backed by KaWe research and development which has been a globally active medical technology company for over 130 years.

Available in the following colors:

Night 01.13500.021

Sky 01.13500.232

Stone 01.13500.262

BV Medical® Chrome Penlight

The BV Medical Chrome Penlight comes with a pupil gauge ranging from 2 to 9 mm as well as a 5 cm ruler for your needs. The penlight has a LED bulb that illuminates brightly for you to see. With a metal pocket clip and easy push button on/off switch, you can bring your penlight anywhere. The chrome body is durable for continued use.

40-403-020





BV Medical® LED Penlight

The BV Medical penlight has a LED bulb that illuminates brightly (5 Lumens) for you to see. With a metal pocket clip and easy push button.

60-405-250

Available in the following colors:

Jet Black

60-405-010

Magenta

60-405-160

Cherry Red

60-405-240

Silver

60-405-260



BV Medical® Penlight

The BV Medical standard illumination penlight comes with a pupil gauge ranging from 2 to 9 mm. The penlight has a LED bulb.

60-402-000

Available in the following colors:

Pink/White

60-402-210

Available with no pupil gauge:

White

60-401-000



BV Medical® Standard Illumination Penlight

This penlight comes with a pupil gauge ranging from 2 to 9 mm, 5 cm ruler, LED bulb, metal pocket clip, and an easy push on/off button.

60-460-250

Available in the following colors:

Pink

60-460-210

Purple

60-460-230

Silver

60-460-260



BV Medical® Lister Bandage Scissors

The BV Medical® Lister Bandage Scissors are made of durable stainless steel. These pairs do not come with a clip, and the blades measure 5.5".

80-150-000



BV Medical® Dressing Scissors

The BV Medical® Dressing Scissors are made of durable stainless steel and have a sharp and blunt blade. These are straight scissors with blades measuring 5.5".

80-113-000



BV Medical® Lister Bandage Scissors - Clip

The BV Medical® Lister Bandage Scissors are made of durable stainless steel. These pairs do come with a clip, and the blades measure 5.5".

80-153-000



BV Medical® Utility Scissors

Durable stainless steel utility scissors with a plastic handle and plastic tip. These are suitable for both left or right handed use. These universal shears may be used by a Emergency Medical Technician (EMT), Paramedic, Nurse, or other medical professional. The plastic tip of the utility scissors are meant to minimize risk of injury to individual, and the rugged construction allows the ability to cut through strong materials such as car seat belts, leather, denim, etc.

80-180-SAM

BV Medical® Stainless Steel IV Pole


This IV Pole is made of durable stainless steel. You can quickly and easily adjust the height of the IV pole. IV pole measures 42 inches at the lowest height and 76 inches at its highest. The quad (4) ram hook can hold up to 44 lbs, and it has coverings at the end of the hooks for protection. The 5-legged polymer base provides stability and includes 5 wheels with locking casters if you need to make the IV Pole stay in place. Simple assembly is required and all tools are provided for assembly.

80-335-020









BV Medical® Sundry Jars

Each sundry jar is made of acrylic and clear for you to see the contents within it. The height of each jar is 6.5 inches (approximately 7 inches with lid), 4 inches in diameter, and a circumference of 13.25 inches. These medical sundry jars are able to be used to store medical necessities and reused.

Unlabeled Sundry Jars

80-275-000

Also available in the following model:

Labeled Sundry Jars

80-265-000





Hemorrhage Tourniquet

Hemorrhage tourniquets can quickly control lifethreatening extremity bleeding. This tourniquet allows for single-hand use. Includes adjustable sizing that extends to 39 inches at full length. Suitable for most hands and robust thighs.

Available in the following colors:

Black

80-700-000

Orange

80-700-190

Blue

80-700-250

Single Patient Tourniquets

This blue tourniquet may be used for single patient use and disposed of. Made of strong and durable elastic material for optimal performance. Tourniquets are also tear resistant which is essential for use on extremities. Measures 1" by 18".


Available in the following packs:

10 Pack

80-679-250-10

1.000 Pack

80-679-250-L



Conditions of Sale

Sales Policy

All prices are subject to change without notice.

Terms

Customers within Continental US: Terms are Net 30 days after credit approval, Prepaid, or Credit Card. An unpaid invoice is considered late if not received within 30 days from invoice date. If you are approved as a Credit Card customer then the invoice must be paid within 7 business days or else we will cancel the order. A service charge of 1½% per month will be charged on all overdue accounts. There is no minimum order requirement. For international customers, please contact Customer Service for terms.

Shipping

We reserve the right to use our own discretion in making shipments, unless specific shipping instructions are given on an order. Shipping charges for customers in the 50 states, free shipping may be offered (customers in a U.S. territory, please call for charges). Shipping charges on our website are subject to additional fees from chosen delivery services. BV Medical is not liable for shipping delays from carriers. Customers who wish to export BV devices assume all responsibility to meet destination country import requirements. Please contact our Customer Service Department regarding shipping at (888) 822-8293.

Claims

If you have received a damaged package or refuse shipments, immediately notify our Customer Service Department at (888) 822-8293, and we will notify the freight delivery company. Please note there is a 3 business day grace period to report any damaged merchandise from the time the shipment is received. After this time period, any damages become the responsibility of the customer. BV Medical cannot be responsible for any loss, damage, detention or delay caused by the elements, fire, strike, insurrection, riot, shortage of goods or transportation, civil or military authority, or by any other cause which is unavoidable or beyond our reasonable control.

Substitution Policy

BV Medical reserves the right to substitute equal or higher value/quality parts or products to fulfill an order, within reason, considering comparable functionality and appearance. Minor cosmetic differences sometimes occur on product packaging or in some products, compared to advertising literature, due to continuous improvement process.

Return Policy

BV Medical is committed to ensuring that every order you receive arrives to your satisfaction. If your product is the wrong size, fit, or damaged please call (888) 822-8293 and a Customer Service Representative will assist you with your questions and concerns. A Return Authorization number is required for all returns. Once an RA number has been issued, carefully repackage the item(s) you wish to return and ship to:

BV Medical Attention: Return Department/RA Number 28W206 Commercial Avenue Unit B Lake Barrington, IL 60010

Freight collect returns will not be accepted. All returns must be shipped prepaid by customer and clearly marked with the Return Authorization number issued by our Customer Service Department. Please retain all tracking information and a copy of your packing slip for your records. All authorized returns not received by BV Medical within 14 days after being issued will be voided. All authorized returns are subject to a twenty-five percent (25%) restocking charge. Returned goods must be in an unused, resalable condition, and may not be older than 30 days. Damaged merchandise or items requiring re-packaging will be billed at full cost. In addition, BV Medical reserves the right to refuse such shipments. Parts are not returnable for credit. BV Medical is not liable for shipping delays from carriers and are not subject to refund.

Limited Warranty

All BV Medical products which prove to be defective in material or workmanship within the warranty period will be repaired or replaced at our discretion when returned prepaid to:

BV Medical 28W206 Commercial Ave., Unit B Lake Barrington, IL 60010

We reserve the right to make the final determination as to whether the specific warranty is applicable to the returned item. If you have received a defective item from BV Medical, please call our Customer Service Department at (888) 822-8293 and a Customer Service Representative will assist you with the return process

Battery Powered Products

Where batteries are supplied with the original equipment, BV Medical does not provide any warranty as to the life that those batteries will provide. See the individual instruction manual for proper information as to the replacement of batteries.

Intended Use

Products sold by BV Medical are for their intended use only. BV Medical is not liable for any misuse of products.

Resal

By purchasing products from BV Medical, you also abide by our associated suppliers' conditions of sale. Omron®, Littman®, and Kinsa do not allow for marketplace competition meaning if you purchase these branded products from BV Medical then you are not allowed to sell on online marketplaces such as Amazon, Walmart, eBay, etc.

BV Medical® does not allow for online marketplace competition. If you purchase BV Medical product then you agree not to sell on marketplaces such as Amazon, Walmart, eBay, etc. without prior authorization. For questions or concerns regarding this policy then please contact us.

Fraud

All e-commerce payments are authorized before charged. As soon as you place an order online, your credit card/payment processor will not be charged until authorized. We have zero tolerance for credit card fraud and theft of our merchandise. Inconsistent/unusual information on your invoice details or transaction information will flag fraud detectors. As an extra step to protect our customers from internet fraud, we may require the customer to submit a government issued photo ID or valid license to authorize their transaction. Any fraudulent purchase attempts will be investigated and prosecuted to the fullest extent of the law.

Agreement to Conditions

Purchase of products sold by BV Medical implies customer agreement to Conditions of Sale. Updated Conditions of Sale can be found on website and are subject to change.

