

OMRON® Blood Pressure Comparison Chart



Model Number	BP7100	BP7200	BP7250	BP7350	BP7450	BP7000	BP7900	BP6100	BP6350	HEM-907XL
Type of Blood Pressure Monitor	Upper Arm	Upper Arm	Upper Arm	Upper Arm	Upper Arm	Upper Arm	Upper Arm	Wrist	Wrist	<p>This is a professional series upper arm blood pressure monitor. Request a Spec Sheet for details.</p>
Advanced Accuracy	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			
Hypertension Indicator		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Bluetooth®			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	
TruRead					<input checked="" type="checkbox"/>					
Extra-Large Digits with Backlight					<input checked="" type="checkbox"/>					
One Touch Use	<input checked="" type="checkbox"/>					<input checked="" type="checkbox"/>				
Two User Mode			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				
Body Movement Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				
Advanced Averaging		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Irregular Heartbeat Detector	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Memory Per User	14	60	60	60	100	100	90	60	90	
UltraSilent Inflation									<input checked="" type="checkbox"/>	
Heart Zone Guidance									<input checked="" type="checkbox"/>	
EKG							<input checked="" type="checkbox"/>			
Side-by-Side Comparison				<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>					
Cuff	Wide Range D-Ring 9"-17"	Wide Range D-Ring 9"-17"	Wide Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Wide Range D-Ring 9"-17"			
Compatible Cuffs	CD-WR17 and CD-CS9	HEM-RML31-B and HEM-CS24-B	HEM-RML31-B and HEM-CS24-B	HEM-FL31-B	HEM-FL31-B		HEM-RML31-B and HEM-CS24-B			
Storage Case							<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Power Options	4 x AA Batteries Not Included / Not Compatible with AC Adapter	4 x AA Batteries Included / AC Adapter Not Included	4 x AA Batteries Included / AC Adapter Not Included	4 x AA Batteries Included / AC Adapter Not Included	4 x AA Batteries Included / AC Adapter Not Included	4 x AAA Batteries Included / Not Compatible with AC Adapter	4 x AA Batteries Included / AC Adapter Not Included	2 x AAA Batteries Included	2 x AAA Batteries Included	