

Stair-Trac SA-S Wheelchair requirements

A: Wheelchair height

760mm to 950mm from the floor

B: Space under the X-frame

Center height: 210mm min.

End height: 170mm min. Crawler under-frame limits

C: Inside width between base tubes


320 mm or wider

D: Main Wheel diameter

22 inch or larger

The following Wheelchair types cannot be mounted

- Reclining type Wheelchair
- Foldable or detachable handle type Wheelchair
- Detachable armrest type Wheelchair
- 4. Caliper brake installed Wheelchair
- 5. Any Wheelchair that has machine part trouble



The dimensions shown above are standard values.
Some Wheelchairs may not fit.

