

# **GAM**

## GAM SERIES STAINLESS STEEL TECHNICAL PRECISION SCALE



IP65 high precision scales, suitable for laboratory and industrial use. The completely stainaless steel structure makes this scale suitable for the environments with strict hygiene requirements. Available also CE-M APPROVED for legal for trade use.

DIVISION (g): from 0.1 to 1. CAPACITY (g): from 320 to 15000.

## **TECHNICAL FEATURES**

- High contrast backlit display.
- Splashproof membrane keypad.
- Splashproof constuction, upper structure completely in stainless steel.
- Overall IP65 protection.
- Stainless steel plate, size:
  - GAM320 and GAM820 models Ø 140mm
  - GAM2200, GAM3200, GAM6200, GAM8200, GAM15000 models 190x190mm.
- Weight detection through tuning fork system.
- Level and adjustable foot assembly.
- Power supply: standard fitted with 230 Vac power adapter or optional kit for power supply through 4 AA batteries.
- RS232/c port for connection to PC, 3590E series indicator or printer.
- Upon request:
  - ACCREDIA calibration with issuing of relative certificate;
  - CE-M approval (Initial Verification Included, validity 3 years);

## STANDARD FUNCTIONS

- totalisation;
- · weight checking;
- counting;
- % weight;
- suspended weighing (for gravimetric measurements);
- animals weighing;
- zeroing;
- semi-automatic tare;
- digital calibration.

# ADVANCED APPLICATION WITH 3590E TOUCH SCREEN

- This configuration allows you to use the GAM scale as a digital platform for 3590E TOUCH SCREEN indicators.
- To create a weighing station with 3590E TOUCH SCREEN, you need to order together with the scale, the indicator, the RSCBSHG cable, and the S1 advanced configuration service.

#### ADVANCED FUNCTIONS via serial port connection:

- Statistical sampling of prepackaged goods
- Microdosing
- Multiscale counting
- Tolerance control
- Wizard formulation
- · Receipts, labels printing, storage and data transmission

# DETAIL 1







# VERSIONS

Versioni disponibili			
	Max	d	M
Codice	(g)	(g)	(g)*
GAM320	320	0,01	0,01
GAM820	820	0,01	0,01
GAM2200	2200	0,1	0,1
GAM3200	3200	0,1	0,1
GAM6200	6200	0,1	0,1
GAM8200	8200	0,1	0,1
GAM15000	15000	1	1

<sup>(\*)</sup> Divisioni CE-M ottenibili solo con le relative opzioni.



GRAVITATION LTD.

UNIT 2B SECTION 5 DUNSHAUGHLIN - - CO. MEATH - IRELAND

Tel. 00353469078900 Fax. 00353469078911 E-mail info@gravitation.ie Web www.gravitation.ie

