

★ Storage

Store at RT

☞ Contents

- Product Manual
- 20X TBS Buffer(Tris Buffered Saline), pH 7.6

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

★ Introduction

20X TBS (Tris Buffered Saline) solution consisting of 500mM Tris, 2.8 M NaCl in high purity RO water, pH 7.6. Solution is 0.2 μ m filtered and dispensed into durable, PET bottles. Dilute to 1X with high purity water for a solution consisting of 25mM Tris, pH 7.6, 140 mM NaCl. TBS buffer is used as wash buffer and diluent for ELISA, Western blotting and other immuno assays.

★ Formulation

Components	Conc (20X)
TRIS, pH 7.4	500 mM
NaCl	2.8 M