

## nanoComposix 2024 Standard Product Catalog Overview

Rev 3.0, 9 January 2024

### LATERAL FLOW ASSAY GOLD NANOPARTICLES & KITS

Type	Description
<b>BioReady Gold Nanoparticles for Lateral Flow:</b> <i>Nanospheres &amp; Nanoshells</i>	High quality monodisperse Colloidal Gold and high sensitivity Gold Nanoshells produced in bulk quantities. Citrate surface options available for passive adsorption and carboxyl surface options for covalent conjugation. See Gold Nanoparticles for details.
<b>Gold Conjugates:</b> <i>Nanospheres &amp; Nanoshells</i>	Streptavidin, Protein G, or Anti-Human IgG functionalized 40 nm gold nanospheres and 150 nm gold nanoshells. Covalently conjugated for high stability and provided at 10 OD. Available in bulk quantities.
<b>Covalent Conjugation Kit</b>	Choose 40 or 80 nm gold nanospheres or high sensitivity gold nanoshells. Includes EDC & Sulfo-NHS reagents plus reaction buffers and diluent for covalent conjugation to antibodies and some proteins.
<b>Passive Conjugation Kit</b>	Includes 40 nm and 80 nm colloidal gold plus a range of buffers and diluent for passive adsorption to proteins and antibodies
<b>Lateral Flow Material Starter Kit</b>	Quickly screen different raw materials to build and optimize your lateral flow assay. Includes a selection of nitrocellulose membranes, samples pads, conjugate pads, wick pad, and backing card.
<b>NHS-Dried Gold Kit</b>	Choose either 40 nm gold nanospheres or 150 nm gold nanoshells provided in a NHS-activated powder form. Includes antibody purification buffer and reaction buffer for quickly screening and evaluating conjugates.
<b>AmpliDetect NALF Kit</b>	The AmpliDetect assay is a universal tool for simple and quick detection of nucleic acid sequences and other molecules labeled with fluorescein (FITC/FAM) and biotin.
<b>BioReady Conjugation Buffers</b>	The same high quality covalent conjugation reaction buffers, block buffer, and conjugate diluent found in our High Sensitivity Gold Conjugation Kit, available in catalog and bulk volumes.

### GOLD NANOPARTICLES

Material	Types & Surfaces	Particle Sizes (nm)	Concentrations (OD or Au mass)
Lateral Flow Gold Reagents	Gold Nanospheres	Citrate (Passive Adsorption)	20 OD
		Carboxyl (Covalent Conjugation)	20 OD
	Gold Nanoshells	Citrate (Passive Adsorption)	20 OD
		Carboxyl (Covalent Conjugation)	20 OD
Magnetic Gold Nanoshells	Carboxyl (Covalent Conjugation)	200 nm	10 OD
Gold Conjugates (Nanospheres & Nanoshells)	Streptavidin	40, 150	10 OD
	Protein G	40	10 OD
	Anti-Human IgG (H+L)	40	10 OD
Gold Nanospheres	Citrate, PEG, Carboxyl, PVP	5 - 200	1 mg/mL (Low Endotoxin) 0.05 mg/mL (Citrate Only)
Ultra Uniform Gold Nanospheres	Ultra-low CV (<5%)	10, 20, 30, 50, 100	0.05 mg/mL
Ultra Uniform Gold Nanoparticle Number Standard	Ultra-low CV (<6%)	30, 50, 100	1 x 10 <sup>8</sup> particles/mL
Gold Nanorods	Citrate, PEG, Carboxyl	660, 800, 980 nm resonant	1 OD, 50 OD
Silica-Shelled Gold Nanospheres	Silanol Surface, Aminated Surface	60, 90, 140	1 mg/mL
Dodecanethiol-Stabilized Gold Nanospheres	Organic Solvent Compatible	2, 5	Dried
Polystyrene-Coated Gold Nanospheres	Organic Solvent Compatible	10, 50, 100	1 mg/mL
Gold Nanoshells	PEG, PVP, Carboxyl	660, 800, & 980 nm resonant	1 mg/mL, 0.05 mg/mL

### SILVER NANOPARTICLES

Material	Types & Surfaces	Particle Sizes (nm)	Concentrations (Ag mass)
Silver Nanospheres	Citrate, PEG, Carboxyl, PVP	5 - 200	1 mg/mL (Low Endotoxin) 0.02 mg/mL (Citrate Only)
Silver Nanocubes	Ultra-low CV (<5%)	10, 20, 30, 50, 100	1 mg/mL
Silver Nanoplates	Silanol Surface, Aminated Surface	60, 90, 140	1 mg/mL
Silica-Shelled Silver Nanospheres	Organic Solvent Compatible	2, 5	Dried
Dodecanethiol-Stabilized Silver Nanospheres	Organic Solvent Compatible	10, 50, 100	1 mg/mL

### SILICA NANOPARTICLES

Material	Types & Surfaces	Particle Sizes (nm)	Concentrations (SiO <sub>2</sub> mass)
Silica Nanospheres	Silanol Surface, Aminated Surface	20-1000	10 mg/mL or dried (5 mg/mL for 20 nm size)
Mesoporous Silica Nanospheres	Silanol Surface, Aminated Surface	100	5 mg/mL or dried

### PLATINUM & MAGNETITE NANOPARTICLES

Material	Types & Surfaces	Particle Sizes (nm)	Concentrations (metal mass)
Platinum Nanospheres	Citrate	5, 30, 50, 70	1 mg/mL (Low Endotoxin) 0.05 mg/mL (Citrate Only)
Magnetite Nanospheres	PVP Stabilized	20	1 mg/mL (Low Endotoxin) 20 mg/mL, dried Dried

Catalog Number (SKU)	Full Item Description	2024 List Price
<b>Gold Nanospheres for Conjugation</b>		
<b>Bare Gold Nanospheres for Passive Adsorption, BioReady, 20 OD</b>		
AUCR40-5M	BioReady Gold Nanospheres - Bare (Citrate) - 40 nm, 20 OD in aq. 0.02 mM sodium citrate, 5 mL	\$80.00
AUCR40-30M	BioReady Gold Nanospheres - Bare (Citrate) - 40 nm, 20 OD in aq. 0.02 mM sodium citrate, 30 mL	\$425.00
AUCR40-100M	BioReady Gold Nanospheres - Bare (Citrate) - 40 nm, 20 OD in aq. 0.02 mM sodium citrate, 100 mL	\$1,295.00
AUCR80-5M	BioReady Gold Nanospheres - Bare (Citrate) - 80 nm, 20 OD in aq. 0.02 mM sodium citrate, 5 mL	\$90.00
AUCR80-30M	BioReady Gold Nanospheres - Bare (Citrate) - 80 nm, 20 OD in aq. 0.02 mM sodium citrate, 30 mL	\$490.00
AUCR80-100M	BioReady Gold Nanospheres - Bare (Citrate) - 80 nm, 20 OD in aq. 0.02 mM sodium citrate, 100 mL	\$1,490.00
<b>Gold Nanospheres for Covalent Conjugation, BioReady, 20 OD</b>		
AUXR40-5M	BioReady Gold Nanospheres - Carboxyl - 40 nm, 20 OD in water, 5 mL	\$150.00
AUXR40-30M	BioReady Gold Nanospheres - Carboxyl - 40 nm, 20 OD in water, 30 mL	\$645.00
AUXR40-100M	BioReady Gold Nanospheres - Carboxyl - 40 nm, 20 OD in water, 100 mL	\$1,800.00
AUXR80-5M	BioReady Gold Nanospheres - Carboxyl - 80 nm, 20 OD in water, 5 mL	\$150.00
AUXR80-30M	BioReady Gold Nanospheres - Carboxyl - 80 nm, 20 OD in water, 30 mL	\$645.00
AUXR80-100M	BioReady Gold Nanospheres - Carboxyl - 80 nm, 20 OD in water, 100 mL	\$1,800.00
<b>BioReady Covalent Conjugation Kit for Lateral Flow, Gold Nanospheres, 20 OD</b>		
AUZR40-10M	BioReady Covalent Gold Conjugation Kit for Lateral Flow - 40 nm Gold Nanospheres	\$1,130.00
AUZR80-10M	BioReady Covalent Gold Conjugation Kit for Lateral Flow - 80 nm Gold Nanospheres	\$1,130.00
<b>Gold Nanospheres for Covalent Conjugation, BioReady, Dried (Conjugation Kit)</b>		
AUNR40-3SRF	BioReady Gold Nanospheres - NHS (Dried Kit) - 40 nm, Small (3 Units)	\$260.00
AUNR40-10SRF	BioReady Gold Nanospheres - NHS (Dried Kit) - 40 nm, Small (10 Units)	\$635.00
AUNR40-1LRF	BioReady Gold Nanospheres - NHS (Dried Kit) - 40 nm, Large (1 Unit)	\$635.00
<b>BioReady Passive Conjugation Kit for Lateral Flow, Gold Nanospheres, 20 OD</b>		
PCKR	BioReady Passive Gold Conjugation Kit - 40 nm, and 80 nm, Citrate Gold Nanospheres	\$800.00
<b>Gold Nanoshells for Conjugation</b>		
<b>Gold Nanoshells for Conjugation, BioReady, 20 OD</b>		
GSXR150-5M	BioReady Gold Nanoshells - Carboxyl - 150 nm, 20 OD in aqueous 0.02 mM potassium carbonate, 5 mL	\$185.00
GSXR150-30M	BioReady Gold Nanoshells - Carboxyl - 150 nm, 20 OD in aqueous 0.02 mM potassium carbonate, 30 mL	\$755.00
GSXR150-100M	BioReady Gold Nanoshells - Carboxyl - 150 nm, 20 OD in aqueous 0.02 mM potassium carbonate, 100 mL	\$2,130.00
GSCRI150-5M	BioReady Gold Nanoshells - Citrate - 150 nm, 20 OD in aqueous 0.2 mM sodium citrate, 5 mL	\$175.00
GSCRI150-30M	BioReady Gold Nanoshells - Citrate - 150 nm, 20 OD in aqueous 0.2 mM sodium citrate, 30 mL	\$720.00
GSCRI150-100M	BioReady Gold Nanoshells - Citrate - 150 nm, 20 OD in aqueous 0.2 mM sodium citrate, 100 mL	\$2,030.00
<b>BioReady High Sensitivity Conjugation Kit for Lateral Flow, Gold Nanoshells, 20 OD</b>		
GSZR150-10M	BioReady High Sensitivity Gold Conjugation Kit for Lateral Flow - 150 nm Gold Nanoshells	\$1,130.00
<b>Gold Nanoshells for Conjugation, BioReady, Dried</b>		
GSNR150-3SRF	BioReady Gold Nanoshells - NHS (Dried Kit) - 150 nm, Small (3 Units)	\$290.00
GSNR150-10SRF	BioReady Gold Nanoshells - NHS (Dried Kit) - 150 nm, Small (10 Units)	\$710.00
GSNR150-1LRF	BioReady Gold Nanoshells - NHS (Dried Kit) - 150 nm, Large (1 Unit)	\$710.00
<b>Magnetic Gold Nanoshells</b>		
<b>Magnetic Gold Nanoshells, BioReady, 10 OD</b>		
MSXR200-1M	BioReady Magnetic Gold Nanoshells - Carboxyl - 200 nm, 10 OD in aqueous 0.02 mM potassium carbonate, 1 mL	\$165.00
MSXR200-5M	BioReady Magnetic Gold Nanoshells - Carboxyl - 200 nm, 10 OD in aqueous 0.02 mM potassium carbonate, 5 mL	\$295.00
MSXR200-10M	BioReady Magnetic Gold Nanoshells - Carboxyl - 200 nm, 10 OD in aqueous 0.02 mM potassium carbonate, 10 mL	\$465.00
MSXR200-50M	BioReady Magnetic Gold Nanoshells - Carboxyl - 200 nm, 10 OD in aqueous 0.02 mM potassium carbonate, 50 mL	\$2,095.00
<b>Gold Conjugates (Nanospheres &amp; Nanoshells)</b>		
<b>Gold Nanospheres with Streptavidin, BioReady, 10 OD</b>		
AUIR40-200U	BioReady Gold Nanospheres - Streptavidin - 40 nm, 10 OD in aqueous conjugation buffer, 200 µL	\$150.00
AUIR40-1M	BioReady Gold Nanospheres - Streptavidin - 40 nm, 10 OD in aqueous conjugation buffer, 1 mL	\$475.00
AUIR40-5M	BioReady Gold Nanospheres - Streptavidin - 40 nm, 10 OD in aqueous conjugation buffer, 5 mL	\$1,905.00
AUIR40-10M	BioReady Gold Nanospheres - Streptavidin - 40 nm, 10 OD in aqueous conjugation buffer, 10 mL	\$3,215.00
<b>Gold Nanospheres with Protein G, BioReady, 10 OD</b>		
AUPGR40-200U	BioReady Gold Nanospheres - Protein G - 40 nm, 10 OD in aqueous conjugation buffer, 200 µL	\$150.00
AUPGR40-1M	BioReady Gold Nanospheres - Protein G - 40 nm, 10 OD in aqueous conjugation buffer, 1 mL	\$400.00
AUPGR40-5M	BioReady Gold Nanospheres - Protein G - 40 nm, 10 OD in aqueous conjugation buffer, 5 mL	\$1,750.00
AUPGR40-10M	BioReady Gold Nanospheres - Protein G - 40 nm, 10 OD in aqueous conjugation buffer, 10 mL	\$3,200.00
<b>Gold Nanospheres with Anti-Human IgG (H+L), BioReady, 10 OD</b>		
AUAHR40-200U	BioReady Gold Nanospheres - Anti-Human IgG (H+L) - 40 nm, 10 OD in aq. conjugation buffer, 200 µL	\$150.00
AUAHR40-1M	BioReady Gold Nanospheres - Anti-Human IgG (H+L) - 40 nm, 10 OD in aq. conjugation buffer, 1 mL	\$400.00
AUAHR40-5M	BioReady Gold Nanospheres - Anti-Human IgG (H+L) - 40 nm, 10 OD in aq. conjugation buffer, 5 mL	\$1,750.00

UAHR40-10M	BioReady Gold Nanospheres - Anti-Human IgG (H+L) - 40 nm, 10 OD in aq. conjugation buffer, 10 mL	\$3,200.00
------------	--	------------

Gold Nanoshells with Streptavidin, BioReady, 10 OD		
GSIR150-200U	BioReady Gold Nanoshells - Streptavidin - 150 nm, 10 OD in aqueous conjugation buffer, 200 µL	\$170.00
GSIR150-1M	BioReady Gold Nanoshells - Streptavidin - 150 nm, 10 OD in aqueous conjugation buffer, 1 mL	\$510.00
GSIR150-5M	BioReady Gold Nanoshells - Streptavidin - 150 nm, 10 OD in aqueous conjugation buffer, 5 mL	\$2,130.00
GSIR150-10M	BioReady Gold Nanoshells - Streptavidin - 150 nm, 10 OD in aqueous conjugation buffer, 10 mL	\$3,445.00

Other Gold Nanospheres		
Gold Nanospheres, BioPure, 1 mg/mL		
AUCB5-1M	BioPure Gold Nanospheres - Bare (Citrate) - 5 nm, 1 mg/mL in water, 1 mL	\$285.00
AUCB5-5M	BioPure Gold Nanospheres - Bare (Citrate) - 5 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUCB5-10M	BioPure Gold Nanospheres - Bare (Citrate) - 5 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUCB5-30M	BioPure Gold Nanospheres - Bare (Citrate) - 5 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUCB10-1M	BioPure Gold Nanospheres - Bare (Citrate) - 10 nm, 1 mg/mL in water, 1 mL	\$285.00
AUCB10-5M	BioPure Gold Nanospheres - Bare (Citrate) - 10 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUCB10-10M	BioPure Gold Nanospheres - Bare (Citrate) - 10 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUCB10-30M	BioPure Gold Nanospheres - Bare (Citrate) - 10 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUCB20-1M	BioPure Gold Nanospheres - Bare (Citrate) - 20 nm, 1 mg/mL in water, 1 mL	\$285.00
AUCB20-5M	BioPure Gold Nanospheres - Bare (Citrate) - 20 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUCB20-10M	BioPure Gold Nanospheres - Bare (Citrate) - 20 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUCB20-30M	BioPure Gold Nanospheres - Bare (Citrate) - 20 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUCB30-1M	BioPure Gold Nanospheres - Bare (Citrate) - 30 nm, 1 mg/mL in water, 1 mL	\$285.00
AUCB30-5M	BioPure Gold Nanospheres - Bare (Citrate) - 30 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUCB30-10M	BioPure Gold Nanospheres - Bare (Citrate) - 30 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUCB30-30M	BioPure Gold Nanospheres - Bare (Citrate) - 30 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUCB40-1M	BioPure Gold Nanospheres - Bare (Citrate) - 40 nm, 1 mg/mL in water, 1 mL	\$285.00
AUCB40-5M	BioPure Gold Nanospheres - Bare (Citrate) - 40 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUCB40-10M	BioPure Gold Nanospheres - Bare (Citrate) - 40 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUCB40-30M	BioPure Gold Nanospheres - Bare (Citrate) - 40 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUCB50-1M	BioPure Gold Nanospheres - Bare (Citrate) - 50 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 1 mL	\$285.00
AUCB50-5M	BioPure Gold Nanospheres - Bare (Citrate) - 50 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 5 mL	\$1,130.00
AUCB50-10M	BioPure Gold Nanospheres - Bare (Citrate) - 50 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 10 mL	\$1,865.00
AUCB50-30M	BioPure Gold Nanospheres - Bare (Citrate) - 50 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 30 mL	\$3,740.00
AUCB60-1M	BioPure Gold Nanospheres - Bare (Citrate) - 60 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 1 mL	\$285.00
AUCB60-5M	BioPure Gold Nanospheres - Bare (Citrate) - 60 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 5 mL	\$1,130.00
AUCB60-10M	BioPure Gold Nanospheres - Bare (Citrate) - 60 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 10 mL	\$1,865.00
AUCB60-30M	BioPure Gold Nanospheres - Bare (Citrate) - 60 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 30 mL	\$3,740.00
AUCB80-1M	BioPure Gold Nanospheres - Bare (Citrate) - 80 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 1 mL	\$285.00
AUCB80-5M	BioPure Gold Nanospheres - Bare (Citrate) - 80 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 5 mL	\$1,130.00
AUCB80-10M	BioPure Gold Nanospheres - Bare (Citrate) - 80 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 10 mL	\$1,865.00
AUCB80-30M	BioPure Gold Nanospheres - Bare (Citrate) - 80 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 30 mL	\$3,740.00
AUCB100-1M	BioPure Gold Nanospheres - Bare (Citrate) - 100 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 1 mL	\$285.00
AUCB100-5M	BioPure Gold Nanospheres - Bare (Citrate) - 100 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 5 mL	\$1,130.00
AUCB100-10M	BioPure Gold Nanospheres - Bare (Citrate) - 100 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 10 mL	\$1,865.00
AUCB100-30M	BioPure Gold Nanospheres - Bare (Citrate) - 100 nm, 1 mg/mL in aqueous 0.2 mM sodium citrate, 30 mL	\$3,740.00
AUPB10-1M	BioPure Gold Nanospheres - PVP - 10 nm, 1 mg/mL in water, 1 mL	\$285.00
AUPB10-5M	BioPure Gold Nanospheres - PVP - 10 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUPB10-10M	BioPure Gold Nanospheres - PVP - 10 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUPB10-30M	BioPure Gold Nanospheres - PVP - 10 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUPB50-1M	BioPure Gold Nanospheres - PVP - 50 nm, 1 mg/mL in water, 1 mL	\$285.00
AUPB50-5M	BioPure Gold Nanospheres - PVP - 50 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUPB50-10M	BioPure Gold Nanospheres - PVP - 50 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUPB50-30M	BioPure Gold Nanospheres - PVP - 50 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUPB100-1M	BioPure Gold Nanospheres - PVP - 100 nm, 1 mg/mL in water, 1 mL	\$285.00
AUPB100-5M	BioPure Gold Nanospheres - PVP - 100 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AUPB100-10M	BioPure Gold Nanospheres - PVP - 100 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AUPB100-30M	BioPure Gold Nanospheres - PVP - 100 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AULB10-1M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 10 nm, 1 mg/mL in water, 1 mL	\$285.00
AULB10-5M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 10 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AULB10-10M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 10 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AULB10-30M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 10 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AULB50-1M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 50 nm, 1 mg/mL in water, 1 mL	\$285.00
AULB50-5M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 50 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AULB50-10M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 50 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AULB50-30M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 50 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AULB100-1M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 100 nm, 1 mg/mL in water, 1 mL	\$285.00
AULB100-5M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 100 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AULB100-10M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 100 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AULB100-30M	BioPure Gold Nanospheres - Carboxyl (Lipoic Acid)- 100 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AUGB10-1M	BioPure Gold Nanospheres - PEG - 10 nm, 1 mg/mL in water, 1 mL	\$285.00



AUPD100-10MG	NanoXact Gold Nanospheres - PVP (Dried) - 100 nm, 10 mg (by gold mass)	\$800.00
AUPD100-30MG	NanoXact Gold Nanospheres - PVP (Dried) - 100 nm, 30 mg (by gold mass)	\$1,745.00

## Ultra Uniform Gold Nanospheres

### Gold Nanospheres, Ultra Uniform, 0.05 mg/mL

AUXU10-1M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 10 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$125.00
AUXU10-5M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 10 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$445.00
AUXU10-10M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 10 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$720.00
AUXU10-30M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 10 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$1,565.00
AUXU20-1M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 20 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$125.00
AUXU20-5M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 20 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$445.00
AUXU20-10M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 20 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$720.00
AUXU20-30M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 20 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$1,565.00
AUXU30-1M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 30 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$125.00
AUXU30-5M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 30 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$445.00
AUXU30-10M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 30 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$720.00
AUXU30-30M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 30 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$1,565.00
AUXU50-1M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 50 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$125.00
AUXU50-5M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 50 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$445.00
AUXU50-10M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 50 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$720.00
AUXU50-30M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 50 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$1,565.00
AUXU100-1M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 100 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$125.00
AUXU100-5M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 100 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$445.00
AUXU100-10M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 100 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$720.00
AUXU100-30M	Ultra Uniform Gold Nanospheres - PEG-Carboxyl - 100 nm, 0.05 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$1,565.00

### Ultra Uniform Gold Nanoparticle Number Standard, 1 x 10<sup>8</sup> particles/mL

AUXS30-25M	Ultra Uniform Gold Nanoparticle Number Standard - 30 nm, 1E8 particles/mL, 25 mL	\$835.00
AUXS50-25M	Ultra Uniform Gold Nanoparticle Number Standard - 50 nm, 1E8 particles/mL, 25 mL	\$835.00
AUXS100-25M	Ultra Uniform Gold Nanoparticle Number Standard - 100 nm, 1E8 particles/mL, 25 mL	\$835.00

## Silica-Shelled Gold Nanospheres

### Silica-Shelled Gold Nanospheres, NanoXact, 1 mg/mL

AUSH20-1M	NanoXact Gold Nanospheres - Silica Shelled - 60 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AUSH20-5M	NanoXact Gold Nanospheres - Silica Shelled - 60 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AUSH20-10M	NanoXact Gold Nanospheres - Silica Shelled - 60 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AUSH20-30M	NanoXact Gold Nanospheres - Silica Shelled - 60 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AUSH50-1M	NanoXact Gold Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AUSH50-5M	NanoXact Gold Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AUSH50-10M	NanoXact Gold Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AUSH50-30M	NanoXact Gold Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AUSH100-1M	NanoXact Gold Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AUSH100-5M	NanoXact Gold Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AUSH100-10M	NanoXact Gold Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AUSH100-30M	NanoXact Gold Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AUAH20-1M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 60 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AUAH20-5M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 60 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AUAH20-10M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 60 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AUAH20-30M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 60 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AUAH50-1M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AUAH50-5M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AUAH50-10M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AUAH50-30M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AUAH100-1M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AUAH100-5M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AUAH100-10M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AUAH100-30M	NanoXact Gold Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00

## Dodecanethiol-Stabilized Gold Nanospheres

### Dodecanethiol-Stabilized Gold Nanospheres, NanoXact, Dried

AUDD2-1MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 2 nm, 1 mg	\$220.00
AUDD2-5MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 2 nm, 5 mg	\$870.00
AUDD2-25MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 2 nm, 25 mg	\$2,555.00
AUDD2-50MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 2 nm, 50 mg	\$4,070.00
AUDD5-1MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 5 nm, 1 mg	\$220.00
AUDD5-5MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 5 nm, 5 mg	\$870.00
AUDD5-25MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 5 nm, 25 mg	\$2,555.00
AUDD5-50MG	NanoXact Gold Nanospheres - Dodecanethiol (Dried) - 5 nm, 50 mg	\$4,070.00

## Polystyrene-Coated Gold Nanospheres

Polystyrene-Coated Gold Nanospheres, NanoXact, 1 mg/mL

AUYH10-1M	NanoXact Gold Nanospheres - Polystyrene - 10 nm, 1 mg/mL, 1 mL	\$430.00
AUYH10-5M	NanoXact Gold Nanospheres - Polystyrene - 10 nm, 1 mg/mL, 5 mL	\$1,640.00
AUYH10-10M	NanoXact Gold Nanospheres - Polystyrene - 10 nm, 1 mg/mL, 10 mL	\$2,650.00
AUYH10-30M	NanoXact Gold Nanospheres - Polystyrene - 10 nm, 1 mg/mL, 30 mL	\$5,625.00
AUYH50-1M	NanoXact Gold Nanospheres - Polystyrene - 50 nm, 1 mg/mL, 1 mL	\$430.00
AUYH50-5M	NanoXact Gold Nanospheres - Polystyrene - 50 nm, 1 mg/mL, 5 mL	\$1,640.00
AUYH50-10M	NanoXact Gold Nanospheres - Polystyrene - 50 nm, 1 mg/mL, 10 mL	\$2,650.00
AUYH50-30M	NanoXact Gold Nanospheres - Polystyrene - 50 nm, 1 mg/mL, 30 mL	\$5,625.00
AUYH100-1M	NanoXact Gold Nanospheres - Polystyrene - 100 nm, 1 mg/mL, 1 mL	\$430.00
AUYH100-5M	NanoXact Gold Nanospheres - Polystyrene - 100 nm, 1 mg/mL, 5 mL	\$1,640.00
AUYH100-10M	NanoXact Gold Nanospheres - Polystyrene - 100 nm, 1 mg/mL, 10 mL	\$2,650.00
AUYH100-30M	NanoXact Gold Nanospheres - Polystyrene - 100 nm, 1 mg/mL, 30 mL	\$5,625.00

## Silver-Shelled Gold Nanospheres

Silver-Shelled Gold Nanospheres, NanoXact, 1 mg/mL

BMCH20-1M	NanoXact Silver Shelled Gold - Bare (Citrate) - 20 nm, 1 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$285.00
BMCH20-5M	NanoXact Silver Shelled Gold - Bare (Citrate) - 20 nm, 1 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$1,130.00
BMCH20-10M	NanoXact Silver Shelled Gold - Bare (Citrate) - 20 nm, 1 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$1,865.00
BMCH20-30M	NanoXact Silver Shelled Gold - Bare (Citrate) - 20 nm, 1 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$3,740.00
BMCH60-1M	NanoXact Silver Shelled Gold - Bare (Citrate) - 60 nm, 1 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$285.00
BMCH60-5M	NanoXact Silver Shelled Gold - Bare (Citrate) - 60 nm, 1 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$1,130.00
BMCH60-10M	NanoXact Silver Shelled Gold - Bare (Citrate) - 60 nm, 1 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$1,865.00
BMCH60-30M	NanoXact Silver Shelled Gold - Bare (Citrate) - 60 nm, 1 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$3,740.00
BMCH80-1M	NanoXact Silver Shelled Gold - Bare (Citrate) - 80 nm, 1 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$285.00
BMCH80-5M	NanoXact Silver Shelled Gold - Bare (Citrate) - 80 nm, 1 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$1,130.00
BMCH80-10M	NanoXact Silver Shelled Gold - Bare (Citrate) - 80 nm, 1 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$1,865.00
BMCH80-30M	NanoXact Silver Shelled Gold - Bare (Citrate) - 80 nm, 1 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$3,740.00

## Other Gold Nanoshells

Gold Nanoshells, NanoXact, 1 mg/mL

GSPH660-1M	NanoXact Gold Nanoshells - PVP - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 1 mL	\$285.00
GSPH660-5M	NanoXact Gold Nanoshells - PVP - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 5 mL	\$1,130.00
GSPH660-10M	NanoXact Gold Nanoshells - PVP - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 10 mL	\$1,805.00
GSPH660-30M	NanoXact Gold Nanoshells - PVP - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 30 mL	\$3,615.00
GSPH800-1M	NanoXact Gold Nanoshells - PVP - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 1 mL	\$285.00
GSPH800-5M	NanoXact Gold Nanoshells - PVP - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 5 mL	\$1,130.00
GSPH800-10M	NanoXact Gold Nanoshells - PVP - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 10 mL	\$1,805.00
GSPH800-30M	NanoXact Gold Nanoshells - PVP - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 30 mL	\$3,615.00
GSPH980-1M	NanoXact Gold Nanoshells - PVP - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 1 mL	\$285.00
GSPH980-5M	NanoXact Gold Nanoshells - PVP - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 5 mL	\$1,130.00
GSPH980-10M	NanoXact Gold Nanoshells - PVP - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 10 mL	\$1,805.00
GSPH980-30M	NanoXact Gold Nanoshells - PVP - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 30 mL	\$3,615.00
GSLH660-1M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 1 mL	\$285.00
GSLH660-5M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 5 mL	\$1,130.00
GSLH660-10M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 10 mL	\$1,805.00
GSLH660-30M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 30 mL	\$3,615.00
GSLH800-1M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 1 mL	\$285.00
GSLH800-5M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 5 mL	\$1,130.00
GSLH800-10M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 10 mL	\$1,805.00
GSLH800-30M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 30 mL	\$3,615.00
GSLH980-1M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 1 mL	\$285.00
GSLH980-5M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 5 mL	\$1,130.00
GSLH980-10M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 10 mL	\$1,805.00
GSLH980-30M	NanoXact Gold Nanoshells - Carboxyl - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in 0.02 mM potassium carbonate, 30 mL	\$3,615.00
GSGH660-1M	NanoXact Gold Nanoshells - PEG - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 1 mL	\$285.00
GSGH660-5M	NanoXact Gold Nanoshells - PEG - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 5 mL	\$1,130.00
GSGH660-10M	NanoXact Gold Nanoshells - PEG - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 10 mL	\$1,805.00
GSGH660-30M	NanoXact Gold Nanoshells - PEG - Peak @ 660 nm · 120 nm diameter, 1 mg/mL in water, 30 mL	\$3,615.00
GSGH800-1M	NanoXact Gold Nanoshells - PEG - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 1 mL	\$285.00
GSGH800-5M	NanoXact Gold Nanoshells - PEG - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 5 mL	\$1,130.00
GSGH800-10M	NanoXact Gold Nanoshells - PEG - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 10 mL	\$1,805.00
GSGH800-30M	NanoXact Gold Nanoshells - PEG - Peak @ 800 nm · 160 nm diameter, 1 mg/mL in water, 30 mL	\$3,615.00
GSGH980-1M	NanoXact Gold Nanoshells - PEG - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 1 mL	\$285.00
GSGH980-5M	NanoXact Gold Nanoshells - PEG - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 5 mL	\$1,130.00
GSGH980-10M	NanoXact Gold Nanoshells - PEG - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 10 mL	\$1,805.00
GSGH980-30M	NanoXact Gold Nanoshells - PEG - Peak @ 980 nm · 240 nm diameter, 1 mg/mL in water, 30 mL	\$3,615.00

Gold Nanoshells, NanoXact, 0.05 mg/mL







AGPB100-30M	BioPure Silver Nanospheres - PVP - 100 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AGLB50-1M	BioPure Silver Nanospheres - Carboxyl - 50 nm, 1 mg/mL in 0.02 mM potassium carbonate, 1 mL	\$285.00
AGLB50-5M	BioPure Silver Nanospheres - Carboxyl - 50 nm, 1 mg/mL in 0.02 mM potassium carbonate, 5 mL	\$1,130.00
AGLB50-10M	BioPure Silver Nanospheres - Carboxyl - 50 nm, 1 mg/mL in 0.02 mM potassium carbonate, 10 mL	\$1,865.00
AGLB50-30M	BioPure Silver Nanospheres - Carboxyl - 50 nm, 1 mg/mL in 0.02 mM potassium carbonate, 30 mL	\$3,740.00
AGLB100-1M	BioPure Silver Nanospheres - Carboxyl - 100 nm, 1 mg/mL in 0.02 mM potassium carbonate, 1 mL	\$285.00
AGLB100-5M	BioPure Silver Nanospheres - Carboxyl - 100 nm, 1 mg/mL in 0.02 mM potassium carbonate, 5 mL	\$1,130.00
AGLB100-10M	BioPure Silver Nanospheres - Carboxyl - 100 nm, 1 mg/mL in 0.02 mM potassium carbonate, 10 mL	\$1,865.00
AGLB100-30M	BioPure Silver Nanospheres - Carboxyl - 100 nm, 1 mg/mL in 0.02 mM potassium carbonate, 30 mL	\$3,740.00
AGGB50-1M	BioPure Silver Nanospheres - PEG - 50 nm, 1 mg/mL in water, 1 mL	\$285.00
AGGB50-5M	BioPure Silver Nanospheres - PEG - 50 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AGGB50-10M	BioPure Silver Nanospheres - PEG - 50 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AGGB50-30M	BioPure Silver Nanospheres - PEG - 50 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AGGB100-1M	BioPure Silver Nanospheres - PEG - 100 nm, 1 mg/mL in water, 1 mL	\$285.00
AGGB100-5M	BioPure Silver Nanospheres - PEG - 100 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AGGB100-10M	BioPure Silver Nanospheres - PEG - 100 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AGGB100-30M	BioPure Silver Nanospheres - PEG - 100 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AGBB50-1M	BioPure Silver Nanospheres - BPEI - 50 nm, 1 mg/mL in water, 1 mL	\$285.00
AGBB50-5M	BioPure Silver Nanospheres - BPEI - 50 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AGBB50-10M	BioPure Silver Nanospheres - BPEI - 50 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AGBB50-30M	BioPure Silver Nanospheres - BPEI - 50 nm, 1 mg/mL in water, 30 mL	\$3,740.00
AGBB100-1M	BioPure Silver Nanospheres - BPEI - 100 nm, 1 mg/mL in water, 1 mL	\$285.00
AGBB100-5M	BioPure Silver Nanospheres - BPEI - 100 nm, 1 mg/mL in water, 5 mL	\$1,130.00
AGBB100-10M	BioPure Silver Nanospheres - BPEI - 100 nm, 1 mg/mL in water, 10 mL	\$1,865.00
AGBB100-30M	BioPure Silver Nanospheres - BPEI - 100 nm, 1 mg/mL in water, 30 mL	\$3,740.00
<b>Silver Nanospheres, NanoXact, 0.02 mg/mL</b>		
AGCN10-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 10 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN10-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 10 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN10-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 10 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN10-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 10 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN20-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 20 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN20-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 20 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN20-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 20 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN20-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 20 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN30-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 30 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN30-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 30 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN30-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 30 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN30-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 30 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN40-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 40 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN40-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 40 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN40-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 40 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN40-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 40 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN50-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 50 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN50-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 50 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN50-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 50 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN50-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 50 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN60-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 60 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN60-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 60 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN60-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 60 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN60-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 60 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN80-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 80 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN80-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 80 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN80-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 80 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN80-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 80 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN100-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 100 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN100-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 100 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN100-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 100 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN100-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 100 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
AGCN200-50M	NanoXact Silver Nanospheres - Bare (Citrate) - 200 nm, 0.02 mg/mL in 2 mM sodium citrate, 50 mL	\$125.00
AGCN200-100M	NanoXact Silver Nanospheres - Bare (Citrate) - 200 nm, 0.02 mg/mL in 2 mM sodium citrate, 100 mL	\$190.00
AGCN200-500M	NanoXact Silver Nanospheres - Bare (Citrate) - 200 nm, 0.02 mg/mL in 2 mM sodium citrate, 500 mL	\$720.00
AGCN200-1L	NanoXact Silver Nanospheres - Bare (Citrate) - 200 nm, 0.02 mg/mL in 2 mM sodium citrate, 1 L	\$1,200.00
<b>Silver Nanospheres, NanoXact, Dried</b>		
AGPD10-1MG	NanoXact Silver Nanospheres - PVP (Dried) - 10 nm, 1 mg (by silver mass)	\$135.00
AGPD10-5MG	NanoXact Silver Nanospheres - PVP (Dried) - 10 nm, 5 mg (by silver mass)	\$475.00
AGPD10-10MG	NanoXact Silver Nanospheres - PVP (Dried) - 10 nm, 10 mg (by silver mass)	\$800.00
AGPD10-30MG	NanoXact Silver Nanospheres - PVP (Dried) - 10 nm, 30 mg (by silver mass)	\$1,745.00
AGPD50-1MG	NanoXact Silver Nanospheres - PVP (Dried) - 50 nm, 1 mg (by silver mass)	\$135.00
AGPD50-5MG	NanoXact Silver Nanospheres - PVP (Dried) - 50 nm, 5 mg (by silver mass)	\$475.00

AGPD50-10MG	NanoXact Silver Nanospheres - PVP (Dried) - 50 nm, 10 mg (by silver mass)	\$800.00
AGPD50-30MG	NanoXact Silver Nanospheres - PVP (Dried) - 50 nm, 30 mg (by silver mass)	\$1,745.00
AGPD75-1MG	NanoXact Silver Nanospheres - PVP (Dried) - 75 nm, 1 mg (by silver mass)	\$135.00
AGPD75-5MG	NanoXact Silver Nanospheres - PVP (Dried) - 75 nm, 5 mg (by silver mass)	\$475.00
AGPD75-10MG	NanoXact Silver Nanospheres - PVP (Dried) - 75 nm, 10 mg (by silver mass)	\$800.00
AGPD75-30MG	NanoXact Silver Nanospheres - PVP (Dried) - 75 nm, 30 mg (by silver mass)	\$1,745.00
AGPD100-1MG	NanoXact Silver Nanospheres - PVP (Dried) - 100 nm, 1 mg (by silver mass)	\$135.00
AGPD100-5MG	NanoXact Silver Nanospheres - PVP (Dried) - 100 nm, 5 mg (by silver mass)	\$475.00
AGPD100-10MG	NanoXact Silver Nanospheres - PVP (Dried) - 100 nm, 10 mg (by silver mass)	\$800.00
AGPD100-30MG	NanoXact Silver Nanospheres - PVP (Dried) - 100 nm, 30 mg (by silver mass)	\$1,745.00

## Dodecanethiol-Stabilized Silver Nanospheres

Dodecanethiol-Stabilized Silver Nanospheres, NanoXact, Dried

AGDD4-1MG	NanoXact Silver Nanospheres - Dodecanethiol (Dried) - 4 nm, 1 mg	\$190.00
AGDD4-5MG	NanoXact Silver Nanospheres - Dodecanethiol (Dried) - 4 nm, 5 mg	\$720.00
AGDD4-25MG	NanoXact Silver Nanospheres - Dodecanethiol (Dried) - 4 nm, 25 mg	\$2,120.00
AGDD4-50MG	NanoXact Silver Nanospheres - Dodecanethiol (Dried) - 4 nm, 50 mg	\$3,375.00

## Silica-Shelled Silver Nanospheres

Silica-Shelled Silver Nanospheres, NanoXact, 1 mg/mL

AGSH50-1M	NanoXact Silver Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AGSH50-5M	NanoXact Silver Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AGSH50-10M	NanoXact Silver Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AGSH50-30M	NanoXact Silver Nanospheres - Silica Shelled - 90 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AGSH100-1M	NanoXact Silver Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AGSH100-5M	NanoXact Silver Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AGSH100-10M	NanoXact Silver Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AGSH100-30M	NanoXact Silver Nanospheres - Silica Shelled - 140 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AGAH50-1M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AGAH50-5M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AGAH50-10M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AGAH50-30M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 90 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00
AGAH100-1M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 1 mL	\$285.00
AGAH100-5M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 5 mL	\$1,130.00
AGAH100-10M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 10 mL	\$1,865.00
AGAH100-30M	NanoXact Silver Nanospheres - Silica Shelled (Aminated) - 140 nm, 1 mg/mL in ethanol, 30 mL	\$3,740.00

## Silver Nanoplates

Silver Nanoplates, NanoXact, 1 mg/mL

SPPH660-1M	NanoXact Silver Nanoplates - PVP - Peak @ 660 nm · 40 nm diam, 1 mg/mL in 5 mM sodium borate, 1 mL	\$285.00
SPPH660-5M	NanoXact Silver Nanoplates - PVP - Peak @ 660 nm · 40 nm diam, 1 mg/mL in 5 mM sodium borate, 5 mL	\$1,130.00
SPPH660-10M	NanoXact Silver Nanoplates - PVP - Peak @ 660 nm · 40 nm diam, 1 mg/mL in 5 mM sodium borate, 10 mL	\$1,805.00
SPPH660-30M	NanoXact Silver Nanoplates - PVP - Peak @ 660 nm · 40 nm diam, 1 mg/mL in 5 mM sodium borate, 30 mL	\$3,615.00
SPPH800-1M	NanoXact Silver Nanoplates - PVP - Peak @ 800 nm · 70 nm diam, 1 mg/mL in 5 mM sodium borate, 1 mL	\$285.00
SPPH800-5M	NanoXact Silver Nanoplates - PVP - Peak @ 800 nm · 70 nm diam, 1 mg/mL in 5 mM sodium borate, 5 mL	\$1,130.00
SPPH800-10M	NanoXact Silver Nanoplates - PVP - Peak @ 800 nm · 70 nm diam, 1 mg/mL in 5 mM sodium borate, 10 mL	\$1,805.00
SPPH800-30M	NanoXact Silver Nanoplates - PVP - Peak @ 800 nm · 70 nm diam, 1 mg/mL in 5 mM sodium borate, 30 mL	\$3,615.00
SPPH980-1M	NanoXact Silver Nanoplates - PVP - Peak @ 980 nm · 100 nm diam, 1 mg/mL in 5 mM sodium borate, 1 mL	\$285.00
SPPH980-5M	NanoXact Silver Nanoplates - PVP - Peak @ 980 nm · 100 nm diam, 1 mg/mL in 5 mM sodium borate, 5 mL	\$1,130.00
SPPH980-10M	NanoXact Silver Nanoplates - PVP - Peak @ 980 nm · 100 nm diam, 1 mg/mL in 5 mM sodium borate, 10 mL	\$1,805.00
SPPH980-30M	NanoXact Silver Nanoplates - PVP - Peak @ 980 nm · 100 nm diam, 1 mg/mL in 5 mM sodium borate, 30 mL	\$3,615.00
SPPH1064-1M	NanoXact Silver Nanoplates - PVP - Peak @ 1064 nm · 150 nm diam, 1 mg/mL in 5 mM sodium borate, 1 mL	\$285.00
SPPH1064-5M	NanoXact Silver Nanoplates - PVP - Peak @ 1064 nm · 150 nm diam, 1 mg/mL in 5 mM sodium borate, 5 mL	\$1,130.00
SPPH1064-10M	NanoXact Silver Nanoplates - PVP - Peak @ 1064 nm · 150 nm diam, 1 mg/mL in 5 mM sodium borate, 10 mL	\$1,805.00
SPPH1064-30M	NanoXact Silver Nanoplates - PVP - Peak @ 1064 nm · 150 nm diam, 1 mg/mL in 5 mM sodium borate, 30 mL	\$3,615.00

## Silica-Shelled Silver Nanoplates

Silica-Shelled Silver Nanoplates, NanoXact, 1 mg/mL

SPSH660-1M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 660 nm · 80 nm diam, 1 mg/mL, 1 mL	\$285.00
SPSH660-5M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 660 nm · 80 nm diam, 1 mg/mL, 5 mL	\$1,130.00
SPSH660-10M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 660 nm · 80 nm diam, 1 mg/mL, 10 mL	\$1,805.00
SPSH660-30M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 660 nm · 80 nm diam, 1 mg/mL, 30 mL	\$3,615.00
SPSH800-1M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 800 nm · 110 nm diam, 1 mg/mL, 1 mL	\$285.00
SPSH800-5M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 800 nm · 110 nm diam, 1 mg/mL, 5 mL	\$1,130.00
SPSH800-10M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 800 nm · 110 nm diam, 1 mg/mL, 10 mL	\$1,805.00
SPSH800-30M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 800 nm · 110 nm diam, 1 mg/mL, 30 mL	\$3,615.00
SPSH980-1M	NanoXact Silver Nanoplates - SiO <sub>2</sub> Shelled - Peak @ 980 nm · 140 nm diam, 1 mg/mL, 1 mL	\$285.00

SPSH980-5M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 980 nm · 140 nm diam, 1 mg/mL, 5 mL	\$1,130.00
SPSH980-10M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 980 nm · 140 nm diam, 1 mg/mL, 10 mL	\$1,805.00
SPSH980-30M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 980 nm · 140 nm diam, 1 mg/mL, 30 mL	\$3,615.00
SPSH1064-1M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 1064 nm · 190 nm diam, 1 mg/mL, 1 mL	\$285.00
SPSH1064-5M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 1064 nm · 190 nm diam, 1 mg/mL, 5 mL	\$1,130.00
SPSH1064-10M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 1064 nm · 190 nm diam, 1 mg/mL, 10 mL	\$1,805.00
SPSH1064-30M	NanoXact Silver Nanoplates - SiO2 Shelled - Peak @ 1064 nm · 190 nm diam, 1 mg/mL, 30 mL	\$3,615.00

## Silver Nanocubes

### Silver Nanocubes, NanoXact, 1 mg/mL

SCPH55-1M	NanoXact Silver Nanocubes - PVP - 55 nm, 1 mg/mL in ethanol, 1 mL	\$445.00
SCPH55-5M	NanoXact Silver Nanocubes - PVP - 55 nm, 1 mg/mL in ethanol, 5 mL	\$1,705.00
SCPH55-10M	NanoXact Silver Nanocubes - PVP - 55 nm, 1 mg/mL in ethanol, 10 mL	\$2,765.00
SCPH55-30M	NanoXact Silver Nanocubes - PVP - 55 nm, 1 mg/mL in ethanol, 30 mL	\$5,875.00
SCPH75-1M	NanoXact Silver Nanocubes - PVP - 75 nm, 1 mg/mL in ethanol, 1 mL	\$445.00
SCPH75-5M	NanoXact Silver Nanocubes - PVP - 75 nm, 1 mg/mL in ethanol, 5 mL	\$1,705.00
SCPH75-10M	NanoXact Silver Nanocubes - PVP - 75 nm, 1 mg/mL in ethanol, 10 mL	\$2,765.00
SCPH75-30M	NanoXact Silver Nanocubes - PVP - 75 nm, 1 mg/mL in ethanol, 30 mL	\$5,875.00
SCPH100-1M	NanoXact Silver Nanocubes - PVP - 100 nm, 1 mg/mL in ethanol, 1 mL	\$445.00
SCPH100-5M	NanoXact Silver Nanocubes - PVP - 100 nm, 1 mg/mL in ethanol, 5 mL	\$1,705.00
SCPH100-10M	NanoXact Silver Nanocubes - PVP - 100 nm, 1 mg/mL in ethanol, 10 mL	\$2,765.00
SCPH100-30M	NanoXact Silver Nanocubes - PVP - 100 nm, 1 mg/mL in ethanol, 30 mL	\$5,875.00

## Silica Nanospheres

### Silica Nanospheres, NanoXact, 5 mg/mL

SISN20-25M	NanoXact Silica Nanospheres - 20 nm, 5 mg/mL in water, 25 mL	\$235.00
SISN20-100M	NanoXact Silica Nanospheres - 20 nm, 5 mg/mL in water, 100 mL	\$875.00
SISN20-500M	NanoXact Silica Nanospheres - 20 nm, 5 mg/mL in water, 500 mL	\$3,980.00
SISN20-1L	NanoXact Silica Nanospheres - 20 nm, 5 mg/mL in water, 1 L	\$7,120.00

SIAN20-25M	NanoXact Silica Nanospheres - Aminated - 20 nm, 5 mg/mL in ethanol, 25 mL	\$235.00
SIAN20-100M	NanoXact Silica Nanospheres - Aminated - 20 nm, 5 mg/mL in ethanol, 100 mL	\$875.00
SIAN20-500M	NanoXact Silica Nanospheres - Aminated - 20 nm, 5 mg/mL in ethanol, 500 mL	\$3,980.00
SIAN20-1L	NanoXact Silica Nanospheres - Aminated - 20 nm, 5 mg/mL in ethanol, 1 L	\$7,120.00

### Silica Nanospheres, NanoXact, 10 mg/mL

SISN50-25M	NanoXact Silica Nanospheres - 50 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN50-100M	NanoXact Silica Nanospheres - 50 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN50-500M	NanoXact Silica Nanospheres - 50 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN50-1L	NanoXact Silica Nanospheres - 50 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN60-25M	NanoXact Silica Nanospheres - 60 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN60-100M	NanoXact Silica Nanospheres - 60 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN60-500M	NanoXact Silica Nanospheres - 60 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN60-1L	NanoXact Silica Nanospheres - 60 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN80-25M	NanoXact Silica Nanospheres - 80 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN80-100M	NanoXact Silica Nanospheres - 80 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN80-500M	NanoXact Silica Nanospheres - 80 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN80-1L	NanoXact Silica Nanospheres - 80 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN100-25M	NanoXact Silica Nanospheres - 100 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN100-100M	NanoXact Silica Nanospheres - 100 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN100-500M	NanoXact Silica Nanospheres - 100 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN100-1L	NanoXact Silica Nanospheres - 100 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN120-25M	NanoXact Silica Nanospheres - 120 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN120-100M	NanoXact Silica Nanospheres - 120 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN120-500M	NanoXact Silica Nanospheres - 120 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN120-1L	NanoXact Silica Nanospheres - 120 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN200-25M	NanoXact Silica Nanospheres - 200 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN200-100M	NanoXact Silica Nanospheres - 200 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN200-500M	NanoXact Silica Nanospheres - 200 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN200-1L	NanoXact Silica Nanospheres - 200 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN300-25M	NanoXact Silica Nanospheres - 300 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN300-100M	NanoXact Silica Nanospheres - 300 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN300-500M	NanoXact Silica Nanospheres - 300 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN300-1L	NanoXact Silica Nanospheres - 300 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN500-25M	NanoXact Silica Nanospheres - 500 nm, 10 mg/mL in water, 25 mL	\$210.00
SISN500-100M	NanoXact Silica Nanospheres - 500 nm, 10 mg/mL in water, 100 mL	\$755.00
SISN500-500M	NanoXact Silica Nanospheres - 500 nm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN500-1L	NanoXact Silica Nanospheres - 500 nm, 10 mg/mL in water, 1 L	\$6,030.00
SISN1000-25M	NanoXact Silica Nanospheres - 1 µm, 10 mg/mL in water, 25 mL	\$210.00
SISN1000-100M	NanoXact Silica Nanospheres - 1 µm, 10 mg/mL in water, 100 mL	\$755.00
SISN1000-500M	NanoXact Silica Nanospheres - 1 µm, 10 mg/mL in water, 500 mL	\$3,450.00
SISN1000-1L	NanoXact Silica Nanospheres - 1 µm, 10 mg/mL in water, 1 L	\$6,030.00

SIAN50-25M	NanoXact Silica Nanospheres - Aminated - 50 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIAN50-100M	NanoXact Silica Nanospheres - Aminated - 50 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN50-500M	NanoXact Silica Nanospheres - Aminated - 50 nm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIAN50-1L	NanoXact Silica Nanospheres - Aminated - 50 nm, 10 mg/mL in ethanol, 1 L	\$6,030.00
SIAN80-25M	NanoXact Silica Nanospheres - Aminated - 80 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIAN80-100M	NanoXact Silica Nanospheres - Aminated - 80 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN80-500M	NanoXact Silica Nanospheres - Aminated - 80 nm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIAN80-1L	NanoXact Silica Nanospheres - Aminated - 80 nm, 10 mg/mL in ethanol, 1 L	\$6,030.00
SIAN100-25M	NanoXact Silica Nanospheres - Aminated - 100 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIAN100-100M	NanoXact Silica Nanospheres - Aminated - 100 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN100-500M	NanoXact Silica Nanospheres - Aminated - 100 nm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIAN100-1L	NanoXact Silica Nanospheres - Aminated - 100 nm, 10 mg/mL in ethanol, 1 L	\$6,030.00
SIANI20-25M	NanoXact Silica Nanospheres - Aminated - 120 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIANI20-100M	NanoXact Silica Nanospheres - Aminated - 120 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN200-25M	NanoXact Silica Nanospheres - Aminated - 200 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIAN200-100M	NanoXact Silica Nanospheres - Aminated - 200 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN200-500M	NanoXact Silica Nanospheres - Aminated - 200 nm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIAN200-1L	NanoXact Silica Nanospheres - Aminated - 200 nm, 10 mg/mL in ethanol, 1 L	\$6,030.00
SIAN300-25M	NanoXact Silica Nanospheres - Aminated - 300 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIAN300-100M	NanoXact Silica Nanospheres - Aminated - 300 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN300-500M	NanoXact Silica Nanospheres - Aminated - 300 nm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIAN300-1L	NanoXact Silica Nanospheres - Aminated - 300 nm, 10 mg/mL in ethanol, 1 L	\$6,030.00
SIAN500-25M	NanoXact Silica Nanospheres - Aminated - 500 nm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIAN500-100M	NanoXact Silica Nanospheres - Aminated - 500 nm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIAN500-500M	NanoXact Silica Nanospheres - Aminated - 500 nm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIAN500-1L	NanoXact Silica Nanospheres - Aminated - 500 nm, 10 mg/mL in ethanol, 1 L	\$6,030.00
SIANI000-25M	NanoXact Silica Nanospheres - Aminated - 1 µm, 10 mg/mL in ethanol, 25 mL	\$210.00
SIANI000-100M	NanoXact Silica Nanospheres - Aminated - 1 µm, 10 mg/mL in ethanol, 100 mL	\$755.00
SIANI000-500M	NanoXact Silica Nanospheres - Aminated - 1 µm, 10 mg/mL in ethanol, 500 mL	\$3,450.00
SIANI000-1L	NanoXact Silica Nanospheres - Aminated - 1 µm, 10 mg/mL in ethanol, 1 L	\$6,030.00

#### Silica Nanospheres, NanoXact, Dried

SISD80-250MG	NanoXact Silica Nanospheres - (Dried) - 80 nm, 250 mg	\$210.00
SISD80-1G	NanoXact Silica Nanospheres - (Dried) - 80 nm, 1 g	\$755.00
SISD80-5G	NanoXact Silica Nanospheres - (Dried) - 80 nm, 5 g	\$3,450.00
SISD80-10G	NanoXact Silica Nanospheres - (Dried) - 80 nm, 10 g	\$6,030.00
SISD100-250MG	NanoXact Silica Nanospheres - (Dried) - 100 nm, 250 mg	\$210.00
SISD100-1G	NanoXact Silica Nanospheres - (Dried) - 100 nm, 1 g	\$755.00
SISD100-5G	NanoXact Silica Nanospheres - (Dried) - 100 nm, 5 g	\$3,450.00
SISD100-10G	NanoXact Silica Nanospheres - (Dried) - 100 nm, 10 g	\$6,030.00
SISD120-250MG	NanoXact Silica Nanospheres - (Dried) - 120 nm, 250 mg	\$210.00
SISD120-1G	NanoXact Silica Nanospheres - (Dried) - 120 nm, 1 g	\$755.00
SISD120-5G	NanoXact Silica Nanospheres - (Dried) - 120 nm, 5 g	\$3,450.00
SISD120-10G	NanoXact Silica Nanospheres - (Dried) - 120 nm, 10 g	\$6,030.00
SISD200-250MG	NanoXact Silica Nanospheres - (Dried) - 200 nm, 250 mg	\$210.00
SISD200-1G	NanoXact Silica Nanospheres - (Dried) - 200 nm, 1 g	\$755.00
SISD200-5G	NanoXact Silica Nanospheres - (Dried) - 200 nm, 5 g	\$3,450.00
SISD200-10G	NanoXact Silica Nanospheres - (Dried) - 200 nm, 10 g	\$6,030.00
SISD300-250MG	NanoXact Silica Nanospheres - (Dried) - 300 nm, 250 mg	\$210.00
SISD300-1G	NanoXact Silica Nanospheres - (Dried) - 300 nm, 1 g	\$755.00
SISD300-5G	NanoXact Silica Nanospheres - (Dried) - 300 nm, 5 g	\$3,450.00
SISD300-10G	NanoXact Silica Nanospheres - (Dried) - 300 nm, 10 g	\$6,030.00
SISD500-250MG	NanoXact Silica Nanospheres - (Dried) - 500 nm, 250 mg	\$210.00
SISD500-1G	NanoXact Silica Nanospheres - (Dried) - 500 nm, 1 g	\$755.00
SISD500-5G	NanoXact Silica Nanospheres - (Dried) - 500 nm, 5 g	\$3,450.00
SISD500-10G	NanoXact Silica Nanospheres - (Dried) - 500 nm, 10 g	\$6,030.00
SISD1000-250MG	NanoXact Silica Nanospheres - (Dried) - 1 µm, 250 mg	\$210.00
SISD1000-1G	NanoXact Silica Nanospheres - (Dried) - 1 µm, 1 g	\$755.00
SISD1000-5G	NanoXact Silica Nanospheres - (Dried) - 1 µm, 5 g	\$3,450.00
SISD1000-10G	NanoXact Silica Nanospheres - (Dried) - 1 µm, 10 g	\$6,030.00
SIAD80-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 80 nm, 250 mg	\$210.00
SIAD80-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 80 nm, 1 g	\$755.00
SIAD80-5G	NanoXact Silica Nanospheres - Aminated (Dried) - 80 nm, 5 g	\$3,450.00
SIAD80-10G	NanoXact Silica Nanospheres - Aminated (Dried) - 80 nm, 10 g	\$6,030.00
SIAD100-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 100 nm, 250 mg	\$210.00
SIAD100-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 100 nm, 1 g	\$755.00
SIAD100-5G	NanoXact Silica Nanospheres - Aminated (Dried) - 100 nm, 5 g	\$3,450.00
SIAD100-10G	NanoXact Silica Nanospheres - Aminated (Dried) - 100 nm, 10 g	\$6,030.00
SIAD120-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 120 nm, 250 mg	\$210.00
SIAD120-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 120 nm, 1 g	\$755.00
SIAD200-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 200 nm, 250 mg	\$210.00
SIAD200-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 200 nm, 1 g	\$755.00
SIAD200-5G	NanoXact Silica Nanospheres - Aminated (Dried) - 200 nm, 5 g	\$3,450.00

SIAD200-10G	NanoXact Silica Nanospheres - Aminated (Dried) - 200 nm, 10 g	\$6,030.00
SIAD300-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 300 nm, 250 mg	\$210.00
SIAD300-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 300 nm, 1 g	\$755.00
SIAD300-5G	NanoXact Silica Nanospheres - Aminated (Dried) - 300 nm, 5 g	\$3,450.00
SIAD300-10G	NanoXact Silica Nanospheres - Aminated (Dried) - 300 nm, 10 g	\$6,030.00
SIAD500-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 500 nm, 250 mg	\$210.00
SIAD500-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 500 nm, 1 g	\$755.00
SIAD500-5G	NanoXact Silica Nanospheres - Aminated (Dried) - 500 nm, 5 g	\$3,450.00
SIAD500-10G	NanoXact Silica Nanospheres - Aminated (Dried) - 500 nm, 10 g	\$6,030.00
SIAD1000-250MG	NanoXact Silica Nanospheres - Aminated (Dried) - 1 µm, 250 mg	\$210.00
SIAD1000-1G	NanoXact Silica Nanospheres - Aminated (Dried) - 1 µm, 1 g	\$755.00
SIAD1000-5G	NanoXact Silica Nanospheres - Aminated (Dried) - 1 µm, 5 g	\$3,450.00
SIAD1000-10G	NanoXact Silica Nanospheres - Aminated (Dried) - 1 µm, 10 g	\$6,030.00

## Mesoporous Silica Nanospheres

### Mesoporous Silica Nanospheres, NanoXact, 5 mg/mL

SHSN100-5M	NanoXact Mesoporous Silica Nanospheres - 100 nm, 5 mg/mL in ethanol, 5 mL	\$280.00
SHSN100-30M	NanoXact Mesoporous Silica Nanospheres - 100 nm, 5 mg/mL in ethanol, 30 mL	\$1,430.00
SHSN100-100M	NanoXact Mesoporous Silica Nanospheres - 100 nm, 5 mg/mL in ethanol, 100 mL	\$4,090.00
SHAN100-5M	NanoXact Mesoporous Silica Nanospheres - Aminated - 100 nm, 5 mg/mL in ethanol, 5 mL	\$280.00
SHAN100-30M	NanoXact Mesoporous Silica Nanospheres - Aminated - 100 nm, 5 mg/mL in ethanol, 30 mL	\$1,430.00
SHAN100-100M	NanoXact Mesoporous Silica Nanospheres - Aminated - 100 nm, 5 mg/mL in ethanol, 100 mL	\$4,090.00

### Mesoporous Silica Nanospheres, NanoXact, 5 mg/mL

SHSD100-25MG	NanoXact Mesoporous Silica Nanospheres - (Dried) - 100 nm, 25 mg	\$280.00
SHSD100-150MG	NanoXact Mesoporous Silica Nanospheres - (Dried) - 100 nm, 150 mg	\$1,430.00
SHSD100-500MG	NanoXact Mesoporous Silica Nanospheres - (Dried) - 100 nm, 500 mg	\$4,090.00
SHAD100-25MG	NanoXact Mesoporous Silica Nanospheres - Aminated (Dried) - 100 nm, 25 mg	\$280.00
SHAD100-150MG	NanoXact Mesoporous Silica Nanospheres - Aminated (Dried) - 100 nm, 150 mg	\$1,430.00
SHAD100-500MG	NanoXact Mesoporous Silica Nanospheres - Aminated (Dried) - 100 nm, 500 mg	\$4,090.00

## Magnetite Nanoparticles

### Magnetite Nanoparticles, BioPure, 1 mg/mL

MGPB20-1M	BioPure Magnetite Nanoparticles - PVP - 20 nm, 1 mg/mL, in aq. 2 mM sodium citrate solution, 1 mL	\$285.00
MGPB20-5M	BioPure Magnetite Nanoparticles - PVP - 20 nm, 1 mg/mL, in aq. 2 mM sodium citrate solution, 5 mL	\$1,130.00
MGPB20-10M	BioPure Magnetite Nanoparticles - PVP - 20 nm, 1 mg/mL, in aq. 2 mM sodium citrate solution, 10 mL	\$1,865.00
MGPB20-30M	BioPure Magnetite Nanoparticles - PVP - 20 nm, 1 mg/mL, in aq. 2 mM sodium citrate solution, 30 mL	\$3,740.00

### Magnetite Nanoparticles, NanoXact, 20 mg/mL

MGPN20-5M	NanoXact Magnetite Nanoparticles - PVP - 20 nm, 20 mg/mL, in aq. 2 mM sodium citrate sol., 5 mL	\$245.00
MGPN20-10M	NanoXact Magnetite Nanoparticles - PVP - 20 nm, 20 mg/mL, in aq. 2 mM sodium citrate sol., 10 mL	\$370.00
MGPN20-25M	NanoXact Magnetite Nanoparticles - PVP - 20 nm, 20 mg/mL, in aq. 2 mM sodium citrate sol., 25 mL	\$795.00

### Magnetite Nanoparticles, NanoXact, Dried

MGPD20-100MG	NanoXact Magnetite Nanoparticles - PVP (Dried) - 20 nm, 100 mg	\$245.00
MGPD20-200MG	NanoXact Magnetite Nanoparticles - PVP (Dried) - 20 nm, 200 mg	\$370.00
MGPD20-500MG	NanoXact Magnetite Nanoparticles - PVP (Dried) - 20 nm, 500 mg	\$795.00

## Platinum Nanoparticles

### Platinum Nanoparticles, BioPure, 1 mg/mL

PTCB5-1M	BioPure Platinum Nanoparticles - Bare (Citrate) - 5 nm, 1 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$285.00
PTCB5-5M	BioPure Platinum Nanoparticles - Bare (Citrate) - 5 nm, 1 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$1,130.00
PTCB5-10M	BioPure Platinum Nanoparticles - Bare (Citrate) - 5 nm, 1 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$1,865.00
PTCB5-30M	BioPure Platinum Nanoparticles - Bare (Citrate) - 5 nm, 1 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$3,740.00
PTCB30-1M	BioPure Platinum Nanoparticles - Bare (Citrate) - 30 nm, 1 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$285.00
PTCB30-5M	BioPure Platinum Nanoparticles - Bare (Citrate) - 30 nm, 1 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$1,130.00
PTCB30-10M	BioPure Platinum Nanoparticles - Bare (Citrate) - 30 nm, 1 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$1,865.00
PTCB30-30M	BioPure Platinum Nanoparticles - Bare (Citrate) - 30 nm, 1 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$3,740.00
PTCB50-1M	BioPure Platinum Nanoparticles - Bare (Citrate) - 50 nm, 1 mg/mL in aq. 2 mM sodium citrate, 1 mL	\$285.00
PTCB50-5M	BioPure Platinum Nanoparticles - Bare (Citrate) - 50 nm, 1 mg/mL in aq. 2 mM sodium citrate, 5 mL	\$1,130.00
PTCB50-10M	BioPure Platinum Nanoparticles - Bare (Citrate) - 50 nm, 1 mg/mL in aq. 2 mM sodium citrate, 10 mL	\$1,865.00
PTCB50-30M	BioPure Platinum Nanoparticles - Bare (Citrate) - 50 nm, 1 mg/mL in aq. 2 mM sodium citrate, 30 mL	\$3,740.00
PTCB70-1M	BioPure Platinum Nanoparticles - Bare (Citrate) - 70 nm, 1 mg/mL in aq. 4 mM sodium citrate, 1 mL	\$285.00
PTCB70-5M	BioPure Platinum Nanoparticles - Bare (Citrate) - 70 nm, 1 mg/mL in aq. 4 mM sodium citrate, 5 mL	\$1,130.00
PTCB70-10M	BioPure Platinum Nanoparticles - Bare (Citrate) - 70 nm, 1 mg/mL in aq. 4 mM sodium citrate, 10 mL	\$1,865.00
PTCB70-30M	BioPure Platinum Nanoparticles - Bare (Citrate) - 70 nm, 1 mg/mL in aq. 4 mM sodium citrate, 30 mL	\$3,740.00

### Platinum Nanoparticles, NanoXact, 0.05 mg/mL

PTCN5-50M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 5 nm, 0.05 mg/mL in 1 mM sodium citrate, 50 mL	\$335.00
PTCN5-100M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 5 nm, 0.05 mg/mL in 2 mM sodium citrate, 100 mL	\$520.00
PTCN5-500M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 5 nm, 0.05 mg/mL in 2 mM sodium citrate, 500 mL	\$1,965.00
PTCN5-1L	NanoXact Platinum Nanoparticles - Bare (Citrate) - 5 nm, 0.05 mg/mL in 2 mM sodium citrate, 1 L	\$3,275.00
PTCN30-50M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 30 nm, 0.05 mg/mL in 2 mM sodium citrate, 50 mL	\$335.00
PTCN30-100M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 30 nm, 0.05 mg/mL in 2 mM sodium citrate, 100 mL	\$520.00
PTCN30-500M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 30 nm, 0.05 mg/mL in 2 mM sodium citrate, 500 mL	\$1,965.00
PTCN30-1L	NanoXact Platinum Nanoparticles - Bare (Citrate) - 30 nm, 0.05 mg/mL in 2 mM sodium citrate, 1 L	\$3,275.00
PTCN50-50M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 50 nm, 0.05 mg/mL in 2 mM sodium citrate, 50 mL	\$335.00
PTCN50-100M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 50 nm, 0.05 mg/mL in 2 mM sodium citrate, 100 mL	\$520.00
PTCN50-500M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 50 nm, 0.05 mg/mL in 2 mM sodium citrate, 500 mL	\$1,965.00
PTCN50-1L	NanoXact Platinum Nanoparticles - Bare (Citrate) - 50 nm, 0.05 mg/mL in 2 mM sodium citrate, 1 L	\$3,275.00
PTCN70-50M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 70 nm, 0.05 mg/mL in 4 mM sodium citrate, 50 mL	\$335.00
PTCN70-100M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 70 nm, 0.05 mg/mL in 4 mM sodium citrate, 100 mL	\$520.00
PTCN70-500M	NanoXact Platinum Nanoparticles - Bare (Citrate) - 70 nm, 0.05 mg/mL in 4 mM sodium citrate, 500 mL	\$1,965.00
PTCN70-1L	NanoXact Platinum Nanoparticles - Bare (Citrate) - 70 nm, 0.05 mg/mL in 4 mM sodium citrate, 1 L	\$3,275.00

## Lateral Flow Assay Kits & Consumables

### Lateral Flow Material Starter Kit

MSKR	BioReady Material Starter Kit for Lateral Flow - Full Kit	\$2,410.00
NCMR	Lateral Flow Material Starter Kit - Nitrocellulose Membranes Only	\$900.00

### Assay Kits

ADLF-25T	AmpliDetect - Nucleic Acid Lateral Flow (NALF) assay	\$185.00
----------	--	----------

### Buffers

RXKR-100M	BioReady Potassium Phosphate Reaction Buffer - 5mM Potassium Phosphate, 0.5% PEG20K, pH 7.4, 100 mL	\$185.00
RXKR-500M	BioReady Potassium Phosphate Reaction Buffer - 5mM Potassium Phosphate, 0.5% PEG20K, pH 7.4, 500 mL	\$545.00
RXKR-1L	BioReady Potassium Phosphate Reaction Buffer - 5mM Potassium Phosphate, 0.5% PEG20K, pH 7.4, 1 L	\$820.00
RXNR-100M	BioReady Sodium Phosphate Reaction Buffer - 5mM Sodium Phosphate, 0.5% PEG20K, pH 7.4, 100 mL	\$185.00
RXNR-500M	BioReady Sodium Phosphate Reaction Buffer - 5mM Sodium Phosphate, 0.5% PEG20K, pH 7.4, 500 mL	\$545.00
RXNR-1L	BioReady Sodium Phosphate Reaction Buffer - 5mM Sodium Phosphate, 0.5% PEG20K, pH 7.4, 1 L	\$820.00
RXPR-100M	BioReady PBS Reaction Buffer - 0.01X PBS, 0.5% PEG20K, pH 7.4, 100 mL	\$185.00