

Small Universal Centrifuge

LABEC offer compact benchtop **Small** and **Digital Universal Centrifuges** with the capacity for rotors ranging from 1.5ml to 50ml and speeds of 500-6000rpm's.

M-SCEN-207 / Small Universal Centrifuge

Features:

- Analog speed tachometer.
- Optional rotor selections for various tube sizes.
- Regressive timer: 30min, or continuous operation.
- Angle rotor, stepless speed setting.
- Auto accelerating & decelerating functions.
- Auto-current restraining function.
- Power pilot light & fuse safety device.
- Manual cover lock.
- Speed: 5,500rpm; Timer: 0~30 or "ON".
- Max. capacity: 300ml.
- Weight: 9kg.
- Size: W29 x D32 x H25cm.
- Auto power off when open & auto-balance device.
- ABS hard shell & see-through SCEN-207 window.
- Auto-cooling air circulation.



M-SCEN-207



M-RA-1508
15ml x 8pcs
Max. 5,500rpm

M-RA-1512S
15ml x 12pcs
Max. 5,500rpm

M-RA-5004
50ml x 4pcs
Max. 5,500rpm

M-RA-5006
50ml x 6pcs
Max. 5,500rpm

M-REA-1520
1.5ml x 20pcs
Max. 4,000/5,500rpm

Optional Rotors	Max. Speed (rpm)	Max. RCF (g)	Capacity (ml x pcs)	Tube size (Ø x L) mm
M-RA-1508	0~5,500	3,140	15 x 8	15~16.7 x 90~122
M-RA-1521S		3,280	15 x 12	15~16.7 x 90~122
M-RA-5004		3,240	50 x 4	27.5~29.5 x 90~117
M-RA-5006		3,240	50 x 6	27.5~29.5 x 90~117
M-REA-1520		2,600	1.5 x 20	10~11 x 36~56

Digital Universal Centrifuge

M-SCEN-206 / Digital Universal Centrifuge

Features:

- LED speed and time display.
- Microprocessor control.
- 10 memory programs for preset parameters.
- Auto accelerating and decelerating functions.
- Auto cover lock with motor running.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- Speed: 500~6000rpm;
- Timer: 99min 59 sec
- Max. Capacity: 300ml.
- Weight: 10kg.
- Size: W29 x D32 x H25cm.
- ABS hard shell and see-through window.
- Auto-cooling air circulation.
- Optional rotor selections for various tube size.



M-SCEN-206



M-RA-1508
15ml x 8pcs
Max. 6,000rpm

M-RA-1512S
15ml x 12pcs
Max. 6,000rpm

M-RA-5004
50ml x 4pcs
Max. 6,000rpm

M-RA-5006
50ml x 6pcs
Max. 6,000rpm

M-REA-1520
1.5ml x 20pcs
Max. 6,000rpm

Optional Rotors	Max. Speed (rpm)	Max. RCF (g)	Capacity (ml x pcs)	Tube size (Ø x L) mm
M-RA-1508	6,000	3,740	15 x 8	15~16.7 x 90~122
M-RA-1521S		3,900	15 x 12	15~16.7 x 90~122
M-RA-5004		3,860	50 x 4	27.5~29.5 x 90~117
M-RA-5006		3,860	50 x 6	27.5~29.5 x 90~117
M-REA-1520		3,090	1.5 x 20	10~11 x 36~56