

Clima-Flex[™] Knee Brace More Elexibility

More Flexibility, Better Mobility



+* MOTION IS **MEDICINE***



Clima-Flex[™] Knee Brace More Flexibility, **Better Mobility**


clima = Temperature Control flex = Mobility & Comfort

Clima-Flex[™] OA offers a new level of comfortable, cooling relief through compression and offloading for patients with mild to moderate OA knee pain. New DonJoy C-6 Fabric Technology delivers breathable, moisture-wicking comfort that's soft to the touch, antimicrobial and incorporates temperature control. Keyless offloading enables on-the-fly hinge adjustment to empower patients with convenient, comforting relief from OA knee pain. Clima-Flex[™] OA is sleek and ideal for light activities such as walking, golf, or hiking.

NEW

C-6 Technology[™] Advanced Cooling Material

Pure carbon threads are woven into this moisture-wicking, highly breathable and cooling antimicrobial fabric to naturally repel odor and provide temperature control.



E-Z Wraparound Application

E-Z Grip handles and intuitive wraparound design ensure form-fitting comfort; brace is easily adjustable for more or less compression based on patient needs and comfort.

Perfect Fit & Control

Semi-rigid cuffs concealed in the thigh and calf provide comfortable support and confidence.

Keyless, Dual Dial Off Loading

New, easy-to-use integrated dual-action dial allows on-the-fly hinge adjustment for comfortable, convenient offloading of the knee anytime, anywhere.





ORDERING INFORMATION

Part No	Description
11-8810- <mark>X</mark>	Clima-Flex OA, Right Medial/Left Lateral
11-8811- <mark>X</mark>	Clima-Flex OA, Left medial/Right Lateral

Sizes: XS-XXXL (1-7)

(Measurements taken 6" above mid-patella) X = 1 = XS 13" - 15½" (33-39 CM) X = 2 = S 15½" - 18½" (39-47 CM) X = 3 = M 18½" - 21" (47-53 CM) X = 4 = L 21" - 23½" (53-60 CM) X = 5 = XL 23½" - 26½" (60-67 CM) X = 6 = XXL 26½" - 29½" (67-75 CM) X = 7 = XXXL 29½" - 32" (75-81 CM)



90 Day Test Drive - Simply fit your patients in any DonJoy knee osteoarthritis brace. If they are not 100% satisfied with the brace, you may return it for a full credit or refund within 90 days of brace fitting.



Products with the Arthritis Foundation's Ease-of-Use Commendation are tested and proven to make daily tasks easier. Learn more at arthritis.org/easeofuse.



 DJO, LLC
 I
 A DJO Global Company

 T
 800.336.6569
 D
 760.727.1280
 F
 800.936.6569

 1430
 Decision Street
 I
 Vista, CA
 92081-8553
 I
 U.S.A.

 DJOglobal.com
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D
 D





DR