TTTI I III dina

Procedure Trolleys & Carts



Supplied proudly in New Zealand by



Caretray Trolley - Anaesthesia

Features

- Trolleys suitable for use with Anaesthesia
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
 - High Security Bolt Lock Supplied with two keys
 - · Electronic Push Button Lock Programmable four-digit code
- Trolley pre-configured with the following accessories...
 - Overbridge, Mild Steel
 - Four Compartment Tilt Bin, Mild Steel
 - · Five Compartment Tilt Bin, Mild Steel
 - Waste Bin & Bracket, 8 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Grey White

| Description | High Security Bolt Lock | Electronic Push Button Lock | Weight |
|----------------------------|-------------------------|-----------------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTANAEM1/10/3S2D/HSB | CTANAEM1/10/3S2D/EPB | 41.4kg |
| 4 x 100mm, 1 x 150mm Trays | CTANAEM1/10/4S1D/HSB | CTANAEM1/10/4S1D/EPB | 41.1kg |
| 6 x 100mm Trays | CTANAEM1/10/6S/HSB | CTANAEM1/10/6S/EPB | 42.2kg |

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Carerail Mounted page 62
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71





CTANAEM1/10/3S2D/HSB

Carts

లర

Trolleys

Procedure

Caretray Trolley - Cast

Features

- Trolleys suitable for use within A&E / Plaster Room
- · Quick & easy disassembly, assists thorough & deep cleaning
- · Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
 - · CAM Lock with Wing Handle Supplied with two keys
 - · Electronic Push Button Lock Programmable four-digit code
- Trolley pre-configured with the following accessories...
 - Overbridge, Mild Steel
 - Four Compartment Tilt Bin, Mild Steel
 - · Bowl Support & Bracket
 - Stainless Steel Bowl, 7 Litre
 - Waste Bin & Bracket, 8 Litre
 - · Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- · Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Grey White

| Description | CAM Lock with Wing Handle | Electronic Push Button Lock | Weight |
|----------------------------|---------------------------|-----------------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTCASTM1/10/3S2D/CA | CTCASTM1/10/3S2D/EPB | 43.5kg |
| 4 x 100mm, 1 x 150mm Trays | CTCASTM1/10/4S1D/CA | CTCASTM1/10/4S1D/EPB | 43.2kg |
| 6 x 100mm Trays | CTCASTM1/10/6S/CA | CTCASTM1/10/6S/EPB | 44.3kg |

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Carerail Mounted page 62
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71



Caretray Trolley - Critical Care

Features

- Trolleys suitable for use with Critical Care
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
 - · CAM Lock with Winged Handle Supplied with two keys
 - · Electronic Push Button Lock Programmable four-digit code
- Trolley pre-configured with the following accessories...
 - Overbridge, Mild Steel
 - Five Compartment Tilt Bin, Mild Steel
 - Six Compartment Tilt Bin, Mild Steel
 - Alcohol Gel Holder
 - Infusion Pole & Bracket
 - Waste Bin & Bracket, 8 Litre
 - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

| Description | CAM Lock with Wing Handle | Electronic Push Button Lock | Weight |
|----------------------------|---------------------------|-----------------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTCRITM1/10/3S2D/CA | CTCRITM1/10/3S2D/EPB | 44.6kg |
| 4 x 100mm, 1 x 150mm Trays | CTCRITM1/10/4S1D/CA | CTCRITM1/10/4S1D/EPB | 44.3kg |
| 6 x 100mm Trays | CTCRITM1/10/6S/CA | CTCRITM1/10/6S/EPB | 45.4kg |

Standard Colours

· Paint - Grey White

Options / Accessories

- · Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Carerail Mounted page 62
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71



Carts

න්

Trolleys

Procedure

Caretray Trolley - Infection

Features

- Trolleys designed for use when treating patients that have been infected with multi-resistant pathogens
- Quick & easy disassembly, assists thorough & deep cleaning
- · Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
 Two lacking actions
- Two locking options...
 - High Security Bolt Lock Supplied with two keys
 - Electronic Push Button Lock Programmable four-digit code
- Trolley pre-configured with the following accessories...
 - Waste Bin & Bracket, 8 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied - white, blue, green, yellow & red)
- · Each tray supplied with 2 x long & 2 x short dividers
- · Max load per tray 5kg
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Grey White

| Description | High Security Bolt Lock | Electronic Push Button Lock | Weight |
|----------------------------|-------------------------|-----------------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTINFEM1/10/3S2D/HSB | CTINFEM1/10/3S2D/EPB | 35.5kg |
| 4 x 100mm, 1 x 150mm Trays | CTINFEM1/10/4S1D/HSB | CTINFEM1/10/4S1D/EPB | 35.2kg |
| 6 x 100mm Trays | CTINFEM1/10/6S/HSB | CTINFEM1/10/6S/EPB | 36.3kg |

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Defibrillator Platform page 61
- Carerail Mounted page 62
- Laptop Holders page 68
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71

CTINFEM1/10/3S2D/HSE

Caretray Trolley - IV Preparation

Features

- Trolleys designed for use when preparing & administering intravenous infusions
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
 - High Security Bolt Lock Supplied with two keys
 - · Electronic Push Button Lock Programmable four-digit code
- Trolley pre-configured with the following accessories...
 - Infusion Pole & Bracket
 - Waste Bin & Bracket, 8 Litre
 - · Glove / Mask Box Holder, One Box
 - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Grey White

| Description | High Security Bolt Lock | Electronic Push Button Lock | Weight |
|----------------------------|-------------------------|-----------------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTIVPRM1/10/3S2D/HSB | CTIVPRM1/10/3S2D/EPB | 39.8kg |
| 4 x 100mm, 1 x 150mm Trays | CTIVPRM1/10/4S1D/HSB | CTIVPRM1/10/4S1D/EPB | 39.5kg |
| 6 x 100mm Trays | CTIVPRM1/10/6S/HSB | CTIVPRM1/10/6S/EPB | 40.6kg |

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Defibrillator Platform page 61
- Carerail Mounted page 62
- Laptop Holders page 68
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71





Carts

ංඊ

Trolleys

Procedure

Caretray Trolley - Medical / Surgical

Features

- Trolleys designed for use within the Medical / Surgical Environment
- Quick & easy disassembly, assists thorough & deep cleaning
- · Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
 - · CAM Lock with Wing Handle Supplied with two keys
 - · Electronic Push Button Lock Programmable four-digit code
- Trolley pre-configured with the following accessories...
 - Overbridge, Mild Steel
 - Four Compartment Tilt Bin, Mild Steel
 - Infusion Pole & Bracket
 - Waste Bin & Bracket, 8 Litre
 - · Glove / Mask Box Holder, One Box
 - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- · Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

Paint - Grey White

| Description | CAM Lock with Wing Handle | Electronic Push Button Lock | Weight |
|----------------------------|---------------------------|-----------------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTMEDSM1/10/3S2D/CA | CTMEDSM1/10/3S2D/EPB | 43.2kg |
| 4 x 100mm, 1 x 150mm Trays | CTMEDSM1/10/4S1D/CA | CTMEDSM1/10/4S1D/EPB | 42.9kg |
| 6 x 100mm Trays | CTMEDSM1/10/6S/CA | CTMEDSM1/10/6S/EPB | 44.0kg |

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Carerail Mounted page 62
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71



+44 (0)1258 484455

Caretray Trolley - Plaster

Features

- Trollevs suitable for use within A&E / Plaster Room
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing travs from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Trolley pre-configured with the following accessories...
 - Bucket Support & Bracket
 - Stainless Steel Graduated Bucket, 15 Litre
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied - white, blue, green, yellow & red)
- Each trav supplied with 2 x long & 2 x short dividers
- Max load per trav 5kg
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 1005 x 550 x 993mm
- Plastic Top 915 x 455 (internal) at height 968mm
- Travs

Carts

ංඊ

Trolleys

Procedure

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Grey White

Options / Accessories

- Carerail Mounted page 62
- Laptop Holders page 68
- · Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71



| Description | Item | Weight |
|----------------------------|------------------|--------|
| 6 x 100mm, 4 x 150mm Trays | CTPLASM2/10/6S4D | 50.6kg |
| 8 x 100mm, 2 x 150mm Trays | CTPLASM2/10/8S2D | 49.9kg |

Caretray Trolley - Vaccination

Features

- Trolleys suitable for use with vaccination
- · Quick & easy disassembly, assists thorough & deep cleaning
- · Fully welded mild steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Trolley pre-configured with the following accessories...
 - Alcohol Gel Holder
 - Debris Bin & Bracket
 - · Glove, Mask Box Holder, One Box
 - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- · Each tray supplied with 2 x long & 2 x short dividers
- · Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm polymer castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Grey White

| Description | Item | Weight |
|----------------------------|------------------|--------|
| 3 x 100mm, 2 x 150mm Trays | CTVACIM1/10/3S2D | 32.0kg |
| 4 x 100mm, 1 x 150mm Trays | CTVACIM1/10/4S1D | 31.7kg |
| 6 x 100mm Trays | CTVACIM1/10/6S | 32.8kg |

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Defibrillator Platform page 61
- Carerail Mounted page 62
- Laptop Holders page 68
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71



Caretray Trolley - Basics

Features

- Standard level single column easy clean procedure trolleys
- Design allows for the trolley to be upgraded with side panels & trays
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame
- Recyclable polypropylene one-piece top...
 - Smooth surfaces & radiused edges to assist cleaning (compatible with steam cleaning)
 - Raised edges on three sides, providing defined grab points to act as push handles whilst an integral rim at the front contains spillages (max load 30kg)
- Solid grade laminate full width shelves, removable & supported on the trolley frame (max load per shelf 5kg)
- Carerails fitted to either side support a wide range of accessories
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Full Width Shelves 420 x 420 (internal) at heights 193 & 554mm

Standard Colours

Paint - Spearmint

Options / Accessories

- Overbridge & Tilt Bins page 60
- Defibrillator Platform page 61
- · Carerail Mounted page 62
- · Laptop Holders page 68
- Castors & Ring Buffers page 70
- Trays, Shelves, Dividers & Labels page 71





CTBASIM1/10/1FWS

CTBASIM1/10

| Description | Item | Weight |
|------------------------|------------------|--------|
| Top Panel, No Shelves | CTBASIM1/10 | 14.2kg |
| 1 x Full Width Shelf | CTBASIM1/10/1FWS | 15.5kg |
| 2 x Full Width Shelves | CTBASIM1/10/2FWS | 16.8kg |



Carts

Procedure Trolleys &

Procedure Trolleys & Carts

Mobile Workstation - Mild Steel, Low Level

Features

- Recessed liquid retaining top
- Single point locking using a CAM lock (supplied with two keys)
- Anti-tilt mechanism, allows only one drawer to be opened at any one time
- 85, 170 & 255mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Low level corner bumpers protect the fabric of the building
- 80mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall 615 x 500 x 810mm
- Drawers 480 x 375 x 85 / 170 / 255mm (internal)

Standard Colours

• Paint - Grey White



 Item
 Item

 MW005
 2 x 85mm, 1 x 170mm, 1 x 255mm Drawers

Weight 49.0kg



Manufacturing Medical Furniture & Equipment Since 1953

With over 70 years of experience, Bristol Maid is a leading manufacturer and supplier of medical furniture & equipment.

From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money.

We take pride in working in close partnership with our customers applying innovative design and technology to deliver effective solutions for today while anticipating the needs of tomorrow.



Hospital Metalcraft Limited Blandford Heights,Blandford Forum, Dorset. DT11 7TG, UK

Tel: +44 (0) 1258 451338 Fax: +44 (0) 1258 455056 sales@bristolmaid.com

www.bristolmaid.com

Supplied proudly in New Zealand by

TTTTT

WS Medical

Email: admin@wsmedical.co.nz Freephone: 0800 330 331

www.wsmedical.co.nz