


MODULAR
SYSTEM

Optima


Comfort and versatility
to suit every need



Semirigid Sacrolumbar
Back Support



Dorsolumbar
Back Support



Sacrolumbar Insert



Optima





PRIM



Discover



Optima

Semirigid Sacrolumbar Back Support

Created for greater comfort

High-frequency manufacture. Promotes discretion, providing perfection and modernity. Its soft breathable fabric and lumbar stays, together with a modern pattern, optimise the relationship between the comfort and containment that the patient needs.



HIGH-FREQUENCY TECHNOLOGY

Seamless manufacture


Dual straps for optimum immobilization and comfort.

Micro velcro for secure fastening.


Flexible, lightweight design: promotes comfort even when sitting down.



Hand loops: for easier fitting and adjustment



4 steel lumbar stays to ensure containment.



Includes lumbar plate for greater immobilization.

COLOURS



LOW BACK SUPPORT
HIGH BACK SUPPORT

	Front width	Back width	Ref.	XS	S	M	L	XL	XXL
LOW BACK SUPPORT	16 6.3"	26 10.2"	995/995G	75-82 27.5"-32.3"	82-90 32.3"-35.4"	90-100 35.4"-39.4"	100-110 39.4"-43.3"	110-120 43.3"-47.2"	120-130 47.2"-51.2"
HIGH BACK SUPPORT	20 7.9"	30 11.8"	996/996G	75-82 27.5"-32.3"	82-90 32.3"-35.4"	90-100 35.4"-39.4"	100-110 39.4"-43.3"	110-120 43.3"-47.2"	120-130 47.2"-51.2"

Dimensions: cm/inch. - Hip circumference. - Unisex.



Optima

Dorsolumbar Back Support

The evolution
of adaptability

The OPTIMA back support can be attached to the OPTIMA semirigid sacrolumbar back support to make a dorsolumbar orthosis. The quality of its 100% breathable fabric, together with its reinforcement stays, provides stability, compression and uniform support.

Padded shoulder straps that can be adjusted for a comfortable fit.

Central containment reinforcement with 2 adjustable steel stays.

Micro velcro for secure fastening in any position.

*Step 1:
trim the strap*



*Step 2:
insert the strap
into the closure*



Crocodile-type closures that enable the height of the strap to be adjusted.



COLOURS



EspOne



EspOneG

one
size




Optima

Sacrolumbar Insert

Maximum
immobilization



By combining the OPTIMA semirigid back support with the OPTIMA sacrolumbar insert, you can obtain a corset that maximises immobilization of the sacrolumbar area. This insert can be removed as treatment progresses, leaving the semirigid back support.




Possibility of moulding using heat for a perfect fit to the anatomy of the patient.

*Inner fabric:
Honeycomb*

*Outer fabric:
Lycra*



Soft 100% breathable fabric lining



Easy and safe fastening to the semirigid back support thanks to its micro velcro.

COLOURS



LOW
MODULE

HIGH
MODULE

Back width	XS	S-M	L	XL-XXL
30 11.8"	MDO-1/1G	MDO-2/2G	MDO-3/3G	MDO-4/4G
	75-82 29.5"-32.3"	82-100 32.3"-39.4"	100-110 39.4"-43.3"	110-130 43.3"-51.2"
37 14.6"	MDO-5/5G	MDO-6/6G	MDO-7/7G	MDO-8/8G
	75-82 29.5"-32.3"	82-100 32.3"-39.4"	100-110 39.4"-43.3"	110-130 43.3"-51.2"

Dimensions: cm/inch. - Hip circumference. - Unisex.



Full view of the **Optima** Modular System

Advanced manufacturing technology for a modern, comfortable and effective design.

Optima



Polígono Industrial, 1
Calle F, nº 15
28938 Móstoles - Madrid - Spain
+34 91 334 25 32

www.orthoprim.es
optima.prim.es



Watch the video of the
OPTIMA Modular System