


HOOD
PROSTHETIC SOCK

...making life comfortable



The company HOOD was founded in 1920 by William Hood, a WWI amputee.

It was a family run business until December of 1988 when
CAMP International Limited purchased the company HOOD.

In the spring of 1989, all production was relocated in Trenton, Ontario.

In 1997, Trulife purchased CAMP International.

Trulife is a manufacturer and market leader of silicone and gel-based products,
based out of Ireland since 1968.

The HOOD production area occupies approximately 5,000 square feet
of the 35,000 square feet expanded Trenton facility,
which now houses all manufacturing of HOOD and Trulife products
for the North American and Global markets.

Hundreds of socks and sheaths are manufactured every day,
using state of the art Shima-Seiki computerized knitting machines.

Our production staff has excellent knowledge and countless years of
experience in producing top quality prosthetic socks, sheaths and accessories.

Tel: 1.800.267.2812

Fax: 1.613.392.4139

Email: customerservice@trulife.com



...making life comfortable

TERMS

- All products are Net 30 days.
- All items are F.O.B. shipping point.
- Please ensure that all checks are made payable to Trulife.
Mail to 39 East Davis Street, Trenton, Ontario, K8V 4K8

DAMAGED SHIPMENT POLICY & PROCEDURE

- Please examine the container(s) upon receipt of goods. Note any damage on the delivery bill and have the driver sign, acknowledging the damage.
Contact the delivery service, and they will file an insurance claim. Then call Trulife to arrange for the shipment of replacement merchandise.
- Inspect the contents within 10 days. If you find shipping damage that was not obvious upon receipt, call the carrier immediately to request an inspection of both the container and merchandise. All claims for damage and shortage (container lost in transit) MUST BE MADE WITHIN TEN DAYS OF RECEIPT. If you have any questions, call Trulife.

RETURN POLICY

All merchandise has been carefully inspected and is manufactured from the highest quality materials; please examine your order thoroughly as soon as you receive it.

- Altered, discontinued, or special items are not returnable.
- All returns must be prepaid (collect shipments will be refused) and in resalable condition.
- A product must be returned within 30 days in order to receive credit.

RETURN PROCEDURE

If you need to return an item, please call Trulife for a return authorization number and state the reason why you wish to return the item.

To ensure fast service when you call, please have at hand the following:

- Date of purchase
- Invoice Number or Order Number
- The reason you wish to return the item

Returns are to be shipped prepaid to: Trulife

39 East Davis Street,
Trenton, Ontario,
K8V 4K8

Collect shipments will not be accepted.

Please clearly mark the outside of the package with the Return Authorization Number.

SPECIAL NOTES

Trulife has an ongoing goal of product improvement and we therefore reserve the right to alter product design, specification and prices without prior notice.

STANDARD SIZING CHART

STUMP SOCKS

Top Width	6" (15 cm)	7" (18 cm)	8" (20 cm)	9" (23 cm)	10" (25 cm)	11" (28 cm)	12" (30 cm)
Length	Toe	Toe	Toe	Toe	Toe	Toe	Toe
8" (20 cm)	3 ½" (9 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
10" (25 cm)	3 ½" (9 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
12" (30 cm)	3 ½" (9 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
14" (35 cm)	3 ½" (9 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
16" (40 cm)	3 ½" (9 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
18" (45 cm)	3 ½" (9 cm)	4" (10 cm)	4 ½" (11 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
20" (51 cm)	3 ½" (9 cm)	4" (10 cm)	4 ½" (11 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
22" (56 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
24" (61 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
26" (66 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
28" (71 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
30" (76 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
32" (81 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)
34" (86 cm)	3 ½" (9 cm)	4" (10 cm)	4" (10 cm)	5" (13 cm)	5" (13 cm)	6" (15 cm)	6" (15 cm)

SIZING GUIDE AND ORDERING INFORMATION

► HOW TO DETERMINE CORRECT HOOD SOCK SIZE

1. Measure limb circumference at the top of the prosthesis socket. Divide by 2; this will be the top width of the sock.*
2. Measure the length of the limb from the distal end to the top of the socket. Add 2" - 4" to allow for fit and/or shrinkage. This will be the length of the sock.*
3. Measure circumference 2" above the bottom of the residual limb. Divide by 2; This will be the toe measurement (Toe width of a HOOD sock is automatically determined by the size ordered).*

HOOD PROSTHETIC SOCKS HAVE THE FOLLOWING CHARACTERISTICS:

- 3 part peel-off label
 - ▶ Easy to file/reorder...One for the patient file, one for the patient to keep, and one for your referral source if required
- Color-coded stitching
 - ▶ Easy to identify ply
- Straight top
 - ▶ Non-restrictive
- Most socks available with DISTALHOLE
 - ▶ More choice for your patient
- MADE IN CANADA
 - ▶ All prosthetic socks and sheaths are made in Canada – support our Canadian workers

Each HOOD Prosthetic Sock is coded with a colored thread stitched at the top of the sock. Each color corresponds to the thickness of the sock or number of plies.

HOOD Prosthetic Sock COLOUR CODING SYSTEM

Beige _____	2 Ply
Yellow _____	3Ply
Pink _____	4Ply
Green _____	5Ply
Grey _____	6Ply

SPECIAL SOCKS

We use Shima Seiki electronic knitting machines for the manufacturing of HOOD prosthetic socks.

SPECIAL SOCKS

All special socks are 3 – 6 weeks delivery.

PRICING

Send us the complete requirements and measurements, for an accurate quote.

RETURN POLICY

By their nature, specials are unique, often requiring special handling. Please verify your measurements and specifications, as we do not allow returns on special socks.



ORDER FORM FOR SPECIAL SOCK

Account # _____ PO # _____ Phone # _____


Company Name _____ Date _____

This section is to be filled out by Trulife	
SPEC # _____	ITEM # _____
PRICE EACH _____	ORDER # _____

MATERIAL _____ PLY _____ QUANTITY _____

**ALL MEASUREMENTS DONE FLAT IN INCHES
TOP - WHOLE NUMBERS ONLY / LENGTH - WHOLE EVEN NUMBERS ONLY
TOE - REFER TO SIZING CHART ON PAGE 4**

- A _____ TOP
- B _____ TOE (2" ABOVE DISTAL END)
- C _____ LENGTH



INSTRUCTIONS _____

FAX # 613 392 4139
 PHONE 1 800 267 2812
 EMAIL: CUSTOMERSERVICE@TRULIFE.COM

WOOL CLASSIC

STUMP SOCKS				
Top Width	7" (18 cm)	8" (20 cm)	9" (23cm)	10" (25 cm)
3 ply	51037	51038	51039	51030
5 ply	51057	51058	51059	51050

Find style number in chart, then add the length in inches.
Available in lengths from 08" to 34" (20 - 86 cm).

NOT AVAILABLE WITH DISTAL HOLE

Order in multiples of six.

Lengths are in 2inch increments using even numbers only



- 100% pure wool
- ▶ Absorbant, resilient, mild elasticity
 - ▶ Precise fitting adjustments
 - ▶ Wide range available

SUPERWASH

STUMP SOCKS			
Top Width	7" (18 cm)	8" (20 cm)	9" (23cm)
3 ply	51737	51738	51739
5 ply	51757	51758	51759
6 ply	51767	51768	51769

Find style number in chart, then add the length in inches.
Available in lengths from 08" to 34" (20 - 86 cm).

DISTAL HOLE ONLY AVAILABLE IN SIZES STATED BELOW

3 ply: H1737, H1738, H1739

6 ply: H1767, H1768, H1769

Order in multiples of six.

Lengths are in 2inch increments using even numbers only



- 70% wool,
30% nylon
- ▶ Machine wash and dry
 - ▶ Precise fitting adjustments
 - ▶ Wide range available

COOL SOCK

95% polyester, 5% LYCRA®

COOLMAX fibers

Soft, lightweight, non-chafing and

breathable 2 and 3 Ply

- ▶ Machine wash and dry
- ▶ COOLMAX is great at wicking away perspiration
- ▶ Great for active wearers
- ▶ Precise fitting adjustments

LYCRA® and COOLMAX® are trademarks of INVISTA.

Order in multiples of six.

Top width	Thickness	6" – 7" (15 – 18 cm)	7" – 8" (18 – 20 cm)
Toe width		3" – 4" (7 – 10 cm)	4" – 5" (10 – 13 cm)
Length			
8" – 10" (20 – 25 cm)	2 ply	51922-01	51923-01
	3 ply	51932-01	51933-01
10" – 14" (25 – 35 cm)	2 ply	51922-02	51923-02
	3 ply	51932-02	51933-02
16" – 20" (40 – 51 cm)	2 ply	51922-03	51923-03
	3 ply	51932-03	51933-03
22" – 26" (55 – 66 cm)	2 ply	51922-04	51923-04
	3 ply	51932-04	51933-04



AVAILABLE WITH DISTAL HOLE

To order, substitute the first digit with an "H" (example: 51922-01 becomes H1922-01)

COOL SOCK II

SOFTER & MORE COMFORMABLE

95% polyester, 2% LYCRA®

COOLMAX® fibers

Soft, lightweight, non-chafing and breathable

3 and 5 ply

- ▶ Machine wash and dry
- ▶ COOLMAX® is great at wicking away perspiration
- ▶ Great for active wearers
- ▶ Wide range available

LYCRA® and COOLMAX® are trademarks of INVISTA.

Order in multiples of six.



Top width	Thickness	6" – 7" (15 – 18 cm)	6" – 7" (15 – 18 cm)	7" – 8" (18 – 20 cm)	8" – 9" (20 – 23 cm)	9" – 10" (22 – 25 cm)
Toe width		3" - 4" (7 - 10 cm)	4" - 5" (10 - 13 cm)	4" - 5" (10 - 13 cm)	5" - 6" (13 - 15 cm)	6" - 7" (15 - 18 cm)
Length						
8" – 10" (20 – 25 cm)	3 ply			51135-08	51136-08	
	5 ply	51153-08	51154-08	51155-08	51156-08	51157-08
10" – 12" (25 – 31 cm)	3 ply			51135-10	51136-10	
	5 ply	51153-10	51154-10	51155-10	51156-10	51157-10
12" – 14" (31 – 36 cm)	3 ply			51135-12	51136-12	
	5 ply	51153-12	51154-12	51155-12	51156-12	51157-12
16" – 18" (40 – 46 cm)	3 ply			51135-16	51136-16	
	5 ply	51153-16	51154-16	51155-16	51156-16	51157-16
20" – 22" (51 – 56 cm)	3 ply			51135-20	51136-20	
	5 ply	51153-20	51154-20	51155-20	51156-20	51157-20
24" – 26" (61 – 66 cm)	3 ply			51135-24	51136-24	
	5 ply	51153-24	51154-24	51155-24	51156-24	51157-24

AVAILABLE WITH DISTAL HOLE

To order, substitute the first digit with an "H" (example: 51135-08 becomes H1135-08)

HOOD SEAMLESS DIABETIC SOCK

A soft, seamless and comfortable diabetic sock for people at risk or suffering from diabetes, arthritis, or who simply have sensitive feet.



- 95% polyester, 5% LYCRA®
 - ▶ Machine wash and dry, no bleach or fabric softener
- COOLMAX® fibers
 - ▶ COOLMAX® is great at wicking away moisture for warmer and dryer feet
- Super soft and lightweight
 - ▶ Comfortable and non-irritating
- Non-chafing and breathable
 - ▶ Great for active wearers
- Treated with Ultra-Fresh
 - ▶ Inhibits the growth of odor-causing bacteria
- Non-restrictive cuff
 - ▶ Reduction of pressure and friction
- Seamless design
 - ▶ Ideal for those with diabetes, arthritis or sensitive feet

LYCRA® and COOLMAX® are trademarks of INVISTA.

SIZING CHART

Fits shoe size	02-Small	03-Medium	04-Large	05-X-Large
Women's	7-9	9.5-11	11.5-13	
Men's	6-8	8.5-10	10.5-13	12.5-14

Model	Description	Length	Color	02	03	04	05
51900	COOLMAX®	Crew	White	51900-02	51900-03	51900-04	51900-05
51901	COOLMAX®	Calf	White	51901-02	51901-03	51901-04	51901-05
51910	COOLMAX®	Crew	Black	51910-02	51910-03	51910-04	51910-05
51911	COOLMAX®	Calf	Black	51911-02	51911-03	51911-04	51911-05
51800	COOLMAX® With X-Static®	Crew	Grey	51800-02	51800-03	51800-04	51800-05
51801	COOLMAX® With X-Static®	Calf	Grey	51801-02	51801-03	51801-04	51801-05

SHEATH

100% nylon

11 sizes

Pocket or bottom seam

Wicks perspiration away from the skin, reduces friction between the limb and the prosthesis

- ▶ Stretches, machine wash and dry
- ▶ One size fits many
- ▶ Added protection from abrasion
- ▶ Amputee with poor skin or heavy perspiration will benefit from the use of a sheath



Order in multiples of six.

SIZING CHART

Length	Top Width			
	4" - 5" 10cm - 13cm	6" - 7" 15cm - 18cm	8" - 9" 20cm - 23cm	10" - 12" 25cm - 30cm
12" - 14" 30 - 35 cm	Extra Narrow Short	Narrow Short	Regular Short	Wide Short
16" - 18" 40 - 45 cm		Narrow Regular	Regular Regular	Wide Regular
20" - 24" 51 - 61 cm		Narrow Long	Regular Long	Wide Long
26" - 30" 66 - 76 cm			Regular Extra Long	

Sheath - Pocket Seam

Size	TOP WIDTH			
	Extra Narrow	Narrow	Regular	Wide
Short	54011-02	54012-02	54013-02	54014-02
Regular		54012-03	54013-03	54014-03
Long		54012-04	54013-04	54014-04
Extra Long			54013-05	

Sheath - Bottom Seam

Size	TOP WIDTH			
	Extra Narrow	Narrow	Regular	Wide
Short	54021-02	54022-02	54023-02	54024-02
Regular		54022-03	54023-03	54024-03
Long		54022-04	54023-04	54024-04
Extra Long			54023-05	

FILLER CAST SOCKS

50/50 polyester/cotton

► Comfortable

Lightweight or Heavy weight

► Superior strength and durability

May be used as a filler sock to compensate for changes in sock weight as required for proper functioning of the prosthesis.



FILLER CAST SOCKS – HEAVYWEIGHT & LIGHTWEIGHT

Model #	Width and Length	
FILLER CAST SOCK, HEAVYWEIGHT (Pink thread)		
54213	3" x 15"	7 x 38 cm
54214-12	3½" x 12"	9 x 30 cm
54214-19	3½" x 19"	9 x 48 cm
54215-12	5" x 12"	13 x 30 cm
54215-15	5" x 15"	13 x 38 cm
54215-19	5" x 19"	13 x 48 cm
54219-15	9" x 15"	23 x 38 cm
FILLER CAST SOCK, LIGHTWEIGHT (Blue thread)		
54224	3½" x 19"	9 x 49 cm
54225	5" x 19"	13 x 48 cm

Order in multiples of 12.

STOCKINETTE

Available in cotton or nylon. Packaged 25 yards per box

Size

Hospital Cotton	
Width	Model
2" (5cm)	56102
3" (7.5cm)	56103
4" (10cm)	56104
10" (25.5cm)	56110





Trulife Canada
39 East Davis Street
Trenton, Ontario, K8V 4K8
Tel: 1.800.267.2812 Fax: 1.613.392.4139
customerservice@trulife.com
www.trulife.com