

Child's Play Energy Lower Limb System ^{K3/K4}

- Composite keel with dynamic toe and cushioned heel response
- Lightweight, 4-bar, polycentric knee design
- Available in above knee and below knee configurations
- Child's Play Energy is featured on page 34


Indications: Medium to high impact levels, K3-K4

Component List

ID	BELOW KNEE	ID	ABOVE KNEE
A	AAAKSMP	A	AAAKSMP
B	AAAK210	B	AAAK211
C	AAAK100	K	SSK610A
D	AAAK135	D	AAAK135
E	AAAK285	E	AAAK285
	SFS24x (optional)		SFS24x (optional)
F	SEF132/144	F	SEF132/144

Weight Limit

BELOW KNEE	ABOVE KNEE
65 kg / 144 lb	55 kg / 121 lb



The System That Grows with Them!

- Build height adjusts from 92-220 mm (3.2-8.7 in)

