

Contracture Management Orthotics



**Comprehensive Selection of
Orthotic Therapy Devices**

Head & Neck

DynaPro™ Torticollis Orthosis

Ideal “static – dynamic” therapy to reverse mild to severe Torticollis. Provides abnormal reflex arc neuro therapy for patients with spasticity.

- Treats lateral flexion of the head and neck.
- Bend to fit shoulder cuff which allows for quick and easy customisation for creating a tailored product for the patient.
- Semi-rigid heat mouldable upright with ‘flex’ technology, which allows for a gentle stretch to reverse tissue shortening or manage lateral cervical flexion spasticity

Indications:

Tissue shortening / Torticollis



Selection		
Model No. Right	Model No. Left	Size
DPTO-RR	DPTO-RL	Adult Regular
DPTO-AR	DPTO-AL	Adult Large
DPTO-YR	DPTO-YL	Youth
DPTO-CR	DPTO-CL	Child

Specifications	
DynaPro™ Torticollis Orthosis	
Forehead to shoulder height	Adult Regular: 15-18cm / 6-7" Adult Large: 18-23cm / 7-9" Youth: 13-15cm / 5-6" Child: 10-13cm / 4-5"
Warranty	90 Days

DynaPro™ Flex Elbow

Ideal “static – dynamic” therapy to reverse mild to moderate tissue shortening contractures of the elbow (< 90°)

- Treats mild/moderate flexion.
- Abnormal reflex arc therapy provided by the device assists with neuro-rehabilitation.
- Axial hinge with legible goniometer facilitates brace fitting.
- Antimicrobial material which helps to prevent odour and infection.
- Gel elbow pad included.

Indications:

Tissue shortening – mild to severe tone, fair to good end feel, neurological hypertonicity/spasticity, neuro contracture management.



Selection	
Model No.	Size
DPFE20724-00	Child
DPFE20724-01	Youth
DPFE20724-02	Small
DPFE20724-03	Medium
DPFE20724-04	Large

Specifications	
DynaPro™ Flex Elbow	
Upper Arm Circumference	Child: 13-18cm / 5-7" Youth: 18-25cm / 7-10" Small: 20-25cm / 8-10" Medium: 24-29cm / 9.5-11.5" Large: 28-35.5cm / 11-14"
Lower Arm Circumference	Child: 9-16.5cm / 3.5-6.5" Youth: 13-18cm / 5-7" Small: 13-20cm / 5-8" Medium: 19-25cm / 7.5-10" Large: 23-30cm / 9-12"
Warranty	90 Days

SoftPro™ Static Gel Elbow

Ideal for providing positioning therapy for mild to moderate (< 90°) flexion contractures of the elbow. Double upright with single setting hinges which are set to provide the desired extension of the elbow. Rigid cuffs on bicep and forearm in soft cover; adjustable elbow cap.

- Static Therapy – Support, Protect and Position.
- Progressive Extension Stretch Therapy.
- Easily adjustable single setting hinge.
- Provides 3 point progressive extension stretch.

Indications:

RA-pain, inflammation, deformity, RSD-pain, with movement or stretch, tissue shortening, mild tone, poor end feel, fair end feel, poor skin integrity.



Selection	
Model No.	Size
SP503-02	Small
SP503-03	Medium
SP503-04	Large

Specifications	
SoftPro™ Static Gel Elbow	
Upper Arm Circumference	Small: 20-25cm / 8-10" Medium: 25-30cm / 10-12" Large: 30-35.5cm / 12-14"
Lower Arm Circumference	Small: 18-23cm / 7-9" Medium: 23-28cm / 9-11" Large: 28-33cm / 11-13"
Warranty	90 Days

Wrist & Hand

AirPro™ Grip

- Orthotic therapy for moderate to severe wrist, hand, finger joint stiffness, or contractures.
- Air Therapy - Low Load Prolonged Stretch.
- Provides gentle air stretch to reverse tissue shortening.
- Gently opens even the most severely non-fixed contracted hands.
- Ideal for spasticity management.
- Easy to apply and remove.

Indications:

Moderate to severe wrist, hand, finger joint stiffness, or contractures, tissue shortening of the MCP/PIP joints.



Selection		
Model No. Right	Model No. Left	Size
AP2-31AG-R	AP2-31AG-L	One Size

Specifications	
AirPro™ Grip	
Warranty	90 Days

AirPro™ Air Graduate WHFO

- Conforms to the shape of the hand, accommodating different degrees of finger contracture.
- Ideal for severely contracted fingers.
- The removable cone shaped palmer air bladder can be used independently of a wrist support.
- Antimicrobial material which helps to prevent odour and infection.
- The colour coded straps facilitate easy application of the brace.
- Extra cover is included.

Indications:

Tissue shortening, mild tone, fair to good end feel.



Selection	
Model No.	Size
AP31AG	Adult L/R
AP31AG-Y	Youth L/R
AP31AG-C	Child L/R

Specifications	
AirPro™ Air Graduate WHFO	
Width of MP Joints	Child: 4.5-5cm / 1.75-2" Youth: 5-6.5cm / 2-2.5" Adult: 9-10cm / 3.5-4"
Warranty	90 Days

AirPro™ Bladder Only

- Orthotic therapy for moderate to severe wrist, hand, finger joint stiffness, or contractures.
- Air Therapy - Low Load Prolonged Stretch.
- Provides gentle air stretch to reverse tissue shortening.
- Gently opens even the most severely non-fixed contracted hands.
- Ideal for spasticity management.
- Easy to apply and remove.

Indications:

Tissue shortening, mild tone, fair to good end feel.



Selection	
Model No.	Size
AP31AGHO	Adult

Specifications	
AirPro™ Bladder Only	
Warranty	90 Days

DynaPro™ Flex Hand

- Reverses moderate/severe flexion and/or spasticity of wrist/hand/fingers.
- Two incremental larger palmar rolls included to gradually increase extension of the MCP and PIP joints.
- Provides gentle ulnar drift correction.
- Heat mouldable/adjustable.
- Antimicrobial material which helps to prevent odour and infection.
- Colour coded straps to facilitate easy application of the brace.
- Extra cover included.

Indications:

Tissue shortening, mild to severe tone, fair to good end feel.



Selection		
Model No. Right	Model No. Left	Size
DPWF20721-R	DPWF20721-L	Adult

Specifications	
DynaPro™ Flex Hand	
Width of MP Joints	Adult: 9-10cm / 3.5-4"
Length	Adult: 33.5cm / 13.25"
Warranty	90 Days

Wrist & Hand

DynaPro™ Finger Flex

- A dynamic resting hand designed with a flexible finger platform.
- Flexible wrist, MCP and PIP joints.
- Thumb socket with two gel inserts to help treat thumb abduction.
- Finger separators included.
- Heat mouldable/adjustable.
- Antimicrobial material which helps to prevent odour and infection.
- Colour coded straps to help with easy application of the brace.
- Extra cover included.

Indications:

Mild to moderate contractures of the wrist, hand and fingers. To treat thumb abduction.



Selection		
Model No. Right	Model No. Left	Size
DPWFFR-00	DPWFFL-00	Child
DPWFFR-01	DPWFFL-01	Youth
DPWFFR-02	DPWFFL-02	Small
DPWFFR-03	DPWFFL-03	Medium
DPWFFR-04	DPWFFL-04	Large

Specifications	
DynaPro™ Finger Flex	
Width of MP Joints	Child: 4.5-5cm / 1.75-2" Youth: 5-6.5cm / 2-2.5" Small: 6.5-7.5cm / 2.5-3" Medium: 7.5-9cm / 3-3.5" Large: 9-10cm / 3.5-4"
Length	Child: 21cm / 8.25" Youth: 26cm / 10.25" Small: 28.5cm / 11.25" Medium: 31cm / 12.25" Large: 33.5cm / 13.25"
Warranty	90 Days

DynaPro™ Resting Hand (Thumb Ease)

- Resting hand designed with a rigid finger platform.
- Thumb pocket with 2 gel inserts to help treat thumb abduction.
- Flexible wrist minimises spastic episodes.
- Maintains functional hand alignment.
- Heat mouldable/adjustable.
- Antimicrobial material which helps prevent odour and infection.
- Colour coded straps facilitate easy brace application.
- Extra covers included.

Indications:

Tissue shortening, mild to severe tone to good end feel. Neurological hypertonicity/spasticity, neuro contracture management.



Selection		
Model No. Right	Model No. Left	Size
DPWR20737-02R	DPWR20737-02L	Small
DPWR20737-03R	DPWR20737-03L	Medium
DPWR20737-04R	DPWR20737-04L	Large

Specifications	
DynaPro™ Resting Hand (Thumb Ease)	
Width of MP Joints	Small: Up to 8cm / 2.5-3" Medium: 8-9cm / 3-3.5" Large: 9-10cm / 3.5-4"
Length	Small: 28.5cm / 11.25" Medium: 31cm / 12.25" Large: 33.5cm / 13.25"
Warranty	90 Days

SoftPro™ Champ Resting Hand

- Severely abducted thumb rests on the finger platform.
- The thumb can be positioned under the finger platform as a second step in abducting the thumb.
- Large finger platforms provide a comfortable ‘pan’ to support and protect the hand and fingers.
- Antimicrobial material which helps to prevent odour and infection.
- Colour coded straps to facilitate easy brace application.
- An extra cover is included with every SoftPro™ wrist/hand/finger Orthosis.

Indications:

RA-Pain, inflammation, deformity, RSD-pain with movement or stretch, tissue shortening, mild tone, poor end feel and poor skin integrity.



Selection		
Model No. Right	Model No. Left	Size
SP31CHAR-01	SP31CHAL-01	X-Small
SP31CHAR-02	SP31CHAL-02	Small
SP31CHAR-03	SP31CHAL-03	Medium
SP31CHAR-04	SP31CHAL-04	Large

Specifications	
SoftPro™ Champ Resting Hand	
Width of MP Joints	X-Small: Up to 6-7cm / 2.5" Small: 7-8cm / 2.5-3" Medium: 8-9cm / 3-3.5" Large: 9-10cm / 3.5-4"
Length	X-Small: 25cm / 10" Small: 28cm / 11" Medium: 30.5cm / 12" Large: 33cm / 13"
Warranty	90 Days

SoftPro™ Functional Resting Hand

- Treats moderate contractures of the wrist/hand/fingers.
- Extra long wrist strap holds the wrist in the correct position and applies a gentle stretch.
- Thumb is “Bend to Fit” to support, position and/or correct alignment of the thumb.
- Antimicrobial material which helps to prevent odour and infection.
- Colour coded straps facilitate brace application.
- Extra cover included with every SoftPro™ wrist/hand/finger Orthosis.

Indications:

RA-pain, inflammation, deformity, RSD- pain with movement or stretch, tissue shortening, mild tone, poor end feel, fair end feel, poor skin integrity.



Selection		
Model No. Right	Model No. Left	Size
SP31FFRR-00	SP31FFRL-00	Child
SP31FFRR-01	SP31FFRL-01	Youth
SP31FFRR-02	SP31FFRL-02	Small
SP31FFRR-03	SP31FFRL-03	Medium
SP31FFRR-04	SP31FFRL-04	Large

Specifications	
SoftPro™ Functional Resting Hand	
Width of MP Joints	Child: 4.5-5cm / 1.75-2" Youth: 5-6.5cm / 2-2.5" Small: 6.5-7.5cm / 2.5-3" Medium: 7.5-9cm / 3-3.5" Large: 9-10cm / 3.5-4"
Length	Child: 18cm / 7" Youth: 20cm / 8" Small: 23cm / 9" Medium: 25cm / 10" Large: 28cm / 11"
Warranty	90 Days

Wrist & Hand

SoftPro™ Grip

- Treats moderate to severe contractures of the wrist/hand/fingers.
- Palmar bar comfortably fits into the palm of a fist.
- Two incrementally larger palmar rolls included to gradually increase extension of the MCP and PIP joints.
- Bend to fit, easily moulded without tools yet sturdy enough to maintain shape during wear.
- Antimicrobial material which helps to prevent odour and infection.
- Colour coded straps to facilitate easy brace application.
- Finger separators and Ulnar Drift Strap available separately to allow custom fitting.
- Extra cover included with every SoftPro™ wrist/hand/finger Orthosis.

Indications:

RA-pain, inflammation, deformity, RSD-pain with movement or stretch, tissue shortening, mild tone, poor end feel, fair end feel, poor skin integrity.



Selection		
Model No. Right	Model No. Left	Size
SP31FGR-00	SP31FGL-00	Child
SP31FGR-01	SP31FGL-01	Youth
SP31FGR-02	SP31FGL-02	Small
SP31FGR-03	SP31FGL-03	Medium
SP31FGR-04	SP31FGL-04	Large

Specifications	
SoftPro™ Grip	
Width of MP Joints	Child: 6-7cm / 1.75-2" Youth: 7-8cm / 2-2.5" Small: Up to 8cm / 2.5-3" Medium: 8-9cm / 3-3.5" Large: 9-10cm / 3.5-4"
Length	Child: 18cm / 7" Youth: 20cm / 8" Small: 23cm / 9" Medium: 25cm / 10" Large: 28cm / 11"
Warranty	90 Days

SoftPro™ Palmar Resting Hand

- Severely abducted thumb rests on the finger platform.
- The thumb can be positioned under the finger platform as a second step in abducting the thumb.
- Large finger platforms provide a comfortable 'pan' to support and protect the hand and fingers.
- Antimicrobial material which helps to prevent odour and infection.
- Colour coded straps to facilitate easy brace application.
- An extra cover is included with every SoftPro™ wrist/hand/finger Orthosis.

Indications:

RA-pain, inflammation, deformity, RSD-pain with movement or stretch, tissue shortening, mild tone, poor end feel and poor skin integrity.



Selection		
Model No. Right	Model No. Left	Size
SP31FPRR-02	SP31FPRL-02	Small
SP31FPRR-03	SP31FPRL-03	Medium
SP31FPRR-04	SP31FPRL-04	Large

Specifications	
SoftPro™ Palmar Resting Hand	
Width of MP Joints	Small: Up to 8cm / 2.5-3" Medium: 8-9cm / 3-3.5" Large: 9-10cm / 3.5-4"
Length	Small: 23cm / 9" Medium: 25cm / 10" Large: 28cm / 11"
Warranty	90 Days

SoftPro™ Finger Separators

- Designed to give finger separation and mild abduction.
- The fingers are held in position with padded finger separators with ample padding for a comfortable fit.
- One size fits all.

Indications:

Finger Separator for SoftPro and DynaPro Wrist Orthoses (Sold Separately)



Selection	
Model No.	Size
SPFSS	One Size

Specifications	
SoftPro™ Finger Separators	
Warranty	90 Days

OrthoPro™ HyperEx Knee Brace

- Inflatable air bladder incorporated being the hyperextension control strap prevents migration and allows for adjustment of the knee position.
- Quick release straps are suitable for patients who require one-handed application.
- Lightweight construction and slimline profile ensure the brace conforms to the leg anatomy and offers superior comfort.

Indications:

Knee hyperextension and genu recurvatum.



Selection			
Model No. Right	Model No. Left	Size	Knee Circumference*
OP403-0-R	OP403-0-L	Small	25-30cm / 10-12"
OP403-1-R	OP403-1-L	Medium	33-35cm / 13-14"
OP403-2-R	OP403-2-L	Large	38-41cm / 15-16"

*Measure 10cm / 4" above knee center

Specifications	
OrthoPro™ HyperEx Knee Brace	
Thigh Circumference	Small: 23-33cm / 9-13" Medium: 35.5-48cm / 14-19" Large: 50-63.5cm / 20-25"
Length	Small: 25.5-35.5cm / 10-14" Medium: 35.5-51cm / 14-20" Large: 51-68.5cm / 20-27"
Warranty	90 Days

OrthoPro™ Passive Night Splint

- Comfortable night time solution for plantar fasciitis.
- Provides passive stretching for plantar fascia, calf muscles and Achilles tendon.
- Non-skid bottom for increased safety.
- Pre-set with 5 degree of dorsiflexion.
- Toe wedge included.
- Antimicrobial material which helps to prevent odour and infection.

Indications:

Post surgery/injury and pain with movement.



Selection		
Model No.	Size	Foot Length
OP800-1	Small	15-20cm / 6-8"
OP800-2	Medium	20-25cm / 8-10"
OP800-3	Large	25-28cm / 10-11"

Specifications	
OrthoPro™ Passive Night Splint	
Max Calf Circumference	Small: 38cm / 15" Medium: 43cm / 17" Large: 48cm / 19"
Warranty	90 Days

OrthoPro™ Stretch Night Splint

- Adjustable dorsiflexion straps provide optimal stretch.
- Non-skid bottom for increased safety.
- Toe wedge included.
- Antimicrobial material which helps prevent odour and infection.

Indications:

Post injury/surgery or pain with movement.



Selection		
Model No.	Size	Shoe Size
OP800-4	Small	4-6 / 37-39
OP800-5	Medium	7-9 / 40.5-43
OP800-6	Large	>10 / 44.5

Specifications	
OrthoPro™ Stretch Night Splint	
Shoe Size	Small: 4-6 / 37-39 Medium: 7-9 / 40.5-43 Large: >10 / 44.5
Warranty	90 Days

DynaPro™ AFO

- Antimicrobial material which helps prevent odour and infection.
- Effectively reverses ankle/foot contractures.
- Elastic pull straps effectively treat plantar flexion, inversion or eversion.
- Non-skid bottom for safer transfers.

Indications:

Tissue shortening, mild to severe tone, fair to good end feel, neurological hypertonicity/spasticity and neuro contracture management.



Selection		
Model No.	Size	Foot Length
DPAFO-00	Child	9-13cm / 3.5-5"
DPAFO-01	Youth	13-15cm / 5-6"
DPAFO-02	Small	15-20cm / 6-8"
DPAFO-03	Medium	20-25cm / 8-10"
DPAFO-04	Large	25-30cm / 10-12"

Specifications	
DynaPro™ AFO	
Max Calf Circumference	Child: 30.5cm / 12" Youth: 35.5cm / 14" Small: 38cm / 15" Medium: 43cm / 17" Large: 48cm / 19"
Warranty	90 Days

SoftPro™ In-Bed AFO

- Adjustable dorsiflexion straps provide optimal stretch.
- Non-skid bottom for increased safety.
- Toe wedge included.
- Antimicrobial material which helps prevent odour and infection.

Indications:

Tissue shortening, mild to severe tone, fair to good end feel, neurological hypertonicity/spasticity and neuro contracture management.




Selection		
Model No.	Size	Foot Length
SP10S-00	Child	9-13cm / 3.5-5"
SP10S-01	Youth	13-15cm / 5-6"
SP10S-02	Small	15-20cm / 6-8"
SP10S-03	Medium	20-25cm / 8-10"
SP10S-04	Large	25-30cm / 10-12"
SP10S-05	X-Large	28-31cm / 11-12"
SP10F-07	Bariatric (Fleece Lining)	>31cm / 12"

Specifications	
SoftPro™ In-Bed AFO	
Max Calf Circumference	Child: 30.5cm / 12" Youth: 35.5cm / 14" Small: 38cm / 15" Medium: 43cm / 17" Large: 48cm / 19" X-Large: 58cm / 23" Bariatric: Over 58cm / 23"
Warranty	90 Days

Head & neck			
Model No	Description	Size	Page
DPTO-RR	DynaPro™ Torticollis Orthosis	Adult Regular Right	2
DPTO-AR	DynaPro™ Torticollis Orthosis	Adult Large Right	2
DPTO-YR	DynaPro™ Torticollis Orthosis	Youth Right	2
DPTO-CR	DynaPro™ Torticollis Orthosis	Child Right	2
DPTO-RL	DynaPro™ Torticollis Orthosis	Adult Regular Left	2
DPTO-AL	DynaPro™ Torticollis Orthosis	Adult Large Left	2
DPTO-YL	DynaPro™ Torticollis Orthosis	Youth Left	2
DPTO-CL	DynaPro™ Torticollis Orthosis	Child Left	2
Elbow			
DPFE20724-00	DynaPro™ Flex Elbow	Child	3
DPFE20724-01	DynaPro™ Flex Elbow	Youth	3
DPFE20724-02	DynaPro™ Flex Elbow	Small	3
DPFE20724-03	DynaPro™ Flex Elbow	Medium	3
DPFE20724-04	DynaPro™ Flex Elbow	Large	3
SP503-02	SoftPro™ Static Gel Elbow	Small	3
SP503-03	SoftPro™ Static Gel Elbow	Medium	3
SP503-04	SoftPro™ Static Gel Elbow	Large	3
Wrist & Hand			
AP2-31AG-R	AirPro™ Grip	One Size Right	4
AP2-31AG-L	AirPro™ Grip	One Size Left	4
AP31AG	AirPro™ Air Graduate WHFO	Adult	4
AP31AG-Y	AirPro™ Air Graduate WHFO	Youth	4
AP31AG-C	AirPro™ Air Graduate WHFO	Child	4
AP31AGHO	AirPro™ Bladder Only	Adult	5
DPWF20721-R	DynaPro™ Flex Hand	Adult Right	5
DPWF20721-L	DynaPro™ Flex Hand	Adult Left	5
DPWFFR-00	DynaPro™ Finger Flex	Child Right	6
DPWFFR-01	DynaPro™ Finger Flex	Youth Right	6
DPWFFR-02	DynaPro™ Finger Flex	Small Right	6
DPWFFR-03	DynaPro™ Finger Flex	Medium Right	6
DPWFFR-04	DynaPro™ Finger Flex	Large Right	6
DPWFFL-00	DynaPro™ Finger Flex	Child Left	6
DPWFFL-01	DynaPro™ Finger Flex	Youth Left	6
DPWFFL-02	DynaPro™ Finger Flex	Small Left	6
DPWFFL-03	DynaPro™ Finger Flex	Medium Left	6
DPWFFL-04	DynaPro™ Finger Flex	Large Left	6
DPWR20737-02R	DynaPro™ Resting Hand (Thumb Ease)	Small Right	6
DPWR20737-03R	DynaPro™ Resting Hand (Thumb Ease)	Medium Right	6
DPWR20737-04R	DynaPro™ Resting Hand (Thumb Ease)	Large Right	6
DPWR20737-02L	DynaPro™ Resting Hand (Thumb Ease)	Small Left	6
DPWR20737-03L	DynaPro™ Resting Hand (Thumb Ease)	Medium Left	6
DPWR20737-04L	DynaPro™ Resting Hand (Thumb Ease)	Large Left	6
SP31CHAR-01	SoftPro™ Champ Resting Hand	X-Small Right	7
SP31CHAR-02	SoftPro™ Champ Resting Hand	Small Right	7
SP31CHAR-03	SoftPro™ Champ Resting Hand	Medium Right	7
SP31CHAR-04	SoftPro™ Champ Resting Hand	Large Right	7
SP31CHAL-01	SoftPro™ Champ Resting Hand	X-Small Left	7
SP31CHAL-02	SoftPro™ Champ Resting Hand	Small Left	7
SP31CHAL-03	SoftPro™ Champ Resting Hand	Medium Left	7
SP31CHAL-04	SoftPro™ Champ Resting Hand	Large left	7

Wrist & Hand Continued			
Model No	Description	Size	Page
SP31FFRR-00	SoftPro™ Functional Resting Hand	Child Right	7
SP31FFRR-01	SoftPro™ Functional Resting Hand	Youth Right	7
SP31FFRR-02	SoftPro™ Functional Resting Hand	Small Right	7
SP31FFRR-03	SoftPro™ Functional Resting Hand	Medium Right	7
SP31FFRR-04	SoftPro™ Functional Resting Hand	Large Right	7
SP31FFRL-00	SoftPro™ Functional Resting Hand	Child Left	7
SP31FFRL-01	SoftPro™ Functional Resting Hand	Youth Left	7
SP31FFRL-02	SoftPro™ Functional Resting Hand	Small Left	7
SP31FFRL-03	SoftPro™ Functional Resting Hand	Medium Left	7
SP31FFRL-04	SoftPro™ Functional Resting Hand	Large left	7
SP31FGR-00	SoftPro™ Grip	Child Right	8
SP31FGR-01	SoftPro™ Grip	Youth Right	8
SP31FGR-02	SoftPro™ Grip	Small Right	8
SP31FGR-03	SoftPro™ Grip	Medium Right	8
SP31FGR-04	SoftPro™ Grip	Large Right	8
SP31FGL-00	SoftPro™ Grip	Child Left	8
SP31FGL-01	SoftPro™ Grip	Youth Left	8
SP31FGL-02	SoftPro™ Grip	Small Left	8
SP31FGL-03	SoftPro™ Grip	Medium Left	8
SP31FGL-04	SoftPro™ Grip	Large Left	8
SP31FPRR-02	SoftPro™ Palmar Resting Hand	Small Right	8
SP31FPRR-03	SoftPro™ Palmar Resting Hand	Medium Right	8
SP31FPRR-04	SoftPro™ Palmar Resting Hand	Large Right	8
SP31FPRL-02	SoftPro™ Palmar Resting Hand	Small Left	8
SP31FPRL-03	SoftPro™ Palmar Resting Hand	Medium Left	8
SP31FPRL-04	SoftPro™ Palmar Resting Hand	Large Left	8
SPFSS	SoftPro™ Finger Separators	One Size	9
Knee			
OP403-0-R	OrthoPro™ HyperEx Knee Brace	Small Right	10
OP403-1-R	OrthoPro™ HyperEx Knee Brace	Medium Right	10
OP403-2-R	OrthoPro™ HyperEx Knee Brace	Large Right	10
OP403-0-L	OrthoPro™ HyperEx Knee Brace	Small Left	10
OP403-1-L	OrthoPro™ HyperEx Knee Brace	Medium Left	10
OP403-2-L	OrthoPro™ HyperEx Knee Brace	Large left	10
Ankle			
OP800-1	OrthoPro™ Passive Night Splint	Small	11
OP800-2	OrthoPro™ Passive Night Splint	Medium	11
OP800-3	OrthoPro™ Passive Night Splint	Large	11
OP800-4	OrthoPro™ Stretch Night Splint	Small	11
OP800-5	OrthoPro™ Stretch Night Splint	Medium	11
OP800-6	OrthoPro™ Stretch Night Splint	Large	11
DPAFO-00	DynaPro™ AFO	Child	12
DPAFO-01	DynaPro™ AFO	Youth	12
DPAFO-02	DynaPro™ AFO	Small	12
DPAFO-03	DynaPro™ AFO	Medium	12
DPAFO-04	DynaPro™ AFO	Large	12
SP10S-00	SoftPro™ In-Bed AFO	Child	12
SP10S-01	SoftPro™ In-Bed AFO	Youth	12
SP10S-02	SoftPro™ In-Bed AFO	Small	12
SP10S-03	SoftPro™ In-Bed AFO	Medium	12
SP10S-04	SoftPro™ In-Bed AFO	Large	12
SP10S-05	SoftPro™ In-Bed AFO	X-Large	12
SP10F-07	SoftPro™ In-Bed AFO	Bariatric	12



Trulife
Tel 0114 2618100
Email: info_UK@trulife.com

www.trulife.com