

Hensinger Head Support



Revolutionary in the world of cervical collars

- Soft foam wrap-around collar designed to support the patient's head
- Used in conjunction with a harness
- Maintains mid-line head position, which can help improve breathing, swallowing and feeding
- Can help aid in social interaction and awareness
- Made from a waterproof foam which is coated in a flexible plastic
- Soft covers included

Unmounted Model	Size	Neck Circumference
60816000010	X-Small	25 cm
60816000011	Small	30 cm
60816000012	Medium	38 cm
60816000013	Large	43 cm