



# KINETIC LIGHT

## NOW AVAILABLE WITHOUT PYLON



### Kinetic Light w/o Pylon K2

Model No. Light SKF110, Dark SKF113

#### Features:

- K2 foot with integrated multi-axial ankle
- Proximal male pyramid
- Full length, composite keel
- Integrated ankle function for all categories is the same as the existing Kinetic
- Sandal toe foot shell (non-removable)
- 160kg weight limit
- Foot shells: Light, Dark



#### Indications:

Moderate impact level, K2, maximum weight of 160 kg /350 lb\*

**Suggested code:** L5972,L5986





# KINETIC LIGHT

## w/o Pylon K2



Model No.	Product Description	Size	Weight Limit
SKF110	Kinetic Light w/o Pylon Light Foot Shell ●	22 - 30cm	160kg / 350lb
SKF113	Kinetic Light w/o Pylon Dark Foot Shell ●	22 - 30cm	160kg / 350lb

Specifications	Product Description
Foot Shell	Sandal Toe (Light, Dark)
Product Weight (Cat 1-2, Size 25 without Pylon)	481g / 1.1lb
Product Weight (Cat 3-5, Size 25 without Pylon)	564g / 1.2lb
Min. Clearance	7.3cm / 2.9in
Heel Rise	9.5mm / 3/8in
Warranty Keel / Foot Shell	2 Years / 1 Year

**Preassembled with Foot Shell**

\*

Foot Length (cm) ▶	22	23	24	25	26	27	28	29	30
Patient Weight ▼	Cat	Cat	Cat	Cat	Cat	Cat	Cat	Cat	Cat
137 - 160kg / 301 - 350lb	-	-	-	-	5	5	5	5	5
101 - 136kg / 221 - 300lb	-	-	4	4	4	4	4	4	4
81 - 100kg / 177 - 220lb	-	2	3	3	3	3	3	3	-
< 80kg / < 176lb	1	1	2	2	2	2	2	-	-