

rce

POSTERIOR LATERAL STRUT AFO

## Matrix L-Pace

New

The NEW Matrix L-Pace offers Comfort, Durability and Style

#### **Features:**

- Posterior lateral strut of intermediate stiffness for best combination of support, flexibility and energy return.
- Made of aerospace grade composite materials for lightweight and durability.
- Ergonomic footplate with curved heel for better foot fitting.
- Flexible high strength band that easily conforms around calf muscle.
- Breathable padded cuff with reverse pull hook and loop compatible strap.

#### **Indications:**

Low to high activity patients, mild to moderate footdrop, with or without mild to moderate foot spasticity, with weak muscles or with loss of calf muscle, for unilateral to bilateral patients.

## Carbon Composite Dynamic AFO's

The Matrix Family has been developed to provide a dynamic response for patients with drop foot, limb proprioception deficit, mild proximal deficit as well as for partial foot amputations.

Trulife's newest inclusion to the Matrix Family, **Matrix L-Pace** is the perfect **posterior lateral** AFO addition.

- Stylish design
- Fits a wide range of footwear
- Lightweight
- Durable
- 2 year warranty





### Matrix L-Pace

Selection						
Model No. Right	Model No. Left	Size	Height	Foot Length	Shoe Size US	Shoe Size EU
FH450R02	FH450L02	Small	36.5cm / 14.4in	<23.6cm / <9.3in	5/6	36/38
FH450R03	FH450L03	Medium	36.5cm / 14.4in	25cm / 9.8in	7/8	39/40
FH450R04	FH450L04	Large	38.5cm / 15.2in	26.4cm / 10.4in	9/10	41/43
FH450R05	FH450L05	X-Large	38.5cm / 15.2in	27.4cm / 10.8in	10>	43>

Specifications				
Matrix L-Pace				
Material	Carbon Composite			
Warranty	Hardware: 2 years Soft goods: 90 days			

# **Matrix Family Rigidity Chart**

