

New

MATRIX

L-Pace



**Posterior Lateral
Strut AFO**

MATRIX

Pace



**Posterior Medial
Strut AFO**

Carbon Composite Dynamic AFO's

The Matrix Family has been developed to provide a dynamic response for patients with drop foot, limb proprioception deficit, mild proximal deficit as well as for partial foot amputations.

Trulife's newest inclusions to the Matrix Family, **Matrix Pace & Matrix L-Pace** are the perfect **posterior lateral & medial** AFO additions.

Matrix Pace & L-Pace

Features:

- Posterior **Lateral (L-Pace)** / Posterior **Medial (Pace)** strut of intermediate stiffness for best combination of support, flexibility and energy return.
- Made of aerospace grade composite materials for lightweight and durability.
- Ergonomic footplate with curved heel for better foot fitting.
- Narrower footplate for greater variety of compatible footwear.
- Flexible high strength band that easily conforms around calf muscle.
- Breathable padded cuff with reverse pull hook and loop compatible strap.
- Flexible high strength band that easily conforms around calf muscles

Indications:

Low to high activity patients, mild to moderate footdrop, with or without mild to moderate foot spasticity, with weak muscles or with loss of calf muscle, for unilateral to bilateral patients.



- *Stylish design*
- *Fits a wide range of footwear*
- *Lightweight*
- *Durable*

New
MATRIX
L-Pace



Posterior Lateral
Strut AFO

Selection						
Model No. Right	Model No. Left	Size	Height	Foot Length	Shoe Size US	Shoe Size EU
FH450R02	FH450L02	Small	36.5cm / 14.4in	<23.6cm / <9.3in	5/6	36/38
FH450R03	FH450L03	Medium	36.5cm / 14.4in	25cm / 9.8in	7/8	39/40
FH450R04	FH450L04	Large	38.5cm / 15.2in	26.4cm / 10.4in	9/10	41/43
FH450R05	FH450L05	X-Large	38.5cm / 15.2in	27.4cm / 10.8in	10>	43>

MATRIX
Pace



Posterior Medial
Strut AFO

Selection						
Model No. Right	Model No. Left	Size	Height	Foot Length	Shoe Size US	Shoe Size EU
FH400R02	FH400L02	Small	36.5cm / 14.4in	<23.6cm / <9.3in	5/6	36/38
FH400R03	FH400L03	Medium	36.5cm / 14.4in	25cm / 9.8in	7/8	39/40
FH400R04	FH400L04	Large	38.5cm / 15.2in	26.4cm / 10.4in	9/10	41/43
FH400R05	FH400L05	X-Large	38.5cm / 15.2in	27.4cm / 10.8in	10>	43>

Specifications	
Matrix L-Pace/Matrix Pace	
Material	Carbon Composite
Warranty	Hardware: 2 years Soft goods: 90 days

Matrix Family
Rigidity Chart

