

**NEW ARRIVAL**

**Matrix Pace**

**Your Perfect**

**Fit AFO for**

**Drop Foot**

**Intermediate  
Stiffness**



The New Matrix Pace offers  
Comfort, Durability and Style  
with a 2 year warranty and great price.



[www.trulife.com](http://www.trulife.com)



# Matrix Pace

- 2 year warranty on hardware
- Shaped carbon and glass strut of intermediate stiffness
- Ergonomic footplate with curved heel for better foot fitting
- Narrower footplate for greater variety of compatible footwear
- Breathable padded cuff with reverse pull hook and loop compatible strap
- Flexible high strength Hytrel® band that easily conforms around calf muscles




**Indications:** Low to high activity patients with unilateral or bilateral drop foot

## Matrix Pace

Model No. Left	Model No. Right	Size	Height and Footplate Length (cm / in)	Shoe Size US	Shoe Size EU
FH400L02	FH400R02	Small	36.5 / 14.323.5 / 9.3	5 / 6	36 / 38
FH400L03	FH400R03	Medium	36.5 / 14.325.0 / 9.8	7 / 8	39 / 40
FH400L04	FH400R04	Large	38.5 / 15.126.5 / 10.4	9 / 10	41 / 43
FH400L05	FH400R05	X-Large	38.5 / 15.127.5 / 10.8	10+	43+

## Matrix Family Rigidity Chart



# Matrix

[www.trulife.com](http://www.trulife.com)