

New PECTUS PRO

Pectus Carinatum Orthosis

↓ New and improved design.



Wider size range to accommodate larger chests.



↓ Optional, soft, shoulder straps included for additional stabilization during wear



↑ ~2.5" of lateral adjustment for asymmetric Pectus Carinatum

Pectus Pro

Pectus Carinatum Orthosis

- Provides increased anterior/posterior directed corrective compression forces.
- New and improved design.
- Wider size range (74-104cm/29-41") to accommodate larger chests.
- New ratchet closure for greater durability.
- Optional, soft, shoulder straps included.
- 2.5" of lateral adjustment for asymmetric Pectus Carinatum
- Removable, washable and breathable pads for increased comfort.
- Slide adjustment bars and screws.
- One size fits most

Indications:

Abnormal chest shape, protrusion of both the sternum (breastbone) and the cartilage portion of the ribs, and bowed chest (Pigeon chest).

Selection		
Model No.	Size	Chest Circumference
PECTUS PRO	One Size Fits Most	74-104cm / 29-41in

Specifications	
Pectus Pro	
Material	Low Density Polyethylene/Aluminum/Padding /Latex Free
Warranty	90 Days

