



## ASKIR 30 SUCTION PUMP

NEW ASKIR 30 is a desk type electric suction unit suitable for nasal, oral and endotracheal aspiration of bodily fluids from adults and children. Indicated for tracheotomised patients, minor surgical applications and post operative therapy in home care. Main unit equipped with vacuum-meter (bar and kPa) and vacuum regulator. Complies with EN 60601-1-11 'for use in the home healthcare environment' and carries IP21 degree of protection.

### Specification:

Motor	Oiless and maintenance free piston pump
Power Feeding	230V-50Hz
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (Adjustable)	-0.88 Bar -80kPa -600 mmHG
Noise Level	60,5 dB
Max Free Air Flow Rate	40 l/min
Power Consumption	110 VA
Fuse	1 x F 1,6 A 250V
Duty Cycle	Non Stop Operation
IP Protection Code	IP21
Weight	3,6 KG
Size	35 x 21 x 18 cm
Place of Manufacture	Italy