

Combat CARL Care and Evacuation Manikin




Item no.: T3C

*Experts in medical education
since 1950*

EZ
ERLER ZIMMER

Distributed by;

Q **QUREMED**



Combat CARL offers all functions for TCCC and TECC training.


Life-size (183cm), 70 kg, with highly mobile extremities, no sensitive electronics, no gimmicks: Ideal for treatment and transport in extreme conditions. Suitable for practicing the standard algorithms.

CABCDE

- C** Three included deep wounds can be applied to 7 sites and connected to bleeding modules. Wound packing is possible for all wounds, as is the application of tourniquets.
- A** Airway management with oro- and nasopharyngeal as well as supraglottic devices.
- B** Thoracic relief in monaldi and Bülau position
- C** Hemostasis and access to the vascular system: i.v. in the forearm, i.o. sternal, humerus head and proximal tibia with multiple-use pads.
- D** Voice output from prerecorded files or directly.
- E** Robust and insensitive to water, dirt, cold and heat.

- Thorax relief Monaldi
- Thorax relief Bülau
- IV access possibility

- I.O Access FAST
- I.O Access



- Voice output via loudspeaker in the head by means of sound files or via microphone function
- Airway management oro - and nasopharyngeal

- Openings for wounds, IM puncture pads and blind covers
- Robust and insensitive to water, dirt and other adverse environmental conditions



Airway management



i.V. additions



Patient voice




Wounds in the thoracic region



Relief puncture according to Monaldi




Sternal intraosseous approach



Chest Drainage



Wound-Packing



TQ - placement

Distributed by;

