

PICO arterial syringes

PICO syringes are arterial syringes for pH, blood gas, oximetry, electrolyte and metabolite analysis

HIGH ANALYTICAL QUALITY

- Syringes are preheparinized with dry electrolyte-balanced heparin
- Preheparinized syringes reduce risk of clots, electrolyte bias and sample dilution

REDUCED RISK OF BLOOD SPILLAGE

- TIPCAP attaches securely to the syringes, reducing the risk of blood spillage
- Needle cube available for PICO70 syringes

PICO FEATURES

- The PICO 50 and PICO70 are compatible with an A-line catheter stopcoc
- The TIPCAP attaches to the syringe to reduce the risk of blood spillage
- PICO70 syringes allows for fast filling

PICO SYRINGES ARE AVAILABLE IN TWO MODELS:

- PICO50 aspirator for arterial line sampling
- PICO70 self-fill for arterial puncture sampling



Preheparinized syringe



TIPCAP



Needle comes in three sizes

Specifications	PICO50 Aspirator	PICO70 Self-fill
Heparin	80 IU electrolyte-balanced heparin, coated on cellulose fiber	60 IU electrolyte-balanced heparin, coated on cellulose fiber
Sample volume	0.5-2.0 mL	0.3-1.5 mL
Sterilization	By ethylene oxide	By ethylene oxide

Data subject to change without notice.

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS. © Radiometer Medical ApS, 2700 Brønshøj, Denmark, 2020. All Rights Reserved. MAPSS-000406 R1.