

**3 YEARS** PARTS & LABOUR WARRANTY



**VFT-28** 

## DISPLAY

Top - Current temperature Left - Min temperature Right - Max temperature

## Â,

## DIMENSIONS

Width: 98mm Depth: 22mm Height: 110mm

## **FEATURES**

- Exceptionally easy to view and record hi/lo readings
- Hi, lo and current temperatures always on display
- One button press to reset hi/lo readings
- Audible alarm when outside +2°C to +8°C range
- Alarm has one press deactivation
- Perfect temperature monitor for vaccine fridge or backup cooler
- No setup required
- Includes 2 sensors to measure air or vaccine temperature
- Designed and manufactured specifically for Australian and New Zealand vaccine fridges
- Fold out stand or hang from a hook
- 2x AAA batteries included