

**ADROIT**  
SURGICAL



**VIE SCOPE**  
**bébé Vie™**



Disposable, Single-Use  
Device



Clear Line of Sight



Patented LED Light  
Illumination of the  
Throat



Easy to Use

# One-Size Paediatric bébé Vie™

**Patented one-size paediatric bébé Vie Scope™ is a self-contained, battery-powered, disposable laryngoscope for intubating toddlers to age ten.**



**Best illumination of the vocal cords with minimal chance of light obstruction by secretions or blood**



**Users can intubate children with an age-appropriate endotracheal tube (ETT) directly using a 'sideways' insertion technique.**



**Users can insert the tip of the bébé Vie Scope™ either in the vallecula like a Macintosh laryngoscope**



**You can use the blade tip to elevate the epiglottis just like a Miller.**



**The bébé Vie Scope™ can also be used as a bougie introducer.**

