



HAWKSLEY COVER GLASS



Code	Type
Z5000	Standard
Z20088	Nageotte
Z10088	McMaster

Cover glass, also known as a coverslip, is a thin, glass slide that covers the specimen viewed on a microscope slide or haemocytometer. Hawksley cover glass is rigidly compliant with the BS 748 standard, ensuring the upmost of scientific accuracy. The cover glass support platforms are ground and polished optically flat and co-planar with ruled area to $\pm 0.001\text{mm}$.

Standard and Nageotte cover glass come in packs of 10, and the McMaster cover glass are purchased in single packs. Dimensions are 22x24x0.45mm thick.

www.hawksley.co.uk
enquiries@hawksley.co.uk
+44 (0)1903 752 815