


MatrixDerm™

Regenerative Collagen Dental Membrane






Product Features

- ▶ Highly purified intact collagen
- ▶ High mechanical strength and structural integrity for easy handling
- ▶ Highly biocompatible
- ▶ Permeable to macromolecules and nutrients while creating an effective barrier to epithelial cells
- ▶ Conformable and repositionable for precise adjustment and placement
- ▶ Resorbable in 6 – 9 months, eliminating the need for removal surgery
- ▶ Terminally sterilized and ready for use following brief hydration

Collagen Matrix | Dental
A Division of Collagen Matrix, Inc.

MatrixDerm™ porcine collagen membrane has been precisely enhanced to provide periodontal and dental surgeons with the ideal balance of properties to effectively address a host of clinical indications and surgical procedures.



MatrixDerm strikes the right balance among critical product characteristics to deliver predictable results and easily supports the clinician's preferred techniques.

MatrixDerm has been designed to simultaneously achieve drapability, structural integrity, and the ability to reposition the membrane for ideal placement prior to closure. MatrixDerm conforms to the defect site, but also has sufficient memory to maintain space for bone and tissue ingrowth.

The programmed resorption time of MatrixDerm supports bone or tissue regeneration for 6 to 9 months to allow for ingrowth and remodeling of the defect site. The semi-permeable membrane allows for nutrient exchange while providing a cell barrier to prevent epithelial down growth.



Ordering Information

Ordering Number	Size
PDM1520	15mm x 20mm
PDM2030	20mm x 30mm
PDM3040	30mm x 40mm

MatrixDerm™ Regenerative Collagen Dental Membrane is a nonfriable, conformable, resorbable, membrane matrix consisting of highly purified type I and III collagen derived from intact porcine dermis.

Collagen Matrix | Dental

A Division of Collagen Matrix, Inc.

15 THORNTON ROAD, OAKLAND, NJ 07436

T: 201-405-1477 • TOLL FREE: 888-405-1001

F: 201-405-1355 • WWW.COLLAGENMATRIXDENTAL.COM