

The ATB100 is a comfortable, washable, and reusable face mask made in Korea with unique UV and antibacterial properties. The special cross-section structure quickly releases sweat and dries quickly. The mask blocks harmful UV-B (290-320nm) and UV-A (320-400nm) and is designed with materials to provide all day comfort.



Model Number: ATB100

“CDC recommends that people wear a face covering to cover their nose and mouth in the community setting. This is an additional public health measure people should take to reduce the spread of Covid-19 in addition to (not instead of) social distancing, frequent hand cleaning and other everyday preventive actions.”

General Use Face Mask Information	
User	General Public
Risk Factor	Low Risk Situations
Design	Designed to contain coughs and sneezes and prevent disease transmission to others – NOT to protect the wearer from breathing in aerosols.
When To Wear	Recommended for use when a person can't perform social distancing. <ul style="list-style-type: none"> • When coughing or sneezing. • Using public transportation. • Shopping and working at essential businesses like grocery stores and pharmacies.
Fit Testing Required	No Scarves and bandanas can be used if necessary.
Use Limitations	Homemade/Cloth Masks: <ul style="list-style-type: none"> • Should be washed after each use. • Should not be worn damp or when wet from spit or mucus.



Features

- Washable and reusable (less waste)
- UV protection
- Very comfortable materials and fit
- Individually packed
- Comfortable earloop design
- Fast drying
- Easy to breathe through
- Various sizes
- Made in Korea

Intended Use

For use by the general public. Using a face mask may limit the spread of User's respiratory droplets.

Source: PA Department of Health



MADE IN KOREA

Materials
ATB Yarn 95% Polyester
Span Yarn 5% Polyurethane

Packaging	
Qty/ Box	Boxes/ Case
5	20


Fitting Instructions



Sanitize your hands thoroughly before touching the mask



Inspect the overall condition of the mask to ensure there are no defective parts



Unfold the mask with the inside facing up



Place the mask over your mouth and nose and place straps behind ears



Using both hands, press down on the mask to create the best possible fit



Adjust the mask as needed and fit check one last time



Excellent antibacterial performance



- It has excellent antibacterial effect by using safe silver ingredient approved by FDA.
- Unlike post-processed products, it shows excellent antibacterial performance even after washing.



Quick absorption and drying



- The special cross-section structure quickly releases sweat and dries quickly.
- It gives you the comfort of the summer when you sweat a lot.



UV protection



- It blocks not only UV-B (290 – 320nm) but also UV-A (320 – 400nm). Furthermore it even blocks ultraviolet rays that penetrate deep into the skin.



Warnings and Limitations

- Always confirm the product is suitable for the application
- Keep the mask on during all periods of exposure
- Fit testing should be performed before entering exposed areas
- Inspect all parts of the mask to ensure it is in good condition before use
- Dispose according to applicable regulations
- Do not alter the product in any way
- Masks should be stored in a cool and dry area