



SIZED ROM KNEE BRACE

SIZE	MEASUREMENT A RANGE	MEASUREMENT B RANGE
Small	Less than 18"	Less than 16"
Medium	18" - 21"	16" - 18"
Large	21" - 24"	18" - 21"
XL	24" - 27"	21" - 24"
XXL	27" - 30"	24" - 27"
XXXL	30" - 33"	27" - 30"
4XL	33" - 36"	30" - 33"
5XL	36" - 39"	33" - 36"

HOW TO CHOOSE THE RIGHT SIZE?



Measure the circumference around your leg at

- A** 6" above the middle of your kneecap
- B** 6" below the middle of your kneecap

Stand up while measuring.

Always use a soft/flexible tape measure.

Measurements will be more accurate if you have someone measure for you.

