



Building a Stronger America Together: Veteran-Led PPE, Small Business Allies.



100% made in the USA



Made with USA raw materials

4-PLY ASTM LEVEL 3 SURGICAL DISPOSABLE FACE MASK COMFORTABLE, LIGHTWEIGHT MATERIAL

SKUS:

- USAM350 (Box of 50)
- USAM32000 (Case of 2000)



FEATURES:

- Lightweight, non-woven polypropylene material design is soft on the face
- Elastic loops fit comfortably around ears
- Nose clip molds mask securely to face
- Fluid resistant and disposable, single use
- Not made with Natural Rubber Latex
- Mask Size 7" Wide x 3.625" Tall
- Surgical Grade Mask made with 100% USA Raw Materials, Non-Sterile

ASTM F2100-11	Level 3
Fluid Resistance, mmHg	160
BFE (Bacterial Filtration Efficiency)	>98%
PFE (Particulate Filtration Efficiency)	>98%
Differential Pressure, mm H ₂ O/cm ²	<5.0
Flame Spread	Class 1

*Tested and certified in an independent laboratory

SHIPPING DETAILS:

- Eaches per box: 50 eaches per box
- Eaches per case: 2000 eaches per case
- Weight per case: 20lbs
- Dims per case: 20" x 16" x 16"
- # of Cases per pallet: 20
- Weight per pallet: 445lbs
- Dims per pallet: 48" x 42" x 68"



medava®
medavausa.com